


# OM-A5-01-1229

## Revision history

Rev.	Date	Edited by	Sections affected	Approved by
0.1.1	2008-OCT-07	Csaba Magyar	ALL – initial revision	
1.0.17	2009-FEB-25	Csaba Magyar	1.0.17 – closing revision	Csaba Magyar
1.0.18	2009-MAR-28	Miklos Svarcz	1.0.18-finalised version	Svarcz Miklós
1.0.20	2009-APR-15	Csaba Magyar	Troubleshooting	
1.0.21	2009-MAY-25	Csaba Magyar		
1.1.1	2009-JUN-10	C.Magyar	SW revision - update	
1.2.1	2009-JUN-11	C.Magyar	Reorder	
1.2.2.	2009.JUN-12	C. Magyar	Accommodate changes and review reordering	
1.2.3.	2009.JUN.15	C. Magyar	Edit	
1.2.5	2009-JUN-30	C.Magyar	Limited features removed	
1.2.6	2009-AUG-14	C.Magyar	Menu system review	
1.2.7	2009-SEP-16	C.Magyar	List of accessories	
1.2.8	2009-OCT-14	C.Magyar	SV cleaning (11.5.1), Result interpretation (sect 6), Reag.cons update,	
			Fluidic schematics (14.1-3); Measurement principles detailed;	
1.2.9	2009-OCT-29	C.MAGYAR	Xb description; measurement methods	
1.2.10				
1.2.11	2009-12-12	C.Magyar	Autoloder prep.	
1.2.12	2010-01-18 - 2010-02-17	Z.Borbás	Rewording, parameter clean-up, screen-shots, pictures, drawings, re-structuring, Autosampler, reagent lock	
1.2.13	2010-02-24.	Z.Borbás	Improved flagging description	
1.2.14	2010-03-03	C.Magyar	Troubleshooting added	
1.2.15	2010-04-12	C.Magyar	Update to SW 1.0.339	
1.2.16	2010-04-15	C:Magyar	XB related functions removed	
1.2.17	2010-06-04	C.Magyar	Screenshot updates, Measurement related troubleshooting temporarily removed	
1.2.18	2010-06-11	C.Magyar	SV cleaning revised	
1.2.19	2010-06-29	C. Magyar	Reagent filtering method (specs) ultra→micro filtered	
1.2.20	2011-02-14	C.Magyar	X-B functions reinserted. Based on review by David Lopez; based on 2 <sup>nd</sup> review by GyA; Vienna sales office removed; installation by service → service manual; load qc values added X-B explanations included along with screenshots; sw versions updated in the top; “set date and time” removed from installation section; reagent lock section rephrased; troubleshooting part rearranged; reagent consumption updated; printout samples updated; EMC Class A specific information added to section 1.1.3; Reagent replacement procedure: steps to replace containers added: section 15.1.11; meas related troubleshooting received hard cleaning; 1 <sup>st</sup> (approval page) added; garbage can icon added (section 1.1.3)	Á. Gyetvai A. Galgóczi
1.2.21	2011-03-10	C.Magyar	New SV cleaning added to SV cleaning chapter	
1.2.22	2011-04-29	A.Tremmel	Update screens; to match FDA release	
1.2.23	2011-07-06	D. Lopez de Quintana	Intensive review and update for US release.	
1.2.24	2011-07-14	D. Lopez de Quintana	Added section for Bar Code Symbologies, updated TAB file format, changed recommendation for shear valve cleaning from 300 to 1200, updated External Devices settings to add LIS TCP/IP “listener” port information	
1.2.25	2011.08.31	D. Lopez de Quintana	Added Indications for Use statement Qualified sample handling to include only venous blood Updated performance characteristics Removed references to refrigerated samples Added table of background limits	
1.2.26	2011.12.31	D. Lopez de Quintana	Updated analytical measuring range	
1.2.27	2012.02.11	D. Lopez de Quintana	Updated Indications for Use and removed references to Research Use Only (RUO) parameters	
1.2.28	2012.03.07	M. Switzer	Correct Blank Limit and AMR errors	
1.2.29	2012.11.22	T. Jozsi	General Update, reagent consumption update	

APPROVALS

Role	Name	Date	Signature
Author	Tünde Józsi		
R&D Director	Gábor Farkas		
Quality Systems	Mike Switzer		



**Abacus 5**

**Innovative Solutions in Hematology  
Operator's Manual**

Revision 1.2.29

Issued: November 22, 2012





This user manual is intended to give detailed information for end users of the Diatron 'Abacus 5' optical hematology analyzer.

Descriptions contained herein are relevant to 'Abacus 5', software versions:

- Firmware: 2.66 (1125)
- PIC sw: 2.3
- Laser sw: 3.6
- TCU sw: 3.52

#### Warranty

THIS WARRANTY IS EXCLUSIVE AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE.

## Exemptions

Diatron's obligation or liability under this warranty does not include any transportation or other charges or liability for direct, indirect or consequential damages or delay resulting from the improper use or application of the product or the use of parts or accessories not approved by Diatron or repairs by people other than Diatron authorized personnel.

This warranty shall not extend to:

- Any Diatron product which has been subjected to misuse, negligence or accident.
- Any Diatron product from which Diatron's original serial number tag or product identification markings have been altered or removed.
- Any product of any other manufacturer.

#### Company Contact

Manufacturer DIATRON MI PLC

Address Táblás u. 39, H-10978 Budapest, HUNGARY

Tel: +36-1-436 9800

Fax: +36-1-436 9809

[www.diatron.com](http://www.diatron.com)

[support@diatron.com](mailto:support@diatron.com)





## Table of Contents

<b>1</b>	<b><i>Indications for Use</i></b>	<b>12</b>
<b>2</b>	<b><i>For Your Safety</i></b>	<b>12</b>
2.1	Who Should Use This Manual	12
2.2	Special Symbols Used In This Manual	12
2.3	General Precautions	12
2.4	Environmental Factors	14
2.5	Electrical Requirements	14
2.6	Space Requirements	14
2.7	Weight Requirements	16
2.8	Waste Disposal	17
2.9	Known Limitations	17
2.10	Emergency Situations	17
<b>3</b>	<b><i>Product Support</i></b>	<b>18</b>
<b>4</b>	<b><i>Installation</i></b>	<b>20</b>
4.1	The 'Abacus 5' Analyzer Package Contents	20
4.2	The Autosampler Package Contenta	20
<b>5</b>	<b><i>General Overview and Principles of Operation</i></b>	<b>22</b>
5.1	General Overview	22
5.1.1	Measured Parameters	22
5.1.2	Approved Reagents and Control Materials	23
5.1.3	Genuine Diatron Reagents	24
5.2	Principles of Operation	26
5.2.1	Volumetric Impedance Method	26
5.2.2	Photometric Light Absorbance Method	26
5.2.3	Optical Light Scatter and Diffraction Method	27
5.3	Inside the 'Abacus 5' Analyzer	29
5.3.1	Front Panel	29
5.3.2	Back Panel	30
5.3.3	Left Side Assembly	32
5.3.4	Right Side Assembly	33
5.3.5	Front Assembly	34
<b>6</b>	<b><i>User Interface</i></b>	<b>36</b>
6.1	Using the Touch Screen	36
6.2	Using the 'Start' Button	36



6.3	Using an External Mouse .....	36
6.4	Using an External Keyboard .....	36
6.5	Using the On-Screen Virtual Keypads .....	37
6.6	Using the Graphical User Interface .....	38
6.6.1	Quick Links .....	39
6.6.2	The Interactive Display Area .....	39
6.6.3	The Status Display .....	39
6.6.4	Entering Information.....	40
6.7	The Menu System.....	40
6.7.1	Primary Menu Items .....	40
6.7.2	Starting a Manual Single Tube Measurement.....	40
6.7.3	Start Automated Measurements .....	41
6.7.4	Access the Database .....	41
6.7.5	Initiate Printing .....	41
6.7.6	Main Menu.....	41
6.7.7	Autosampler Control Panel.....	42
6.7.8	Adjust the Time and Date .....	42
6.7.9	Open the Warning Panel.....	42
6.7.10	Menu Tree.....	42
6.7.11	Safety Access Codes .....	45
7	<b>Start Up and Shut Down the 'Abacus 5'</b> .....	46
7.1	Start Up and Shut Down Overview .....	46
7.2	Starting Up the 'Abacus 5' Analyzer .....	46
7.2.1	Visual Inspection .....	46
7.2.2	Power Up the 'Abacus 5' Analyzer .....	46
7.2.3	Start Up the User Interface .....	47
7.2.4	User Logon .....	48
7.2.5	Pneumatic System Start and Blank Measurement .....	48
7.3	Exiting the 'Abacus 5' Analyzer .....	49
7.3.1	Log Off.....	49
7.3.2	Shut Down.....	50
7.3.3	Prepare for Shipment.....	50
7.3.4	Emergency Shut Down.....	52
7.3.5	Repackaging the 'Abacus 5' Analyzer.....	53
8	<b>Sample Measurement</b> .....	54
8.1	Sample Types Supported by the 'Abacus 5' .....	54
8.1.1	Supported Sample Tubes Types .....	54
8.1.2	Sampling Depth .....	54
8.1.3	Open or Closed Sample Tubes .....	55
8.1.4	Sample Collection and Handling .....	55
8.2	Sample Types and Modes .....	56
8.3	Sample Identification.....	56

<b>8.4</b>	<b>Running Samples.....</b>	<b>57</b>
8.4.1	Manual Mode .....	57
8.4.2	Automatic Mode.....	59
<b>8.5</b>	<b>Result Display .....</b>	<b>69</b>
<b>8.6</b>	<b>Printing Reports .....</b>	<b>70</b>
<b>8.7</b>	<b>The Measurement Process.....</b>	<b>70</b>
<b>9</b>	<b>Result Interpretation.....</b>	<b>72</b>
9.1	The Result Screen .....	72
9.2	Sample Identification Information .....	72
9.3	Parameter Information.....	73
9.3.1	Scatter Diagrams and Histograms .....	74
9.3.2	Warnings.....	76
<b>10</b>	<b>Database Functions.....</b>	<b>82</b>
10.1	Database Overview .....	82
10.2	Scrolling the Database View .....	83
10.3	Sorting Database Information.....	83
10.4	Manual Selection of Database Records .....	84
10.5	Automatic Selection of Database Records .....	84
10.6	Viewing Detailed Results .....	84
10.7	Statistics .....	85
10.8	Managing Database Records.....	85
10.8.1	Select By .....	85
10.8.2	Importing.....	86
10.8.3	Export .....	87
10.8.4	Send to LIS .....	87
10.8.5	Save Tab File .....	87
10.8.6	Save Raw Data .....	88
10.8.7	Delete .....	88
<b>11</b>	<b>Calibration .....</b>	<b>90</b>
11.1	Calibrating the 'Abacus 5' .....	90
11.1.1	View Calibrations .....	92
<b>12</b>	<b>Quality Control.....</b>	<b>94</b>
12.1	Set QC Reference .....	94
12.2	QC Measure .....	96
12.3	View QC References .....	96
12.4	View QC Data.....	96

12.5	View QC Diagrams .....	96
12.6	View X-B Data .....	96
12.7	View X-B Diagrams .....	97
<b>13</b>	<b><i>Patients</i> .....</b>	<b>98</b>
<b>14</b>	<b><i>Multi-User Mode</i>.....</b>	<b>100</b>
14.1	Types of Users .....	100
14.2	Creating a New User .....	100
14.3	Deleting or Changing Users .....	102
<b>15</b>	<b><i>Settings</i> .....</b>	<b>104</b>
15.1	Customize Settings .....	104
15.2	Laboratory Settings .....	105
15.3	External Devices .....	105
15.4	System Settings .....	106
15.5	Units .....	108
15.6	Printer Settings .....	108
15.7	Profile Limits Settings .....	109
15.8	X-B Settings.....	110
15.9	User Settings .....	111
15.10	Date and Time Adjustment .....	111
<b>16</b>	<b><i>Instrument Diagnostics</i> .....</b>	<b>114</b>
16.1	Analyzer Self Test .....	114
16.2	Log.....	116
16.3	Reagent Status .....	117
16.4	Statistics .....	118
16.5	Information.....	118
<b>17</b>	<b><i>Maintenance</i> .....</b>	<b>120</b>
17.1	Opening the Front Panel .....	120
17.2	Closing the Front Panel .....	120
17.3	Removing the Side Panels.....	121
17.4	User Maintainable Parts of the Analyzer .....	121
17.5	Software Maintenance Functions.....	122
17.6	Cleaning the Shear Valve .....	123

17.7	Cleaning the Washing Head .....	126
17.8	Software Cleaning Maintenance Functions.....	127
17.8.1	Daily Cleaning .....	128
17.8.2	Cleaning the Measurement System.....	128
17.8.3	Extended Cleaning (Hard Cleaning) .....	128
17.9	Replacing Reagents .....	129
18	<b>Reagent Locking .....</b>	<b>132</b>
19	<b>The Daily Routine .....</b>	<b>133</b>
20	<b>Troubleshooting .....</b>	<b>136</b>
20.1	Software error messages .....	136
20.2	Pneumatic error messages.....	136
20.3	Mechanical Problems .....	137
20.3.1	Sample Rotor (SR) Failures .....	137
20.3.2	Needle Mechanics, Vertical Motor (Mvert) Problems.....	137
20.3.3	Shear Valve (SV) Related Errors .....	138
20.3.4	Dilutor Errors .....	138
20.3.5	Priming Problems.....	139
20.3.6	Electronics Related Problems .....	139
20.3.7	The Analyzer Does Not Power On.....	139
20.3.8	I <sup>2</sup> C Errors Displayed At Startup .....	139
20.4	Measurement results related problems .....	139
20.4.1	Fluctuating PLT background values .....	139
20.4.2	Long, smeared population .....	140
21	<b>Accessories.....</b>	<b>142</b>
22	<b>Appendix.....</b>	<b>142</b>
22.1	Reagent Consumption .....	142
22.2	Display Ranges .....	143
22.3	Fluidic System .....	144
22.4	Printed Report Formats .....	145
22.5	Parameter Calculation .....	147
22.6	Specifications.....	149
22.7	Performance Characteristics .....	151
22.7.1	Precision .....	151
22.7.2	Accuracy.....	151
22.7.3	Linearity .....	151
22.7.4	Carryover .....	152
22.7.5	Sample Stability .....	152
22.7.6	Mode to Mode.....	152
22.7.7	Reference Ranges .....	152

22.7.8	Background Limits.....	153
22.7.9	Analytical Measurement Range.....	153
22.7.10	Interfering Substances.....	154
<b>22.8</b>	<b>Reagent System.....</b>	<b>155</b>
22.8.1	Diluent.....	155
22.8.2	Lyse Reagent 1.....	155
22.8.3	Lyse Reagent 2.....	155
22.8.4	Cleaner.....	155
<b>22.9</b>	<b>Tab File Format .....</b>	<b>156</b>
<b>23</b>	<b>Index .....</b>	<b>157</b>

## Table of Figures

Figure 1.	Abacus 5 with Autosampler Space Requirements.....	15
Figure 2.	Abacus 5 Without Autosampler Space Requirements.....	16
Figure 3.	Volumetric Impedance Method.....	26
Figure 4.	Photometric Light Absorbance Method .....	27
Figure 5.	'Abacus 5' Optical Head Block Diagram.....	27
Figure 6.	Optical Signal Processing System.....	28
Figure 7.	Cellular Light Scatter .....	28
Figure 8.	4DIFF Scatter Diagram .....	29
Figure 9.	Abacus 5 Front Panel .....	29
Figure 10.	Abacus 5 Back Panel .....	30
Figure 11.	Main Board Back Panel I/O Ports.....	31
Figure 12.	Left Side Assembly .....	32
Figure 13.	Left Side Assembly .....	33
Figure 14.	Front Assembly Behind the Front Panel .....	34
Figure 15.	Alphanumeric Virtual On-Screen Keypad .....	37
Figure 16.	'Sticky' Shift and Symbol Buttons .....	37
Figure 17.	Numeric Virtual On-screen Keypad .....	38
Figure 18.	Date Virtual On-Screen Keypad .....	38
Figure 19.	Graphical User Interface Sections.....	39
Figure 20.	The Status Display.....	40
Figure 21.	'Abacus 5' Main Menu .....	42
Figure 22.	'Abacus 5' Analyzer Startup and Log On Screens .....	47
Figure 23.	Abacus 5 Shutdown Options.....	49
Figure 24.	Log Off.....	49
Figure 25.	Shut Down.....	50
Figure 26.	Prepare for Shipment .....	50
Figure 27.	Connect Distilled Water Panel.....	51
Figure 28.	Remove Reagent Tubing Panel .....	52
Figure 29.	Shut Down Ready Panel.....	52
Figure 30.	Sampling Needle .....	55

Figure 31. Sample Measurement Panel .....	58
Figure 32. Sample Processing Result Screen .....	59
Figure 33. Autosampler Info and AS Panel .....	60
Figure 34. Full Scan Automatic Processing Mode .....	62
Figure 35. Full Scan Mode Tray View Progress .....	63
Figure 36. Full Scan Mode Tray and List Views .....	64
Figure 37. Free List Mode Selection .....	65
Figure 38 Figure 39. Preparing a Free List.....	66
Figure 40. Free List Mode Progress: List View .....	66
Figure 41. Free List Mode Progress: Tray View .....	67
Figure 42. Selected Sample Mode Panel.....	68
Figure 43. Control The Autosampler With Info Panel.....	69
Figure 44. Results Display and Magnified Scatter Plot .....	70
Figure 45. 'Abacus 5' Measurement Process .....	71
Figure 46. 'Abacus 5' Result Screen .....	72
Figure 47. Parameter Information Display.....	73
Figure 48. Graphical Normal Range Display.....	74
Figure 49. Result Screen Scatter Diagrams .....	75
Figure 50. Result Screen Histograms .....	75
Figure 51. Warnings Section of Results Screen.....	76
Figure 52. Accessing the Database.....	82
Figure 53. Scrolling and Selecting .....	83
Figure 54. Multiselect and Multiple Selection .....	84
Figure 55. Manage Records Panel.....	85
Figure 56. Select By Panel .....	85
Figure 57. Database Importing Panel.....	86
Figure 58. Directory Panel for Data Storage .....	87
Figure 59. Calibration Options .....	90
Figure 60. Calibration Panel and Calibration Run Result Panel .....	91
Figure 61. Calibration Factor Panel.....	92
Figure 62. View Calibration Panel .....	93
Figure 63. QC Panel.....	94
Figure 64. Set QC Reference Panel .....	95
Figure 65. Load QC Reference Panel.....	95
Figure 66. View QC Diagrams Panel.....	96
Figure 67. Patients Panel .....	98
Figure 68. Edit Patient Panel.....	99
Figure 69. Logon Panel.....	101
Figure 70. Add User Panel.....	101
Figure 71. Add User Complete Panel .....	101
Figure 72. Settings Panel (Administrator User).....	102
Figure 73. Settings Panel (Administrator User).....	103
Figure 74. Confirm User Delete Panel.....	103
Figure 75. Settings Panel.....	104
Figure 76. Customize Settings Panel.....	105

Figure 77. External Devices Settings Panel .....	106
Figure 78. System Settings Panel .....	107
Figure 79. Units Settings Panel .....	108
Figure 80. Printer Settings Panel .....	109
Figure 81. Profile Limits Settings Panel.....	110
Figure 82. X-B Settings Panel .....	111
Figure 83. User Settings Panel .....	111
Figure 84. Time/Date Settings .....	112
Figure 85. Diagnostics Panel .....	114
Figure 86. Self Test Panel.....	115
Figure 87. Log Panel.....	116
Figure 88. Reagent Status Panel .....	117
Figure 89. Statistics Panel .....	118
Figure 90. Information Panel .....	119
Figure 91. Opening the Front Panel.....	120
Figure 92. Closing the Front Panel.....	121
Figure 93. Software Maintenance Functions.....	122
Figure 94. Shearvalve Cleaning Maintenance Function.....	124
Figure 95. Shearvalve Cleaning Maintenance Function Prompt .....	124
Figure 96. Removing the Shear Valve Axis Thumbscrew .....	125
Figure 97. Disassembling the shear valve.....	125
Figure 98. Cleaning the Shear Valve .....	126
Figure 99. Reassembling the Shear Valve.....	126
Figure 100. Removing the Washing Head.....	127
Figure 101. Cleaning the Washing Head.....	127
Figure 102. Reagent Low Indicator and Explanation Panel .....	129
Figure 103. Reagent Replacement Panel.....	130
Figure 104. Location of Reagent Key in Diatro•Lyse-5P Reagent Container .....	132
Figure 105. Reagent Key Connector .....	132
Figure 106. Pneumactical Error Message Example .....	136
Figure 107. Detailed Error Message .....	137
Figure 108. Proper Seating of Shear Valve Thumbscrew.....	138
Figure 109. Smeared Population .....	140
Figure 110. Fluidic System Diagram.....	144
Figure 111. Detailed Printout Format.....	145
Figure 112. Listed Printout Format.....	146

## Table of Tables

Table 1. 'Abacus 5' Analyzer Parameters.....	23
Table 2. Diatro•Dil-5P Diluent Information .....	24
Table 3. Diatro•Lyse-5P Information .....	24
Table 4. Diatro•Diff-5P Information .....	25
Table 5. Diatro•Hypocleaner CC Information.....	25
Table 6. Reagent Connector Color Codes .....	31

Table 7. 'Abacus 5' Analyzer Menu Tree.....	45
Table 8. Supported Autosampler Bar Code Reader Symbolologies.....	61
Table 9. Sample Identification Information .....	73
Table 10. Normal Range Flags.....	74
Table 11. Linearity Range Flag .....	74
Table 12. High Blank Flag .....	74
Table 13. Warning Flags .....	80
Table 14. Interpretive Flags .....	81
Table 15. Morphological Flags .....	81
Table 16. Database Sort Criteria .....	83
Table 17. Restrictions by User Type.....	100
Table 18. Count Unit Parameters.....	108
Table 19. Self Test Item Limits .....	116
Table 20. Log Filters .....	117
Table 21. Software Maintenance Procedures.....	123
Table 22. Reagent Consumption .....	142
Table 23. Display Ranges.....	143
Table 24. Parameter Calculation.....	148
Table 25. Specifications.....	150
Table 26. Precision Performance .....	151
Table 27. Accuracy Performance .....	151
Table 28. Linearity Performance.....	151
Table 29. Carryover Performance .....	152
Table 30. Closed vs. Open Vial Mode Performance.....	152
Table 31. Closed vs. Autosampler Mode Performance.....	152
Table 32. Reference Ranges.....	153
Table 33. Background Limits .....	153
Table 34. Analytical Measuring Range .....	154





## 1 Indications for Use

The Diatron Abacus 5 System is a quantitative multi-parameter automated hematology analyzer designed for in-vitro-diagnostic use in clinical laboratories for enumeration of the following parameters in K<sub>3</sub>EDTA anti-coagulated venous whole blood samples: WBC, LYM%, LYM#, MON%, MON#, NEU%, NEU#, EOS%, EOS#, BAS%, BAS#, RBC, HGB, HCT, MCV, MCH, MCHC, RDWcv, RDWsd, PCT, PDWcv, PDWsd, PLT and MPV. The Diatron Abacus 5 is indicated for use to identify patients with hematologic parameters within and outside of established reference ranges.

## 2 For Your Safety





### 2.1 Who Should Use This Manual

This user manual is intended for clinical laboratory professionals using the Diatron 'Abacus 5' automated hematology analyzer. The manual includes information about the operation and user interface of the 'Abacus 5' analyzer.





It also contains basic steps necessary to perform the setup procedures to customize the operation of the analyzer to the requirements of your laboratory.

















This manual also describes daily routine maintenance requirements to keep your analyzer functioning properly.

### 2.2 Special Symbols Used In This Manual

Label	Meaning	Explanation
	Biohazard	Blood samples and analyzer waste are potentially infectious materials.
	Corrosive	Reagents may cause corrosion or skin irritation.
	Warning	General warning of possible hazard conditions.
	Sharp needle warning	The sampling needle may be a hazard to the operator.

### 2.3 General Precautions

	The sampling needle and other components inside the analyzer may cause injury, or can get damaged if handled incorrectly. Only certified personnel should open the covers. Running measurements with opened cover is not recommended due to the risk of possible injury. Always wear safety gloves while performing maintenance actions.
	The analyzer weighs 35kg (~77lbs). Please do not attempt to move it alone. The analyzer should always be moved by two persons holding the analyzer by its sides in an upright position.
	Always use safe lifting procedures when lifting the analyzer.
	Make sure to retain the original packaging material for safe transportation and storage in the future.

	To prepare the analyzer for shipping, storage or extended periods of inactivity, please drain the reagents and repackage the 'Abacus 5' in its original packaging. Do not expose the 'Abacus 5' to direct sunlight, extreme temperature or humidity (>80%).
 	The analyzer operates with chemically and biologically active reagents. Physical contact with these reagents should be avoided. Please read reagent descriptions carefully for possible emergency actions.
	To ensure reliable operation and reliable results: <ul style="list-style-type: none"> <li>• Only human blood samples should be analyzed</li> <li>• Only genuine Diatron reagents should be used</li> <li>• Required maintenance (user and service level) should be performed as recommended in this manual</li> <li>• Only Diatron certified service personnel should perform service actions</li> <li>• Only genuine Diatron service materials and spare parts should be used</li> </ul>
	Genuine reagents and service materials and spare parts are available from Diatron.
	Only Diatron certified service personnel that have successfully completed the 'Diatron Abacus 5 Service Training' program are qualified to service the 'Abacus 5' analyzer.
	Before operating the 'Abacus 5' analyzer, all operators should complete a 'Diatron Abacus 5 Operator Training' program. This program is offered by Diatron or by Diatron certified service personnel.
 	Replacement materials or spare parts (tubes, valves, etc.) which might have been in contact with human blood or reagents should be handled as a potentially biologically hazardous and chemically dangerous material. All applicable laws and regulations must be observed in the handling and disposal of these materials.
	The 'Abacus 5' is designed for laboratory operation. Mobile operation is not supported. Operate 'Abacus 5' within the ambient temperature range described in section 2.4.
	The IVD equipment complies with the emission and immunity requirements described in relevant part of the IEC 61326 series.
	This equipment has been designed and tested to CISPR 11 Class A. In a domestic environment it may cause radio interference, in which case, you may need to take measures to mitigate the interference.
	Electromagnetic environment should be evaluated prior to operation of the device.
 	This analyzer contains electronic components. Please handle electronic waste adhering to local or federal regulations.
	CAUTION – Use of controls or adjustments or performance of procedures other than those specified herein may result hazardous radiation exposure.

## 2.4 Environmental Factors

Operate the 'Abacus 5' analyzer within the ambient temperature range of 15-30°C (59-86 °F) and a relative humidity range of 10% - 80%. The optimum operating temperature is 25°C (~77°F).

The 'Abacus 5' analyzer should be stored within the temperature range of 5-35°C (41-95 °F). Avoid exposing the 'Abacus 5' analyzer to direct sunlight or to extreme high or low temperatures. If the 'Abacus 5' analyzer was subjected to extreme temperatures during shipment or storage, the analyzer must be placed for at least one hour in a room whose temperature is within the operational range before installation or use.

Reagents should be stored at a temperature range of 15-30°C (59-86 °F). Reagents may experience a temperature range of at most 5-35°C (41-95 °F) for a maximum of 3 days.

The analyzer should be placed in a well-ventilated location.

Operation at an altitude above 3000 meters (9800 ft) is not recommended.

## 2.5 Electrical Requirements

The analyzer should only be operated from a wall outlet meeting these power input requirements:

- 100-127VAC/200-240VAC; 47Hz to 63Hz
- Power Consumption: maximum 400 VA

Please ensure that the wall outlet is also capable of supplying the power consumption of any additional devices (such as a printer).

Use only the power cord supplied with the instrument. Avoid using extension cords. The 'Abacus 5' analyzer comes with a power cord appropriate for your power system. Proper use of the appropriate power cord assures adequate grounding of the system. If the power is not reliable, contact your representative for options such as the installation of an external UPS module.



Failure to properly ground the 'Abacus 5' bypasses important safety features and may result in electrical hazard.

The instrument should not be placed near potentially interfering devices capable of emitting radio frequencies (e.g. radio or television transmitters/receivers, radars, centrifuges, X-ray devices, fans, etc.).

This analyzer is designed to be safe for transient voltages to INSTALLATION CATEGORY II and POLLUTION DEGREE 2.

## 2.6 Space Requirements

It is important to install the instrument in a suitable location. A poor location can adversely affect its performance.

Select a well-ventilated location near a power source and close to a suitable drain.

Place the unit on a clean, level surface. Leave at least 0.5 m (18 inches) space on both sides and above the instrument to access pneumatics. A minimum of 0.2 m (8 inches) must be maintained between the rear panel and the wall to allow for heat dissipation and tubing clearance.

Ensure there is enough clearance in front of the 'Abacus 5' analyzer to open the front panel. Allow enough space if you want to use optional external keyboard, mouse or bar code reader.

Your selected location should allow placement of the reagents in an unobtrusive location below the laboratory bench that the instrument is placed on, or on the same surface. Placement below the laboratory bench also allows for storage of a spare set of reagents. Never place the reagents above the 'Abacus 5' analyzer.



Placing reagents above the 'Abacus 5' analyzer could result in reagent overflow and spilling.

See Figure 1 and Figure 2 for more information about proper analyzer location and clearance.




Figure 1. Abacus 5 with Autosampler Space Requirements






Figure 2. Abacus 5 Without Autosampler Space Requirements

## 2.7 Weight Requirements


The 'Abacus 5' analyzer weighs 35 Kg (77 lb) without the Autosampler. The 'Abacus 5' with optional Autosampler weighs 47Kg (104 lb). Adding an external keyboard, documents etc. can bring the total weight up to 60 Kg (132 lb). If you decide to store the reagents on the same surface, then the combined weight can reach 100 Kg (220 lb).

Please select a table, laboratory shelf, or other location which can support the weight of the 'Abacus 5' with accessories and is free from vibration.

 	<p>To allow reliable operation and to provide a safe working environment, make sure that the table supporting the unit is stable enough to carry the weight of the instrument and accessories.</p>
--	--

## 2.8 Waste Disposal

'Abacus 5' analyzer waste contains human blood and reagents that are chemically and biologically active, and should be considered to be a potential infection and biohazard threat. Safe laboratory practices must be followed including the use of personal protective when operating the 'Abacus 5' and handling blood, reagents, and waste.

	<p>System waste should be handled as a potentially biohazardous material. All applicable laws, regulations, and laboratory practices should be followed in the handling and disposal of waste.</p>
---	--

## 2.9 Known Limitations


The 'Abacus 5' is not intended for analysis of animal blood samples. Anti-coagulated and homogenized (well-mixed) human blood samples must be free from contamination.

Blood samples must be analyzed within 7 hours of venipuncture.

## 2.10 Emergency Situations

Always follow all applicable laws and regulations with regard to emergency situations.

If the 'Abacus 5' needs to be powered off due to an emergency situation (like fire, thunderstorm etc.), follow the procedures in chapter 7.3.4.

	<p>In case of fire, do not use water to extinguish the fire unless the 'Abacus 5' is disconnected from the electrical network!</p>
---	--

### **3 Product Support**

Your sales representative is trained and qualified to answer questions about applications and the operation of the 'Abacus 5' analyzer. If you have additional questions regarding the 'Abacus 5', please ask for manufacturer support through your sales representative.

Diatron makes every effort to provide excellent support and up to date information and services regarding the 'Abacus 5' analyzer. Software updates and application documents are available on the Diatron web site at

<http://www.diatron.com>.





## 4 Installation

Your Diatron certified service engineer will perform the initial installation of the 'Abacus 5' analyzer and train selected laboratory personnel on the proper use and maintenance of the analyzer. The 'Abacus 5' analyzer should only be operated by properly trained laboratory personnel.

### 4.1 The 'Abacus 5' Analyzer Package Contents

After opening the 'Abacus 5' analyzer packaging, you will find an accessory box. The analyzer is in a plastic bag in between cushioning protective material. After opening the bag, remove the desiccant packs. The package should contain the following items:

- The 'Abacus 5' analyzer
- Accessories box on top of the analyzer, including:
  - Power cord
  - User Manual (this document)
  - Sample tube adapter for the individual sample unit (1 pc)
  - Reagent tube kit in a plastic bag, with 3 caps and 5 color-coded tubing
  - Plastic waste container
  - Cardboard box for the waste container
  - "Preparing for shipment" tube kit (1 piece)
  - 'Abacus 5' analyzer final check report

Please inform your Diatron Representative if any item is missing or damaged. Retain the original 'Abacus 5' analyzer packaging for future transportation and storage.

### 4.2 The Autosampler Package Content

The Autosampler is an optional unit that attaches to the Abacus 5. If the optional Autosampler unit was ordered, package should contain the following items:

- Autosampler unit
- A sample tray
- 10 pieces sample racks
- Autosampler final check report

Please inform your sales representative if any item is missing or damaged. Retain the original Autosampler packaging for future transportation and storage.



## 5 General Overview and Principles of Operation

### 5.1 General Overview

The 'Abacus 5' is a fully automated high quality hematology analyzer for *in vitro* diagnostic use in clinical laboratories. It provides precise and accurate 5-part differential measurement using a laser based optical measurement technology.

The 'Abacus 5' analyzer uses the impedance method for measurement of leukocytes (WBC), erythrocytes (RBC) and platelet (PLT) concentrations. Measurement of the hemoglobin (HGB) content of red blood cells is accomplished by photometric measurement technology. Five part leukocyte differential (LYM%, MON%, NEU%, EOS%, BAS%) is accomplished using optical laser-based flow cytometric technology.

A vivid color touch screen display is featured with an intuitive, informative, and attractive user interface. A START button allows one-touch operation for ease of use.

The 'Abacus 5' analyzer's unique software system supports the use of many commonly used external printers with its USB connections. The 'Abacus 5' internal database is capable of storing 100,000 patient, QC, and calibration result records including flags and graphical scatter diagrams and histograms. The system software is field-upgradeable to ensure up-to-date operation.

An automatic optional Autosampler is available (sold separately) for automated processing of up to 100 sample tubes for increased laboratory efficiency.

The 'Abacus 5' features advanced Ethernet LIS connectivity using the HL7 protocol in addition to standard serial interfaces, providing the clinical laboratory with flexible connectivity options.

#### 5.1.1 Measured Parameters

'Abacus 5' determines 21 hematology parameters including 5-part WBC differential. The instrument requires 110  $\mu$ l of whole blood-sample in Closed- and Open-vial mode. Cycle time is 60 seconds.

Parameter	Symbol	Measurement method
White Blood Cell count	WBC	Impedance measurement
Neutrophil absolute count	NEU	Calculated
Lymphocyte absolute count	LYM	Calculated
Monocyte absolute count	MON	Calculated
Eosinophil absolute count	EOS	Calculated
Basophil absolute count	BAS	Calculated
Neutrophil percentage	NEU%	Optical measurement
Lymphocyte percentage	LYM%	Optical measurement
Monocyte percentage	MON%	Optical measurement
Eosinophil percentage	EOS%	Optical measurement
Basophil percentage	BAS%	Optical measurement
Red Blood Cell count	RBC	Impedance measurement
Hemoglobin	HGB	Photometric measurement
Hematocrit	HCT	Calculated
Mean Corpuscular Volume	MCV	Derived
Mean Corpuscular Hemoglobin	MCH	Calculated
Mean Corpuscular Hemoglobin Concentration	MCHC	Calculated
Red Blood cells Distribution Width	RDWsd	Calculated

Parameter	Symbol	Measurement method
Red Blood cells Distribution Width	RDWcv	Derived
Platelet count	PLT	Impedance measurement
Mean Platelet Volume	MPV	Derived
Platelet Distribution Width	PDWsd	Calculated
Platelet Distribution Width	PDWcv	Calculated
Thrombocrit	PCT	Calculated

Table 1. 'Abacus 5' Analyzer Parameters

### 5.1.2 Approved Reagents and Control Materials

Always use reagents, control materials and calibration materials that are recommended and approved by the manufacturer. The 'Abacus 5' analyzer, control material, calibrator material and reagents are part of a system and are carefully designed and selected for optimal performance. Using unapproved reagents or controls may cause false flagging or incorrect, inaccurate results.

Only genuine Diatron reagents may be used in the 'Abacus 5' analyzer. Diatron recommends the use of Diacon 5 control material and Diatrocacal calibration material from Diatron. These materials have been specifically assayed for use in the 'Abacus 5' analyzer.

For additional details about these materials, please contact the manufacturer directly.

### 5.1.3 Genuine Diatron Reagents

Only use the genuine Diatron reagents below to ensure proper performance of the 'Abacus' analyzer. Reagents are manufactured and provided by Diatron MI Plc., and are for *in vitro* diagnostic use only.



If the eyes or skin come into contact with any of the reagents, flush abundantly with water. If ingested, drink water immediately and seek medical attention.

#### 5.1.3.1 Diatro•Dil-5P Diluent

Description	Multiple micro-filtered, particle free buffered isotonic solution, containing stabilizers, special additives and preservatives.
Application	Quantitative and qualitative determination of RBC, WBC, PLT and HGB concentration.
Appearance	Colorless, odorless solution.
Shelf-life	36 months
Open bottle stability	120 days
Storage	Between +15 °C and +30 °C. (~59-86 °F)
Diatron reagent code	D151-X (X: packaging size; 2=20L)
Color code	Green

Table 2. Diatro•Dil-5P Diluent Information

The product is environmentally friendly, and does not contain azide nor cyanide

#### 5.1.3.2 Diatro•Lyse-5P

Description	Multiple micro-filtered, particle free reagent solution, containing lysing detergents, stabilizers, leucoprotective components, special additives and preservatives.
Application	Quantitative and qualitative determination of WBC, 5diff and HGB concentration.
Appearance	Colorless solution, foaming by shaking.
Shelf-life	24 months
Open bottle stability	120 days
Storage	Between +15 °C and +30 °C. (~59-86 °F)
Diatron reagent code	D301-X (X: packaging size 5=5L)
Color code	Yellow

Table 3. Diatro•Lyse-5P Information

The product is environmentally friendly, and does not contain azide nor cyanide

#### 5.1.3.3 Diatro•Diff-5P

Description	Multiple micro-filtered, particle free reagent solution, containing stabilizers, leucoprotective components, special additives and preservatives.
Application	Quantitative determination of WBC, leukocyte four-part

	differentiation (LYM, MON, NEU, EOS).
Appearance	Colorless, odorless solution.
Shelf-life	24 months
Open bottle stability	120 days
Storage	Between +15 °C and +30 °C. (~59-86 °F)
Diatron reagent code	D302-X (X: packaging size 1=1L)
Color code	Orange

Table 4. Diatro•Diff-5P Information

The product is environmentally friendly, and does not contain azide nor cyanide

#### ***5.1.3.4 Diatro•Hypocleaner CC***

Description	Containing alkaline hypochlorite, special additives and preservatives.
Application	Capillaries, tubing and chambers, removing blood component precipitates.
Appearance	Slightly yellow liquid with chlorine odor.
Shelf-life	15 months.
Open bottle stability	120 days
Storage	Between +15 °C and +30 °C (~59-86 °F)
Diatron reagent code	D801-X (X: packaging size 1=100mL)

Table 5. Diatro•Hypocleaner CC Information

The “Hypocleaner CC” is not an online-reagent (not directly connected to the ‘Abacus 5’). This reagent is aspirated from a sample tube into the ‘Abacus 5’ directly.

## 5.2 Principles of Operation

The 'Abacus 5' analyzer uses a combination of methods to provide measurement results:

- Volumetric impedance is used to determine the cellular concentrations and volume distributions of leukocytes (WBC), erythrocytes (RBC), and platelets (PLT).
- Photometric measurement of light absorbance is used to determine hemoglobin (HGB) concentration.
- Optical measurement of light scattering and diffraction is used to determine five part leukocyte (LYM%, MON%, NEU%, EOS%, BAS%) differential parameters.

### 5.2.1 Volumetric Impedance Method

The volumetric impedance method determines cellular concentrations and volume distributions of cells by detecting and measuring changes in electrical impedance when particles suspended in a conductive liquid pass through a small aperture. The method is "volumetric" because a small known volume of blood is precisely diluted with a conductive diluent and forced through the aperture at a fixed rate.




Figure 3. Volumetric Impedance Method

A constant direct current flows between the electrodes on either side of the aperture. Each cell passing through the aperture causes a change in the electrical impedance of the conductive blood cell suspension. This change is sensed by the 'Abacus 5' electronics and converted to an electrical pulse. The quantity of pulses is proportional to the number of particles. The intensity of each pulse is proportional to the volume of the particle. The volume distribution diagrams of the particles result in the WBC, RBC, and PLT histograms that measured in femtoliter volume units.

Electronic discrimination allows separation of erythrocytes (RBC) and platelets (PLT). A lytic reaction lyses erythrocytes to clearly measure leukocytes (WBC).

### 5.2.2 Photometric Light Absorbance Method

A lysed blood sample dilution can be analyzed for hemoglobin (HGB) concentration based on its stable chromogen content. The reagent lyses the red blood cells causing the release of cellular hemoglobin. The hemoglobin



concentration is measured by taking a photometric reading across the 'Abacus 5' WBC chamber. The HGB measurement is calculated as the difference between a blank and a sample measurement with and without illumination to reduce the effect of liquid refraction and incident light.




Figure 4. Photometric Light Absorbance Method

### 5.2.3 Optical Light Scatter and Diffraction Method

Optical measurement of light scattering and diffraction is used to determine five part leukocyte (LYM%, MON%, NEU%, EOS%, BAS%) differential parameters. An optical measuring head contains a focused laser source that is used to illuminate a stream of leukocytes (WBC) suspended in an optically clear diluent moving through a flow cell.

The cells scatter light as they flow through the path of the laser beam. An optical detector senses changes in the intensity of the scattered laser light which are proportional to the cell volume and granularity of the cell's internal structure. The 'Abacus 5' electronics convert these changes to electrical pulses which are gathered and stored for analysis. Five part population discrimination is based on analysis of the two dimensional volume and granularity distribution scatter diagram.




Figure 5. 'Abacus 5' Optical Head Block Diagram

Cells with greater volume or size or more granularity will tend to scatter greater amounts of light. The intensity of scattered light is detected by an optical signal processing system.




Figure 6. Optical Signal Processing System

External structure and the size of the cell cause lower angles of scatter, and internal granularity or complexity causes higher angles of diffraction. Both low and high angles of light scatter are captured by optical sensors, providing the 'Abacus 5' analyzer with two independent measurements for each cell crossing the path of the laser beam.




Figure 7. Cellular Light Scatter

The light scatter data are gathered and plotted as a two dimensional scatter diagram. Similar cells have similar scatter characteristics and tend to group together. This allows the analytical software to differentiate and identify the leukocyte populations and generate the 4DIFF and BASO scatter diagrams.




Figure 8. 4DIFF Scatter Diagram

### 5.3 Inside the 'Abacus 5' Analyzer

#### 5.3.1 Front Panel




Figure 9. Abacus 5 Front Panel

The touch screen display is the primary user interface of the 'Abacus 5' analyzer for user interactions and data entry. The LCD screen can be cleaned damp sponge and gentle detergent.

The front panel is a hinged cover that provides easy access to the shear valve, sampling needle and syringes for cleaning and maintenance. The front panel, START button, and the sample rotor can be cleaned with damp sponge and regular e.g. household cleaning material

Sample tubes are placed into a tube adapter inserted on the sample rotor for closed vial or open vial single-tube manual sample processing. The sample tube is rotated from the front of the unit into the inside of the 'Abacus 5' for aspiration. This safety feature protects the operator from inadvertent contact with the aspiration needle during specimen processing.

Automatic sample processing with the optional Autosampler can be interrupted in order to run one or more stat samples with the sample rotor. Automatic processing can be resumed after stat sample processing is complete.

The START button allows easy one-hand operation of the 'Abacus 5' analyzer. The color of the START button indicates the status of the 'Abacus 5' analyzer.

- Green – ready to process samples
- Red – busy
- Orange – standby

### 5.3.2 Back Panel



Figure 10. Abacus 5 Back Panel

①Power connector: the 'Abacus 5' power connector should be connected to a grounded power outlet that meets the requirements listed in section 2.5.

②Main power switch: this is a small switch located next to the power connector. This switch turns off all electrical power to the 'Abacus 5' analyzer. Leaving this switch on allows the 'Abacus 5' to remain in the standby state. The switch is on when it is in the 'up' position labeled '1' and is off when it is in the down position labeled '0'.

③Standby switch: flipping this switch to the 'up' position takes the 'Abacus 5' main computer out of the standby state and powers up internal analyzer components. Please note that this switch springs back to the down position when released and does not stay in the 'up' position.






	Green	1 connector	Diatro•Dil-DIFF diluent
	Orange	1 connector	Diatro•Diff-5p
	Yellow	1 connector	Diatro•Lyse-5p
	Red	2 connectors	Waste

Table 6. Reagent Connector Color Codes

④Hardware Key connector: this connection is used to plug in the hardware key supplied with every new package of genuine Diatron Diatro•Lyse-5p reagent. A new reagent key provides the analyzer with the ability to run 700 measurements. When the measurement count has been exhausted, the 'Abacus 5' analyzer stops processing samples.

⑤Reagent connectors: five color-coded tubing connectors provide connection to the reagents. There are three reagent connectors and two waste connectors.

⑥Main board back panel I/O Ports: these are the back panel computers of the 'Abacus 5' analyzer main computer board. They provide standard connection to peripherals such as external keyboard and mouse, printers, bar code readers.




1	PS2 Mouse port for optional mouse
1	PS2 Keyboard for optional PS2 keyboard
2	COM port for LIS
1	DVI-D port (unused)
1	VGA port (unused)
2	RJ45 Ethernet port for LIS
4	USB 2.0 ports for peripherals
3	Audio jacks (unused)

Figure 11. Main Board Back Panel I/O Ports

Avoid getting liquids on the electrical connectors. The rear panel of the analyzer should be cleaned with a cloth or sponge with alcohol.

⑦Serial # and manufacturing label:



### 5.3.3 Left Side Assembly

To access the left or right side assemblies, the side panels must be removed. Please refer to section 17.3 for information about removing the side panels.




Figure 12. Left Side Assembly

- Pressure sensors measure the differential pressures of vacuum buffer for impedance base and optical measurement and empty chamber process.
- Small Measure buffer is a vacuum buffer used for impedance measurement.
- Sheath buffer is a reservoir for Diatro•Dil-Diff diluent for optical measurements.



- Diff-5P buffer is a reservoir for Diatro•Diff-5P reagent.
- Lyse-5P buffer is a reservoir for Diatro•Lyse-5P reagent.
- Diluent buffers are reservoirs for Diatro•Dil-5P diluent.
- Valves 1-44 connect pneumatic components to carry out dilution, reaction, measurement, and waste disposal processes.
- Drain buffer controls chamber draining procedures.
- Mix chamber performs the first step of RBC dilution.
- RBC chamber performs the second step of RBC dilution and contains the RBC aperture.
- WBC chamber performs the WBC dilution and holds the WBC aperture and HGB photometric measurement system.

#### 5.3.4 Right Side Assembly




Figure 13. Left Side Assembly

- Aspirating needle aspirates a blood sample from the sample tube.
- Horizontal, vertical motion unit positions the aspiration needle.
- Washing head cleans the aspiration needle to prepare for the next sample.
- Sample rotor holds the sample tube for single tube manual presentation.
- Autosampler connector is the power and control connector of the optional Autosampler.
- Optical head contains the laser source, optical detectors and auto alignment mechanism for the optical measurement.

### 5.3.5 Front Assembly




Figure 14. Front Assembly Behind the Front Panel

- Syringe pumps control the volume and flow rate of reagents through the pneumatic components.
- Shear valve is an integrated multi-functional valve for segmenting precise volumes of blood.





## 6 User Interface

The user interface of the 'Abacus 5' analyzer is designed to be intuitive and easy to use, and provide access to all user level functions. You can interact with the 'Abacus 5' by using only the touch screen and the virtual on-screen keypad, or by using an external keyboard and mouse (PS2 or USB). The Start button is considered an important part of the user interface, allowing one-hand, one-touch processing of sample tubes.

### 6.1 Using the Touch Screen

The touch screen is a pressure sensitive overlay on the screen of the 'Abacus 5'. It allows selection of menu items and pressing of buttons by gently pressing or tapping the surface over the selected item. This method is called 'tapping'. The touch-screen device used in the 'Abacus 5' analyzer is selected for reliability, endurance and compatibility with the medical laboratory environment. Simple operations are supported, but multi-touch "gestures" or multiple tapping are not supported.

The touchable items on the 'Abacus 5' user interface are designed to be large enough to tap with finger. The use of a tapping tool may be useful if the operator has larger or thicker fingers, or when wearing oversize gloves. A tapping tool can be pencil eraser or a specially designed stylus available at computer or consumer electronics stores.

Do not use sharp or heavy objects to operate the touch screen. The surface material is moderately resistant to liquids. Avoid touching the screen with wet fingers, or allow liquids to come in contact with the screen or the surrounding area.

If you find that the 'Abacus 5' touch screen responds in a different area than the one you are tapping, then please perform a touch screen calibration procedure before calling for service or support. See chapter 17.5 for more details on touch screen calibration and other maintenance procedures.

### 6.2 Using the 'Start' Button

The 'Start' button is a physical button on the front cover of the 'Abacus 5'. It has the following functions:

- Start a single tube, closed or open vial measurement using the sample rotor by pressing the Start button.
- Provide color-coded feedback about the status of the measurement system.
- Green means that the 'Abacus 5' is ready to start a new measurement.
- Red means that the system is busy running a measurement or a procedure.
- Orange color means that the analyzer is in standby mode.

### 6.3 Using an External Mouse

An external mouse can be used to perform the same functions as the touch screen. Navigate the mouse pointer over the desired item and click the left mouse button instead of tapping the touch-screen to activate the function. This procedure is referenced as clicking (click) an item. Touch screen operation continues to be available even if an external mouse is connected to the 'Abacus 5' analyzer.

It is possible to move the mouse pointer to a place on the screen where it isn't viewable. Move the mouse pointer up or to the left to make it visible again.

### 6.4 Using an External Keyboard

An external keyboard provides the same data entry features as the on-screen virtual keypad presented by the touch screen user interface with faster typing.

The “Settings/Customize/On screen keyboard active” check-box allows you to choose whether the on-screen virtual keypad is presented or not when data entry is required (name, sample ID, target value etc.). If the check box is left unchecked, then the on-screen virtual keypads are not presented and an external keyboard is required for data entry. If the check box remains checked, the external keyboard can still be used by pressing the ‘ESC’ key on the external keyboard to dismiss the on-screen virtual keyboard and complete the data entry using the external keyboard.

6.5 Using the On-Screen Virtual Keypads

There are three types of virtual on-screen keypads:

- Alpha-numeric with shift and symbol options
- Numeric only
- Date

The particular virtual on-screen keypad presented is tailored to the data entry required. If only numeric values are accepted than a numeric-only keypad is presented. For selection of dates, a date keypad is presented.

The touch-pad of the ‘Abacus 5’ handles single taps only. The Shift and Symbol keys are ‘sticky’ on the virtual on-screen keypads like the ‘caps-lock’ on a PC key-board. Click/tap on the appropriate buttons to type in the required value. Complete the operation by clicking/tapping the ‘Done’ button. The ‘Symbols’ button presents special characters such as the question mark or ampersand.



Figure 15. Alphabetic Virtual On-Screen Keypad




Figure 16. ‘Sticky’ Shift and Symbol Buttons




Figure 17. Numeric Virtual On-screen Keypad




Figure 18. Date Virtual On-Screen Keypad

## 6.6 Using the Graphical User Interface

The graphical user interface (GUI) of the 'Abacus 5' is divided into three main sections. Refer to Figure 19 to see these sections:

1. Quick links to main functions (top section)
2. Interactive display area (middle section)
3. Status display (bottom section)

The main menu shown below is displayed in the interactive display area when the 'Abacus 5' analyzer is initially powered up.



**Figure 19. Graphical User Interface Sections**

### 6.6.1 Quick Links

The upper part of the GUI is always visible during all operations and presents several quick link buttons. Selecting one of the quick link buttons determines what is displayed in the interactive display area. The selected button is framed on three sides for clear visual feedback. The quick link for Autosampler functions is enabled only if the optional Autosampler is connected.

### 6.6.2 The Interactive Display Area

The middle section is used to display the main menu, functions, information, results, and all the panels that contain information about the active mode.

### 6.6.3 The Status Display

The bottom section is called the Status display, and presents the user with the following:

- The Help icon
- Status of the Autosampler unit
- The state of the reagent and waste containers
- Error/ warning icons
- Printer status
- Status of the pneumatic system (standby, manual measure, rinse, clean)
- Progress of the active operation
- Time
- Input language

Clicking or tapping on an error or warning icon navigates directly to the relevant menu where the problem can be addressed, or provides the user with information on how to resolve the problem.



Figure 20. The Status Display

#### 6.6.4 Entering Information

The 'Abacus 5' sometimes requires the operator to enter information such as sample ID or patient data. In order to enter information into a data field simply tap or click on it. The data field will change color to yellow indicating that data entry is active for that field. Information can be entered using the virtual on-screen keypad or with an optional external keyboard.

### 6.7 The Menu System

The graphical user interface of the 'Abacus 5' analyzer has primary functions that are accessible directly by means of icons or controls that are always available to the operator. The 'Abacus 5' analyzer also has secondary functions that are available via two or more taps or clicks through the menu structure, or menu tree.

#### 6.7.1 Primary Menu Items

The following menu items or functions are directly accessible:

- Start a manual single tube measurement
- Start an automated measurement
- Access the database
- Initiate printing
- Access the main menu
- Access the control panel of the Autosampler (opens with double click or tap)
- Access the time adjustment panel (opens with double click/tap)
- Access the warnings panel (only visible if an error or warning was reported)

#### 6.7.2 Starting a Manual Single Tube Measurement

Tapping or clicking on the 'Measure' icon with the single sample tube graphic allows access to the following functions:

- Initiate a primary blank measurement
- Make adjustments and settings about the manual measurement:
  - Enter the Sample ID (unique or auto-increment)
  - Select the associated patient
  - Select the sample type or patient mode
- Start a manual measurement
- View the measurements results in detail:
  - View the measured parameters
  - View the calculated parameters
  - View parameters that are out of the normal range
  - View error, warning, or clinical flags

See section 8.4.1 for more information about manual measurements.

### 6.7.3 Start Automated Measurements

Tapping or clicking the 'AS' icon with the multiple sample tube graphic allows access to the following functions:

- Initiate automated measurements in one of three modes:
  - Full scan mode
  - Free list mode
  - Selected samples mode
  - Remote worklist mode
- Define sample parameters (the parameter list depends which one of the three sample modes are selected)
- Monitor the progress of the automated measurements
- View the details completed measurements

See section 8.4.2 for more information about automated measurements.

### 6.7.4 Access the Database

Tapping or clicking the 'Database' icon with the file cabinet graphic allows access to the following functions:

- Select completed and measurement stored as database records
- Export database records for storage or import results for viewing
- Delete database records
- Send database records to the LIS
- Review the details of selected measurements

See section 1 for more information about accessing and using the 'Abacus 5' database.

### 6.7.5 Initiate Printing

Tapping or clicking the 'Print' icon with the printer graphic will print the active panel or measurement result according to the printer settings.

See section 15.6 for more information about printing on the 'Abacus 5' analyzer.

### 6.7.6 Main Menu

Tapping or clicking the 'Menu' icon with the main menu graphic displays and allows access to the main menu. The main menu is automatically displayed when the 'Abacus 5' analyzer is started up. The main menu icons are arranged in an intuitive and attractive circular design.




Figure 21. 'Abacus 5' Main Menu

### 6.7.7 Autosampler Control Panel

If you have an optional Autosampler installed, the Autosampler control panel can be opened by double clicking/tapping the Autosampler status on the left side of the status display on the bottom of the screen. The Autosampler control panel allows you to see if the Autosampler is operational and perform an Autosampler initialization.

See section 8.4.2.6 for more information about the Autosampler control panel.

### 6.7.8 Adjust the Time and Date

Double click/tap the time display on the right side of the status display on the bottom of the screen to bring up a panel for adjusting the stored time and date of the 'Abacus 5' analyzer. See section 6.7.8 for more information about time and date adjustment.

### 6.7.9 Open the Warning Panel

If an error or warning is currently active, an exclamation point icon appears on the status display. Click the icon to review the active error or warning and jump to the appropriate menu where the error or warning can be addressed (such as reagent replacement).

### 6.7.10 Menu Tree

The following is a listing of the 'Abacus 5' analyzer main menu tree:

<i>Patient</i>
New: Add a new patient
Edit: Modify an individual patient's information
Details: View an individual patient's information
See section 1 for more information about patient data
<i>Exit</i>
Cancel: Cancels shutdown of the 'Abacus 5' analyzer



Prepare for shipment: Performs an extensive drain procedure of the 'Abacus 5' pneumatic system in preparation for shipment or an extended period of inactivity
Log off: Logs off the current user session
Shutdown: Prepares the 'Abacus 5' for a complete shutdown
<i>See section 7.3 to learn how to power down the 'Abacus 5' analyzer correctly</i>
<b>Settings</b>
Customize
Skins: Change the display skin
Patient's displayed data: Patient identifier visible in the database display
Language: Change the active language
Limit Style: View parameter normal range status in numerical or graphical format
Sound Volume: Change the volume of the 'Abacus 5' built-in speaker
On screen keyboard active: Turn off on-screen virtual keypads for use with an optional external keyboard
Demo mode: The analyzer will not perform any pneumatic action
Laboratory: Seven lines of text of lab information be displayed on each printed page
External devices
Sending port baud rate: Select transmission rate for serial LIS
Automatic LIS: Select whether to automatically transmit results each time a sample is processed
LIS: Select Ethernet LIS connection
IP: IP address of Ethernet-based LIS host computer listening port
Port: Port number of Ethernet-based LIS host computer listening port
Bidirectional LIS: Selects whether Autosampler information downloads are allowed
System
Database columns order: Change visibility and position of the database columns
Special flags(G, L): Select whether G, L flags are displayed or printed
Waste container volume: Selects a 10L, 20L, or no waste container (direct to drain)
Database display limit: Selects all or only last month of database record display
Use only Sarstedt-Monovette tub from sample rotor: Change sampling adjustments to accommodate Sarstedt Monovette tubes
Standby time: Time period of inactivity before the 'Abacus 5' automatically enters standby
Offline rinsing frequency: Time period of inactivity before the 'Abacus 5' performs an offline rinse to keep the pneumatic system in peak operating condition
Screen saver: Time period of inactivity before screen saver takes effect
Blank Request Time: Time interval between system requests to process a blank specimen (0 h = no blank request reminders)
Units
HGB unit: Selects the units for HGB and HGB-derived or calculated parameters
Count unit: Selects the units for WBC, RBC, and PLT and derived or calculated parameters
PCT/HCT unit: Selects the units for PCT/HCT parameters
MCV unit: Selects the units for MCV parameters
Printer
Printer: Selects which printer to print to
Printer status: Display the currently selected printer status
Color printing: Selects whether to print in color
Double sided printing: Selects whether to print multi-page printouts on both sides of a page
Items in queue: Number of items in the printer queue
Cancel all jobs: Cancels all items in the selected printer's queue
Pre-printed paper's header length (px) – Reserves space (pixels) at top of page for pre-printed letterhead
Automatic print: Selects whether to print for every new sample processed
Logo visible: Selects whether to print or omit graphical logo to printout
Warning flags are visible: Selects whether to print or omit warning flags to printout
Diagnostic flags are visible: Selects whether to print or omit Diagnostic flags to printout
Refresh printers list: Refreshes the list of displayed printers

Profile limits: Enter normal ranges for Human, Male, Female, and Alternate 1 and 2 profile limits
X-B: Change X-B limits and targets
User: Add and manage users of the 'Abacus 5' analyzer (only available to administrator users)
Delete user: Initiates deletion of a user account
Administrator: Selects or unselects administrator privileges for the currently selected user
<i>See section 1 for more information about 'Abacus 5' analyzer user settings</i>
<b>Maintenance</b>
Cleaning: Perform extended cleaning of various pneumatic components
Draining: Drain various internal system reservoirs
Empty chamber: Empty RBC, WBC, and mix chambers
Prime: Prime one or all reagent reservoirs
Fill: Fills the fluidic components of the system with reagents
Touchscreen: Calibrate the touchscreen location for taps and button presses
<i>See section 1 for more information about 'Abacus 5' maintenance procedures</i>
<b>Calibration</b>
Calibrate
Calibration mode: Select type of calibration procedure to run
Calibration type: Select human blood or calibrator material
Target values: Enter the target ranges for each calibrated parameter
Cancel: Cancels a calibration procedure
Next: Initiates processing of calibration run
View calibrations: View and delete previous calibrations
<i>See section 1 for more information about 'Abacus 5' calibration</i>
<b>QC</b>
QC Measure: Initiate processing of a control sample
QC Reference select: Select which stored QC reference this measurement belongs to
Set QC Reference: Create a new stored QC reference
Load QC reference: load QC reference values from removable media such as a flash drive
Save reference: Save entered QC reference values
View QC references: Browse stored QC references
Delete: Delete a stored QC reference
Details: View detailed results for a stored QC reference
View QC data: Browse individual QC sample measurements
Delete: Delete currently selected QC run
View QC diagrams: View QC diagrams for currently selected
Details: View detailed result screen of a QC run
View QC diagrams: View Levy-Jennings diagrams of QC data for the currently selected QC reference
View QC data: Browse individual QC sample measurements
Next diagrams: Browse the next set of QC diagrams
View X-B data: Browse individual X-B sample measurements
Reject: Reject a run from being included in X-B calculations
View rejected: View previously rejected X-B runs
Undo reject: Restore a previously rejected X-B run and use in X-B calculations
View X-B diagrams: Navigate directly to X-B diagrams
View X-B diagrams: View X-B diagrams of X-B data
<i>See section 1 for more information about quality control on the 'Abacus 5' analyzer</i>
<b>Diagnostics</b>
Selftest: Perform electronic and/or pneumatic system self tests
Load last selftest: Loads the results of the last selftest
Start electronic: Starts the electronic tests
Start both: Starts both electronic and pneumatic tests
Log: Review and obtain additional details about events in the system log
Details: Display details about the selected event in the event log
Selftest: Navigates directly to self test screen

Reagent status: Reset levels for individual reagents and waste (or all) for reagent replacement
Statistics: Provides operating statistics such as cycle counts, errors, etc.
Information: Provides version information for all software items in the 'Abacus 5' analyzer

Table 7. 'Abacus 5' Analyzer Menu Tree

#### 6.7.11 Safety Access Codes

Some functions require password protected confirmation. For example, deleting a record from the database requires entry of a password.

The security password for deleting database entries is: **555**.

## 7 Start Up and Shut Down the 'Abacus 5'

### 7.1 Start Up and Shut Down Overview

The 'Abacus 5' is an advanced instrument with computer and pneumatic subsystems that need to be properly started up and shut down for continuing reliable performance. Not following proper startup and shutdown procedures may lead to increased maintenance and unnecessary service calls.

The 'Abacus 5' analyzer start up sequence requires four steps:

- Power up the electrical system.
- Start up the internal user interface computer and log in if the multi-user mode is activated.
- Start up the pneumatic system components.
- Perform and accept an initial 'blank' measurement.

A normal shut down of the 'Abacus 5' analyzer has the reverse steps:

- Clean and empty the pneumatic system components.
- Log off (multi-user mode is activated) and properly shut down the internal user interface computer.
- Power down the electrical subsystem.

The 'Abacus 5' analyzer also has a 'prepare for shipment' shut down option that performs an extended cleaning and emptying of the pneumatic components. This option can be used to prepare the 'Abacus 5' analyzer for an extended period of inactivity or for shipping to a new location.

There are also emergency shutdown procedures in case of emergency.

### 7.2 Starting Up the 'Abacus 5' Analyzer

#### 7.2.1 Visual Inspection

Perform a visual inspection of the 'Abacus 5' analyzer before starting up. This is particularly important if the 'Abacus 5' is operated by multiple operators.

- Check the reagent tubes by gently pulling on them to make sure that connections are not loose.
- Open the front cover and ensure that there is no leakage around the shear valve or the syringe pumps.
- Check for salt stains around the reagent tubing and the reagent containers.
- Check that the power cord is plugged in both into the back panel of the 'Abacus 5' analyzer and to the wall outlet.
- Empty the waste container. After start up reset the waste level in the "Main menu/ Diagnostics/Reagent status panel".

#### 7.2.2 Power Up the 'Abacus 5' Analyzer

To power up the 'Abacus 5' analyzer, perform the following steps:

- Turn on the main power switch (small switch) on the rear panel of the 'Abacus 5' analyzer located near the power connection to the 'up' position labeled '1' to power on the electrical system of the 'Abacus 5' analyzer.
- If the optional Autosampler is installed, turn the power switch on the right side of the Autosampler to the 'on' position labeled '1' to power on the electrical system of the Autosampler.
- Turn on printers or any peripherals that have their own power switches.

The electrical system of the 'Abacus 5' analyzer can remain powered up indefinitely. Please follow local regulations or laboratory procedures when deciding whether to fully power down the 'Abacus 5' analyzer after work hours.

### 7.2.3 Start Up the User Interface

To start up the internal user interface computer of the 'Abacus 5' analyzer, flip the standby switch near the top of the rear panel of the 'Abacus 5' analyzer to the 'up' position. The standby switch springs back to the 'down' position after activation.




Figure 22. 'Abacus 5' Analyzer Startup and Log On Screens

The following sequence of events takes place after enabling the standby switch:

- The internal computer starts up displaying the 'EPIA' logo in a rotated way.
- The embedded Windows® XP® operating system starts up. As the display initializes the screen changes color.
- The 'Abacus 5' application software starts up displaying the 'Abacus 5' logo screen. Individual component start up status is displayed in the bottom right corner of the logo.
- If the multi-user option is enabled then a user ID and password is required.
- The 'Abacus 5' main menu is displayed.

At this point, the 'Abacus 5' analyzer user interface is ready to use, but the analyzer is not yet ready to measure samples. This requires initialization of the 'Abacus 5' pneumatic system. The following functions are available at this time:

- Use the patient database to add and modify entries.
- Review, archive, send to LIS and print existing measurement records.
- Review calibration and QC results and history.
- Change settings.

### 7.2.4 User Logon

If the multi-user feature is activated, a user must be logged in before the user interface of the 'Abacus 5' can be accessed. You can log in on two different ways:

- Enter the user ID and password assigned to your account.
- Create a new user account.

The multi-user feature of the 'Abacus 5' analyzer tracks which user is logged on when important activities are performed such as sample or QC measurements or calibration. The multi-user feature allows full access to administrator users while restricting regular users from actions that change the way the 'Abacus 5' analyzer operates, such as changing calibration settings.

See section 15.7 for details about user management.

### 7.2.5 Pneumatic System Start and Blank Measurement

Before measuring any blood or control samples, the pneumatic system components of the 'Abacus 5' analyzer need to be initialized and a blank measurement run and accepted. To start the pneumatic system and run a blank measurement, click or tap the 'Measure' icon with the single test tube graphic on the top left corner of the screen. The following sequence takes place:

- The color of the 'Start' button changes to red.
- The reagents are primed.
- The measuring system is flushed.
- The sample-door, syringe, needle x/y movement, and shear valve motors are activated and tested.
- The vacuum system is activated, and pump sensors and vacuum leakage is tested.
- A blank measurement is initiated automatically.
- An empty result screen is displayed.
- The blank measurement is performed.
- The results of the blank measurement are displayed.

The results of a blank measurement must be accepted before blood or control samples can be processed. If the results from the first blank measurement are not acceptable, it does not necessarily mean that something is wrong about your 'Abacus 5' analyzer. Please consider the following:

- Blood proteins or other particles can accumulate in the pneumatic tubing or components.
- Although the reagents contain antimicrobial agents, some microbial growth may still be present.
- An extended period of inactivity without performing a 'prepare for shipment' procedure may result in salt crystal development in the pneumatic components.

Running additional blank measurements will usually bring the blank measurements down to an acceptable level.

Blank measurements will be flagged if they exceed the following thresholds:

- WBC  $>0.20 \times 10^3$  cells/ $\mu$ L
- RBC  $>0.05 \times 10^6$  cells/ $\mu$ L
- PLT  $>15 \times 10^3$  cells/ $\mu$ L
- HGB  $>1.0$  g/dL

If the 'Abacus 5' was not in use for an extended period of time, the start up of the pneumatic system components can take significantly longer and require more blank measurements. To prevent these problems, use the 'prepare for shipment' process when the 'Abacus 5' analyzer will be inactive for an extended time.

### 7.3 Exiting the 'Abacus 5' Analyzer

After you have completed lab operations and you are ready to exit the 'Abacus 5' analyzer, you can:

- Log off
- Shut down the 'Abacus 5'
- Power off the 'Abacus 5'
- Prepare the Abacus 5' for shipment




Figure 23. Abacus 5 Shutdown Options

Click or tap the Exit icon on the main menu to access various exit options. Click or tap the 'Cancel' button to return to normal operation.

#### 7.3.1 Log Off

Clicking or tapping the Exit icon on the main menu brings up the 'Shutting down Analyzer' panel. Select 'Log off' to end the active user session.




Figure 24. Log Off

After the current user is logged off the logon panel is displayed.

### 7.3.2 Shut Down

The shutdown procedure must be used to maintain the reliable operation of the instrument. To properly shut down the 'Abacus 5' analyzer, click or tap the Exit icon on the main menu and select 'Shutdown.'



**Figure 25. Shut Down**


The shutdown procedure performs two key functions:

- If the pneumatic system components of the 'Abacus 5' analyzer were initialized, a shut down procedure is performed that rinses the pneumatic components. This prepares the 'Abacus 5' analyzer for a period of inactivity of up to 5 days.
- The 'Abacus 5' operating software and the Windows® XP® Embedded operating system are shut down. The analyzer saves all settings, closes the file system and database, and automatically shuts down the internal user interface computer.

The 'Abacus 5' analyzer may ask to perform a cleaning cycle based on the analyzer's usage patterns and the state of the pneumatic components. If this occurs, the 'Abacus 5' analyzer will request that a tube with cleaning reagent be introduced at the sample rotor. The cleaning procedure may take a few minutes. Do not turn off the analyzer until it indicates that it is safe to do so.


### 7.3.3 Prepare for Shipment

If the 'Abacus 5' analyzer is not to be used for more than 5 days or is being transported to a new location then the 'prepare for shipment' procedure should be performed instead the regular shutdown procedure. Click or tap the 'Exit' icon on the main menu, then select the 'Prepare for shipment' button.



**Figure 26. Prepare for Shipment**



	<p>If the 'Abacus 5' is powered off without a proper shutdown or prepare for shipment procedure, the following may occur:</p> <ul style="list-style-type: none"> <li>• The pneumatic components may lock in an undefined state. The next startup and pneumatic initialization may take longer with increased reagent consumption.</li> <li>• The 'Abacus 5' application software and/or the Windows® XP® Embedded operating system can suffer damage or data loss.</li> <li>• The 'Abacus 5' analyzer may have to be restarted by a service engineer after an improper shut down or power off sequence.</li> </ul>
---	--

The 'Prepare for Shipment' procedure should be used for following reasons:

- The reagents may spill out if the 'Abacus 5' is tilted, potentially causing electrical short circuit or corrosion.
- Reagents can flow out of the 'Abacus 5' if the reagent intakes are disconnected.
- Evaporation may cause reagents left in the pneumatic components in the 'Abacus 5' analyzer to become concentrated and form crystals that increase blank results, damage pneumatic components, and may require a service engineer to restore proper operation.
- Microbial growth may occur if reagents are left in an inactive analyzer for extended time periods, causing high blank results and may require a service engineer to restore proper operation.
- Accidental storage at very low temperatures may cause freezing that damages the pneumatic components.

It is strongly recommended that the 'Prepare for shipment' procedure be used if the analyzer will be inactive for more than five days. Higher temperature or lower humidity decreases the period of time the 'Abacus 5' can be idle without performing the 'Prepare for shipment' procedure.

The analyzer will guide you through the preparing for shipment process. This process takes about 30 minutes and requires the following items:

- "Prepare for shipment" tube set included with your 'Abacus 5' analyzer
- Distilled water

After initiating the 'Prepare for shipment' procedure the system will prepare itself for the draining process. Do not power off the analyzer during these steps unless the system indicates that is safe to do so.

**Step 1:** The analyzer will ask you to disconnect the reagents from the reagent connectors on the back panel except for the waste connection. Leave the waste container tubing connected to the analyzer. The system will drain all inner reagent buffers throughout the pneumatic system. This step takes approximately 9 minutes. Do not turn the analyzer off during this step.

**Step 2:** The instrument will ask you to connect distilled water to the reagent inputs.




Figure 27. Connect Distilled Water Panel

Use the 'Prepare for shipment' tube set included with your 'Abacus 5' analyzer to perform this step. The system will prime and fully rinse the tubing and pneumatic system with distilled water. This step takes approximately 7 minutes. Do not turn the analyzer off during this step.

**Step 3:** The analyzer will ask you to disconnect the 'Prepare for shipment' tube set from the reagent connectors on the back panel except for the waste connection. Leave the waste container tubing connected to the analyzer.





Figure 28. Remove Reagent Tubing Panel

The analyzer will drain all liquids from the instrument. This step takes approximately 10 minutes. Do not turn the analyzer off during this step. When the process is over, the system will notify you about switching off the instrument. Click or tap the 'OK' button on the notification panel and wait until the Windows® XP® Embedded exits and power off the user interface computer.

	<p>Do not turn off power of the 'Abacus 5' at this stage. First click or tap the 'OK' button, then wait until the Windows® XP® Embedded operating system exits and shuts down the internal 'Abacus 5' user interface computer.</p> <p>After the shut down process is completed, you can turn off the main power switch (small switch) on the rear panel of the 'Abacus 5' analyzer near the power connection to the 'down' position labeled '0'.</p>
---	--





Figure 29. Shut Down Ready Panel

**Step 4:** Disconnect the waste container tubing from the waste reagent connector on the back panel of the 'Abacus 5' analyzer.

Small droplets of distilled water may remain in the tubing after the 'Prepare for shipment' process completes. These droplets will not cause any malfunction during transportation or the next start up procedure.

#### 7.3.4 Emergency Shut Down

 The procedures described in this section can damage the 'Abacus 5' analyzer or cause data loss or damage to critical system computer files. Use these procedures in case of emergency only!

#### 7.3.4.1 Immediate Shut-Down

The immediate shut-down can be initiated by pressing the standby button near the top of the back panel of the 'Abacus 5' analyzer. This is the same button which is used to start up the 'Abacus 5' user interface computer. This shuts down the user interface computer without properly shutting down the pneumatic system components.

Possible side-effects:

- The pneumatic components may be left in an undefined state. The next startup and pneumatic initialization may take longer with increased reagent consumption.
- The last settings changes and measurement results are not saved.

Immediate shut down should be avoided as a routine procedure. It should only be used if loss of electrical power is imminent such as thunderstorms or if loss of uninterruptible power supply (UPS) backup power is imminent.

#### 7.3.4.2 Immediate Switch Off

In case of immediate danger you can power off the 'Abacus 5' analyzer by turning off laboratory power, or by turning off the main power switch (small switch) on the rear panel of the 'Abacus 5' analyzer located near the power connection to the 'down' position labeled '0' or unplugging the power cord from the wall outlet.

#### 7.3.5 Repackaging the 'Abacus 5' Analyzer

Before repackaging your 'Abacus 5' analyzer for transportation or storage, always do the following:

- Perform a 'Prepare for shipment' procedure.
- Fully power down the analyzer.
- Disconnect all reagent tubes, power cords and peripheral connections.
- Always use the original packaging materials to repack the 'Abacus 5' analyzer.

Prepare the following packaging materials were included in your original shipment:

- Mini-pallet
- Upper and lower packaging foams
- Upper and lower carton box
- Plastic bag

Place the pallet on the ground. Place the lower part of the carton box on the pallet, and place the lower packaging foam on the carton bottom.

Put the 'Abacus 5' into the plastic bag, and place it into the bottom foam part. Place the top packaging foam on the analyzer, and put the top carton cover box over the package. Align the lower carton box flaps so that they fit beneath the upper box edges.

Use shipping straps to gently but firmly hold the entire package together.




Figure 30. Repackaging the 'Abacus 5' Analyzer

## 8 Sample Measurement

### 8.1 Sample Types Supported by the 'Abacus 5'


In this chapter you will find information about the samples usable with the 'Abacus 5' instrument.

#### 8.1.1 Supported Sample Tubes Types

The 'Abacus 5' supports the following K3-EDTA, 13x75mm sample-tubes:

- Sarstedt Monovette®
- Becton Dickinson (BD) Vacutainer®
- Terumo Venosafe®

Only use sample tubes with K3-EDTA anticoagulant.

	<p>Please note that the Sarstedt Monovette® sample tubes have a false bottom that is not the same as the bottom of the sample tube. The sampling depth on the 'Abacus 5' analyzer must be changed for a Sarstedt Monovette® sample tube. To change the sampling depth for manually presented Sarstedt Monovette® tubes at the sample rotor, check the 'Use only Sarstedt Monovette tube from sample rotor' check box in the Main menu/Settings/System panel.</p> <p>The Autosampler automatically recognizes Sarstedt Monovette® tubes and automatically changes the sampling depth.</p>
--	--

You may use sample tubes from other vendors if certain conditions are fulfilled. For manually presented sample tubes at the sample rotor measurement check if the new tube type:

- Is mechanically compatible (13x75 mm family)
- Can be pierced or used in open vial mode
- Can be handled by the tube centralizing mechanisms on the 'Abacus 5' analyzer

It is suggested to contact your sales representative before using the new sample tube type.

For automated sample handling with the optional Autosampler, contact your sales representative before using a new sample tube type. The Autosampler may not reliably recognize sample tube types other than the officially supported sample tubes.

Although these sample tube types are designed for multiple piercing, it is recommended to replace the cap after 3-4 piercing cycles. Small rubber bits parts can crumble into the blood and adversely affect the performance of the 'Abacus 5' analyzer.

#### 8.1.2 Sampling Depth

The 'Abacus 5' is equipped with a piercing needle capable of aspirating samples from closed sample tubes. To prevent needle blockage from rubber sample tube cap material and blood particulates, the piercing needle hole is located on the side of the needle rather than at the bottom. Also, the sampling needle is given a safety clearance of two millimeters from the bottom of the tube to prevent accidental breakage of sample tubes. This means that the sampling needle cannot aspirate the last 8 millimeters (0.31 inches) of sample at the bottom of a tube. This results in sample volume that cannot be aspirated of 0.5 ml for a Beckton Dickinson Vacutainer® tube and 0.3 ml for a Sarstedt Monovette® tube.




Figure 30. Sampling Needle

### 8.1.3 Open or Closed Sample Tubes

The 'Abacus 5' can process both open and closed sample tubes in manual mode using the sample rotor. Manually presented sample tubes must be mixed properly before being presented to the 'Abacus 5' analyzer.

Only closed sample tubes will be automatically processed by the optional Autosampler. Open sample tubes placed in the Autosampler will be skipped and not processed.

### 8.1.4 Sample Collection and Handling

Please consider the following when collecting and handling blood samples:

- The 'Abacus 5' analyzer aspirates 110  $\mu$ L of sample in both closed and open vial modes. Under the following conditions, the 'Abacus 5' may aspirate less than 110  $\mu$ L of sample, but will detect and flag any inadequate sample aspirations:
  - The sample tube does not vent properly and contains residual vacuum
  - The viscosity of the blood is very high
- The 'Abacus 5' is unable to aspirate the last 8 mm (0.31 inches) or 0.3 to 0.5 ml of sample volume.
- Make sure that the sample tube is filled to the sample level line indicated by the tube manufacturer. The anticoagulant not only prevents the coagulation of blood, but also dilutes it. If the volume of blood is too low, than this dilution can adversely affect measurement results.
- Ensure that the blood and anticoagulant are properly mixed by inverting it at least 8 times. Do not shake the sample as this could damage blood cells.
- A minimum of 30 minutes should elapse between taking the blood sample and running it on the 'Abacus 5' for analysis to ensure that the interaction of blood and anticoagulant has fully stabilized. Check the data sheet of the sample tube for the exact period of time.
- Use only K3-EDTA anticoagulated fresh whole blood as sample.
- Analyze blood samples within 7 hours of collection.
- Samples should be analyzed at room temperature.
- Do not deliver or store the samples at temperatures above normal body temperature.



All samples should be handled as potentially infectious material and treated as a biohazard.

## 8.2 Sample Types and Modes

The 'Abacus 5' analyzer is designed to process the following sample types:

- Whole human venous blood to determine the hematology parameters of the blood
- Quality control (QC) materials for quality control purposes
- Calibrator materials for calibration of the hematology parameters

QC samples are artificially modified human or human and animal blood samples with known parameter values. QC samples can be stored for a period of time as indicated by the manufacturer.

The 'Abacus 5' analyzer works in the following modes. The sample mode must be selected before a measurement starts:

- Blank mode: No sample need be presented to run a blank. Blank mode checks that the measuring system is operating correctly and that blank measurement results are acceptable.
- Human blood mode: Choose one of five human blood modes (Human, Male, Female, Alternate 1 or Alternate 2). Each mode has a unique set of normal ranges associated with it. NOTE: Alternate 1 and Alternate 2 are sample modes that are different than Human, Male and Female, and have user-definable normal ranges.
- Control (QC) mode: The hematological parameters of QC materials are known prior to analysis. Controls are used to check the long term stability of the 'Abacus 5' analyzer.

Calibration with a calibrator material is not one of the 'Abacus 5' sample modes, but a specific procedure that can be accessed from the Calibration icon on the main menu.

## 8.3 Sample Identification



Use caution and ensure accurate data entry when manually entering sample ID or patient ID information into the user interface of the 'Abacus 5' analyzer to avoid possible misidentification of sample results.

The following sample information items are used to uniquely identify samples:

- Sample Mode
- Measurement time: assigned automatically by the 'Abacus 5' analyzer
- Result ID: unique database identifier automatically assigned by the 'Abacus 5' analyzer
- Sample ID: assigned by the operator
  - The 'Abacus 5' doesn't require the sample ID to be unique
  - The sample ID can submitted by:
    - Manual typing
    - Bar code reading (manual or Autosampler)
    - Auto-increment
  - The sample ID is not applicable on blank measurements
- Patient:
  - A record from the patient database connected to a measurement
  - Not applicable for blank and QC measurements
  - Can be changed by setting the name or the database identifier of the patient displayed with the measurement results
  - The 'default' patient (ID =1) patient is used if no patient is selected

## 8.4 Running Samples


The 'Abacus 5' supports both manual and automatic sample presentation modes. Automatic sample presentation requires the optional Autosampler. The 'Abacus 5' analyzer software automatically detects and initializes the Autosampler.

Before performing any measurement the 'Abacus 5' should be:

- Powered up and started
- The pneumatic system initialized and a blank measurement performed and accepted

### 8.4.1 Manual Mode

To use Sarstedt Monovette® manually presented tubes at the sample rotor, check the 'Use only Sarstedt Monovette tube from sample rotor' check box in the Main menu/Settings/System panel.

	Always ensure the 'Use only Sarstedt Monovette tube from sample rotor' check box in the Main menu/Settings/System panel is checked if you are manually presenting Sarstedt Monovette® tubes in the sample rotor. Needle damage, tube damage and spilled blood may result if Sarstedt Monovette® tubes are presented and this option is not checked.
---	---

Select the 'Measure' quick link at the top left of the screen to initiate a manual measurement. The new measurement panel opens, giving you have the following options:


- Review the results of the previous measurement by using the 'Last Measure' function. This feature can be useful if you were interrupted during the measurement batch, and you want to double-check which sample is the last processed sample.
- Define the sample id of the sample to be processed by clicking or tapping on the "Sample ID" field:
  - You can type it in the Sample ID using the virtual on-screen keypad or an external keyboard.
  - You can use an external (USB) bar-code reader.
  - You can use an auto-generated sample ID.
- Activate or deactivate the automatic generation of sample IDs. If the 'Autoincrement mode' check box is checked, then the 'Abacus 5' will add 1 to the last sample ID if it was fully numeric or offer '0' if the last sample ID had any non-numeric characters.
- Connect a patient to the sample by clicking/ tapping on the Patient field:
  - The name of this field can be 'Name' or 'ID' depending on the 'Main Menu Settings/Customize/Patient's displayed data' setting.
  - If you don't select a patient then the default patient will be used.
  - As you select the patient field, the patient menu is displayed. You can select an existing patient or define a new one. See section 1 or more detail about Patients on the 'Abacus 5' analyzer.
- Select the sample mode:
  - Blank: a blank measurement will be processed. This is the only option available if the pneumatic system is not initialized, or the 'Abacus 5' was inactive for an extended period of time. After performing and accepting a blank measurement the other modes will be available.
  - Control: a QC sample is expected. QC measurements initiated from this menu will not be saved into the QC database, and are not part of the QC procedure. This mode is intended only as a quick system checks. See section 1 for more information about the QC process.

- Human blood modes: Human, Male, Female, Alternate 1, and Alternate 2. Select the mode which has the appropriate normal ranges. NOTE: Alternate 1 and Alternate 2 are sample modes that are different than Human, Male and Female, and have user-definable normal ranges.
- Start the measurement by clicking or tapping the 'Start' button on the screen or by pressing the physical 'Start' button located on the front panel of the 'Abacus 5'.

Figure 31. Sample Measurement Panel

After starting the measurement the sample door rotates and takes the sample tube inside the 'Abacus 5' analyzer. This is a key safety feature of the 'Abacus 5' and protects the operator from the sample needle during routine operation. The Start button changes color to red to indicate that the system is busy processing the sample. The sample is aspirated inside the analyzer and returned to the operator when the sampling process is complete. As soon as the measurement results are available the result screen is displayed. See section 1 for more information about the results screen.





**Figure 32. Sample Processing Result Screen**

The 'Abacus 5' analyzer will be ready to take the next sample after it has finished cleaning and preparing for the next sample. The Start button will change to green indicating it is ready to start the next sample measurement.

On the result screen you have the following options:


- You can print the results by clicking/ tapping on the 'Print' quick link at the top of the screen. Please note that the printing is a feature that is always available. See section 6.6.1 for more details about printing.
- Return to the measurement panel by clicking or tapping the using the 'Back' button.
- Run the next sample in routine mode by:
  - Editing the 'Next sample ID' just below the actual results
  - Accepting the auto-generated sample ID
  - Starting the next measurement by clicking/ tapping the 'Start' button on the screen or by pressing the physical 'Start' button on the front panel of the 'Abacus 5'.
  - Please note that in this mode the patient and sample mode will be carried over from the previous sample unless they are changed.

#### 8.4.2 Automatic Mode

Only human samples in closed sample tubes can be processed in automatic mode from the Autosampler. QC samples and open tubes must be processed in manual mode.

After starting up the 'Abacus 5', select the 'Measure' quick link. Process a new blank must be processed if this is the first time a new measurement panel is opened since starting up, or if an earlier blank has expired due to inactivity. An accepted blank measurement is required before starting automatic sample processing.

Open the AS panel by using the “AS” quick link at the top of the screen.



**Figure 33. Autosampler Info and AS Panel**

From the AS panel, you can initiate automated sample measurement in 3 modes:

- **Full Scan mode:**
  - The Autosampler will scan the whole sample tray
  - The same sample mode (Human, Male, Female, Alternate 1, Alternate 2) will be used for all the samples (NOTE: Alternate 1 and Alternate 2 are sample modes that are different than Human, Male and Female, and have user-definable normal ranges)
  - The sample ID is generated from the sample tube barcode (if available)
  - It is not possible to assign a patient to automatically processed samples. The default patient will be used for all samples.
- **Free List mode:**
  - The samples are defined in a list.
  - The following parameters can be individually defined per sample:
    - Sample ID
    - Patient
    - Sample mode
  - The samples should be placed in the same order in the sample tray as the samples defined in the list.
  - Empty sample positions will be skipped on the sample tray.
- **Selected Samples mode:**
  - A sample can be defined for each position on the sample tray.
  - Only defined samples positions on the sample tray will be processed. Additional samples will be skipped and missing samples will be marked for later review.

- The following parameters can be defined per sample:
  - Sample ID
  - Patient
  - Sample mode
- Remote Worklist (Free list)
  - The samples are defined in a list received from the LIS.
  - The samples should be placed in the same order in the sample tray as the samples defined in the list.
  - Empty sample positions will be skipped on the sample tray.

Samples are automatically scanned and processed in the following order:

- Rack order starting with rack 'A' to rack 'J'
  - Rack 'A' is the rack nearest to the front panel of the 'Abacus 5'
  - Rack 'J' is the rack nearest to the back panel of the 'Abacus 5'
- Sample order for each rack starting with sample '1' to sample '10'
  - Sample '1' is the sample closest to the 'Abacus 5'
  - Sample '10' is the sample closest to rounded edge of the Autosampler cover
- Sample A1 is the first sample processed and sample J10 is the last.


#### 8.4.2.1 Supported Bar Code Symbolologies

Some of the 'Abacus 5' analyzer automatic processing modes use the bar code on the sample tube to determine the sample ID of the sample tube. The following bar code symbolologies are supported by Autosampler bar code reader:

Bar Code Symbolology	Checksum Status
Code 128	Uses check digits
Code 39	No check digits used. If check digit are present they are considered part of the sample ID
Interleave 2 of 5	No check digits used. If check digit are present they are considered part of the sample ID
Codabar	Uses check digits
Code 93	Uses check digits
Code 11	Uses check digits

Table 8. Supported Autosampler Bar Code Reader Symbolologies

The use of bar code symbolologies with check digits will identify incorrect bar code reads in the unlikely event that the Autosampler bar code reader incorrectly reads the sample tube bar code. Diatron recommends the use of bar code symbolologies that uses check digits to help prevent sample misidentification.


	<p>Incorrect reading of a sample tube bar code without check digits may result in sample misidentification. Use bar code symbolologies with check digits to prevent this.</p>
---	---

#### 8.4.2.2 Full Scan Mode

Full Scan mode is suitable for laboratories where:

- There are batches of samples with the same sample mode,
- Barcode labels are in use, and
- An LIS system rather than the 'Abacus 5' is used to associate patients with sample results.


In Full Scan mode it is possible to add more samples to the batch after starting the automated measurement. As usual, an accepted blank is necessary before starting automatic sample processing.



**Figure 34. Full Scan Automatic Processing Mode**

On the Full Scan mode panel you can:

- Return to the AS panel to select a different type of automated measurement by using the 'Back' function.
- Select a common measure mode for all samples in the tray.
- Swap between the tray and list views by using the 'View list' and 'View tray' functions.
- Start the automated measurement by clicking or tapping on the 'Start' button on the GUI.



**Figure 35. Full Scan Mode Tray View Progress**

After starting the automated measurement in Full Scan mode, the Autosampler will scan all the sample positions on the tray. Any rows with missing racks will be skipped. The Autosampler can detect the tube type and whether the tube has a cap. Sample tubes will be skipped if it has no cap or the Autosampler cannot determine the tube type.

The Autosampler then mixes the sample tube and reads the barcode label if one is used. The sample tube is then moved to the aspiration position where the sample is aspirated and measurement is performed. The aspiration and measurement process is the same in automatic mode as in manual mode.

The sample tray progress can be followed in the tray view or the list view. The results for completed samples can be reviewed by clicking or tapping on a completed sample.

The patient ID is not editable in the Full Scan mode as it is automatically set to the default patient. The sample ID is not editable as it is read from the barcode label. If the barcode label is missing or not readable, the 'Abacus 5' analyzer assigns a default value and allows manual entry of the sample ID.



**Figure 36. Full Scan Mode Tray and List Views**

The Autosampler allows the addition of sample tubes to a Full Scan process already in progress. To add more sample tubes, perform the following steps:

- Click or tap the 'Stop' button. The 'Abacus 5' will prompt you to wait until the current sample is completed.
- Acknowledge the message by clicking or tapping the 'OK' button.
- Insert the new samples into a rack.
- Open the cover of the Autosampler.
- Place the new samples behind the last measured sample.
- Close the Autosampler cover.
- Start the automated measurement by clicking or tapping on the 'Start' button on the GUI or by pressing the physical Start button.

The 'Abacus 5' analyzer keeps track of the last analyzed sample on the sample tray of the Autosampler. After restarting the full scan process the 'Abacus 5' moves the Autosampler to the last processed position and continues scanning from there.

### 8.4.2.3 Free List Mode

Free List mode is suitable for laboratories where:

- The barcode labels are not used, or the bar codes are not the primary identifiers.
- Samples are arriving from multiple sources in small batches.
- In a given batch there are samples which should be measured using different modes for different normal ranges.

To start sample measurement in Free List mode, first perform and accept a blank measurement if one is needed. Click or tap the 'AS' quick link at the top of the screen and select 'Free list mode'. The Free List mode only requires that samples be placed in the same sequence as in the free list. It does not require that samples be matched to exact

rack and sample positions. It is not possible to add additional samples to the batch in free list mode after automated measurement has started.




Figure 37. Free List Mode Selection

On the Free List mode panel you can:

- Return to the AS panel to select a different type of automated measurement by using the 'Back' function
- Add or remove list items
- Define the following parameters for the new list items:
  - Measure mode
  - Sample ID
  - Patient ID
- Swap between the tray and list views by using the 'View list' and 'View tray' functions
- Start the automated measurement by clicking or tapping on the 'Start' button on the GUI or by pressing the physical Start button



Figure 38 Figure 39. Preparing a Free List

As with the Full Scan mode, the Autosampler will scan all the sample positions on the tray and rows with missing racks will be skipped. The Autosampler will only process tube types it recognizes with caps attached. The automated mixing, aspiration, and analysis process is the same as the Full Scan Mode, and the progress can be followed in tray or list view.

The main difference between the Free List mode and the Full Scan mode is that the Free List mode identifies samples by matching them sequentially to the list, and comparing their measurement results to the normal range associated with the sample mode selected for that sample.





Figure 40. Free List Mode Progress: List View





**Figure 41. Free List Mode Progress: Tray View**

Sample barcode labels are read and compared with the sample ID that was previously entered into the free list before starting automatic processing. If the two do not match, the barcode part of the tray view sample icon changes to red.

If the automated processing is interrupted for any reason, it can be restarted by clicking or tapping on the 'Start' button on the screen or by pressing the physical Start button.

#### **8.4.2.4 Selected Samples Mode**

Selected Samples mode is suitable for laboratories where:

- Barcode labels are used.
- In a given batch there are samples which should be measured using different modes for different normal ranges.
- Data security is a concern.

Selected Samples mode can also be used to rerun only certain samples in a sample tray already processed in one of the other two modes. Samples cannot be added to the sample tray once automatic processing is started in Selected Samples mode. An accepted blank measurement is required as with all sample processing.

In Selected Samples mode, the 'Abacus 5' analyzer will only look for samples in the sample positions selected for analysis. All other samples in sample positions that are not selected will be ignored. If a sample position is selected for processing but contains no sample tube, the sample position will be marked for review.

On the Selected Sample mode panel you can:

- Return to the AS panel to select a different type of automated measurement by using the 'Back' function.
- Define the samples to be measured:
  - Click or tap the sample tray position icon where the selected sample is located.
  - Select the appropriate measure mode.
  - Define the sample ID.
  - Define the patient ID.
  - Save the changes by clicking or tapping the 'Add' button.
- Remove selected samples by clicking or tapping the "Remove" button.
- Swap between the tray and list views by using the 'View list' and 'View tray' functions.

- Start the automated measurement by clicking or tapping on the 'Start' button on the GUI or by pressing the physical Start button



Figure 42. Selected Sample Mode Panel

#### 8.4.2.5 Interrupting Automatic Processing for Stat Samples

If one or more stat samples need to be processed, you can interrupt automatic processing to analyze the stat sample(s) without having to wait until the entire sample tray has been processed. To interrupt automatic processing for stat specimens, perform the following steps:

- Click or tap the 'Stop' button. The 'Abacus 5' will prompt you to wait until the current sample is completed.
- Acknowledge the message by clicking or tapping the 'OK' button.
- Place the stat sample into the sample rotor.
- As soon as the current measurement completes change to the manual measurement panel by clicking or tapping the 'Measure' quick-link at the upper left corner of the screen.
- Process the stat sample(s) in manual mode.
- Return to the automated measurement panel by clicking or tapping the 'AS' quick link;
- Restart the measurement by clicking or tapping the 'Start' button or by pressing the physical Start button.
- The 'Abacus 5' will continue automated measurement from where it was interrupted.

#### 8.4.2.6 Controlling the Autosampler

Detailed information about the Autosampler is available from the Autosampler info panel. You can also reset the Autosampler from this panel. Double click or tap the Autosampler status in the bottom left corner of the screen to view the Autosampler info panel. With this panel, you can command the Autosampler to repeat the last command, or reset the Autosampler. You can also view Autosampler messages with detailed information about the last action or sample tube processed.



Figure 43. Control The Autosampler With Info Panel

The following is a list of messages from the Autosampler:

- Home message:
  - Response after a successful reset
  - +HOME\_XX.YY (The XX.YY is the software revision number)
- VT zzzzzzz XY:
  - Vacutainer® style tube detected
  - Barcode content: “zzzzzzzzzz”
  - Tray position XY
- MV zzzzzzz XY:
  - Monovette® style tube detected
  - Barcode content: “zzzzzzzzzz”
  - Tray position XY
- +F : confirmation that the tube was moved from the mixing to the sampling position
- The power-on message contains the serial number of the Autosampler unit.

During the reset cycle the Autosampler moves the Autosampler components to their ‘home position’. The mixing stations moves to an upright position, the active rack moves back on to the sample tray, and the entire tray moves to its home position.

## 8.5 Result Display

After the analysis is completed, the result screen displays all measured parameters, flags, histograms and scatter plots. This information is automatically stored in the database, and stored result data can be retrieved at any time from the database.

Histograms and scatter plots can be zoomed for closer inspection. Click or tap the histogram or scatter plot to magnify it. Click or tap the ‘Close’ button on the upper right corner of the image to close it and return to the result display.




Figure 44. Results Display and Magnified Scatter Plot

Comparison of parameters to the normal range can be displayed numerically or graphically. See section 15.1 for more information on display of normal range status. Other kinds of flags are also displayed in the result screen.


See section 1 for detailed information about the 'Abacus 5' result screen.

## 8.6 Printing Reports

Click or tap the Print icon at the top of the screen to send a report (or any displayed screen) to the external printer.

## 8.7 The Measurement Process

The 'Abacus 5' analyzer measures 24 parameters of whole human blood samples. Internal dilutions are made to allow measurement of particles in the sample. Manual or automatic mode samples are processed in the same way.



**Figure 45. 'Abacus 5' Measurement Process**

When you present the sample tube and press the Start button, the 'Abacus 5' rotates the sample rotor and takes the sample tube inside the analyzer. The piercing needle pierces the sample tube cap (if present) and aspirates approximately 110 µL of primary blood sample. Blood sensors ensure that the primary blood sample is continuous and free of major bubbles as it is moved to the shear valve. The shear valve separates the primary sample into further samples of exact volumes of blood that are distributed to provide dilutions for the RBC, WBC, HGB and 5 Diff population measurements. Lysing reagents are used to eliminate RBCs for dilutions that measure WBCs.

RBC/PLT and WBC are measured with impedance technology through 70 and 80 µm apertures respectively. Two portions of the sample are used for determining a five subpopulation differential of WBC in two independent 4DIFF and BASO measurements.

Sample results are displayed on the screen as the process is nearing completion. Results are stored in the internal database and transmitted to an external LIS if configured to do so. The 'Abacus 5' performs a cleaning rinse of the internal pneumatic components and prepares to measure the next sample.

## 9 Result Interpretation

This chapter provides description about interpreting and understanding reports provided by the 'Abacus 5' analyzer.

### 9.1 The Result Screen




Figure 46. 'Abacus 5' Result Screen

The 'Abacus 5' analyzer result screen is divided into 4 major areas:

1. Sample identification information: contains identifying information used to uniquely identify this sample's information.
2. Parameter information: contains parameter names, values, units, and normal range. Normal range, linearity, and high blank flags are also displayed in this area.
3. Histograms and scatter diagrams: Visual elements that provide information about the measured population dispersion. These are active elements that can be individually magnified by clicking or tapping.
4. Flags area: this area displays warning, interpretive, and morphological flags associated with this sample. If flag display overflows its assigned space, an ellipsis will be displayed and the user can click on the information icon to see a detailed display of flags.

### 9.2 Sample Identification Information

This part of the results screen contains information that uniquely identifies a particular sample run.

Sample ID	The Sample ID is the ID entered by the operator from a keyboard or hand-held barcode reader for a manually processed specimen. For automatically processed specimens, the sample ID can be read from a sample tube barcode or entered into a list for one of the automatic processing list modes. Sample IDs can be reused in the 'Abacus 5' analyzer, but other identifying information together with the sample ID uniquely identify the sample run.
Patient ID	This is the patient associated with this sample run. If none was entered or selected, the default patient ID is set to 1.
Run Date/Time	This is the time and date that this sample was run.
Sample Mode	This is the sample mode selected for this specimen. Each sample mode applies a different set of normal ranges to this sample run.

Table 9. Sample Identification Information

### 9.3 Parameter Information

The parameter information area of the 'Abacus 5' analyzer result screen lists information on 24 parameters.

	Name	Value	Flags	Units	Normal Ranges
White Cell Parameters	WBC:	6.85		10 <sup>3</sup> /μL	[ 5.00 - 10.00 ]
	NEU:	3.80		10 <sup>3</sup> /μL	[ 2.00 - 7.50 ]
	LYM:	2.16		10 <sup>3</sup> /μL	[ 1.30 - 4.00 ]
	MON:	0.60		10 <sup>3</sup> /μL	[ 0.15 - 0.70 ]
	EO:	0.27		10 <sup>3</sup> /μL	[ 0.00 - 0.50 ]
	BAS:	0.02		10 <sup>3</sup> /μL	[ 0.00 - 0.15 ]
	NEU%:	55.5		%	[ 40.0 - 75.0 ]
	LYM%:	31.5		%	[ 21.0 - 40.0 ]
	MON%:	8.8	+	%	[ 3.0 - 7.0 ]
	EO%:	3.9		%	[ 0.0 - 5.0 ]
	BAS%:	0.3		%	[ 0.0 - 1.5 ]
Red Cell Parameters	RBC:	4.19		10 <sup>3</sup> /μL	[ 4.00 - 5.50 ]
	HGB:	10.7	-	g/dL	[ 12.0 - 17.4 ]
	HCT:	36.7		%	[ 36.0 - 52.0 ]
	MCV:	87.5		fL	[ 76.0 - 96.0 ]
	MCH:	25.5	-	pg	[ 27.0 - 32.0 ]
	MCHC:	29.2	-	g/dL	[ 30.0 - 35.0 ]
	RDWsd:	49.4	+	fL	[ 20.0 - 42.0 ]
	RDWcv:	11.8		%	[ 0.0 - 16.0 ]
Platelet Parameters	PLT:	378		10 <sup>3</sup> /μL	[ 150 - 400 ]
	PCT:	0.49		%	
	MPV:	12.9		fL	[ 8.0 - 15.0 ]
	PDWsd:	25.2		fL	
	PDWcv:	34.8		%	

Figure 47. Parameter Information Display

Each row shows the parameter name, value, flags (normal range, linearity range, and high blank), units, and normal range. The units area displays the differential percentage for the five leukocyte subpopulations (NEU, LYM, MON, EO, BAS).

Parameter names are listed with their standard abbreviation. Parameter values are always displayed according to their selected unit. In case of error during sample evaluation, the software will display 'E' or " ---" signs to indicate the

problems. Parameter values falling within the normal range are un-flagged and displayed with black text. Parameter values falling above the normal range are marked with a '+' and displayed with red text, and values falling below the normal range are marked with a '-', and displayed in blue text.

Normal range information can also be conveyed graphically as well as numerically. This setting can be changed in the Main menu/Settings/Customize panel.



Figure 48. Graphical Normal Range Display

The flags displayed in the parameter information area are as follows:

Normal range flags			
Flag	Meaning	Hierarchy	Color Code
-	Parameter is under normal range	Lowest	Blue
--	Parameter is under the half of the low limit of normal range	Low	Blue
+	Parameter is over the normal range	Lowest	Red
++	Parameter is over twice the high limit of normal range	Low	Red

Table 10. Normal Range Flags

Linearity Range Flag			
Flag	Meaning	Hierarchy	Color Code
*	The related parameter is out of the linearity range	High	N/A
**	The related parameter is out of the display range	High	N/A

Table 11. Linearity Range Flag

High blank flag			
Flag	Meaning	Hierarchy	Color Code
!	The blank value of the related primary parameter is high	Middle	N/A

Table 12. High Blank Flag

### 9.3.1 Scatter Diagrams and Histograms

The 'Abacus 5' analyzer displays the results of the optical measurements in scatter diagram representation. Scatter diagrams represent two-dimensional data. There are two scatter diagrams in the patient report: the 4-DIFF and BASO scatter diagrams.

The 4-DIFF scatter diagram displays cells identified after the first lysing and optical measurement process. Due to the measurement technology, cells are classified based on their optically detected properties: low and high angle scattered light intensity. The optical detector can measure the intensity of the light scattered or diffracted by each



cell. One portion of this scattered light is proportional to size of the cell; the other is proportional to the complexity of the internal structures in the cell.

Cells in a leukocyte sub-population have similar light scattering properties, allowing them to be grouped together and identified separately from other cell types. Different colors are used to identify various populations of blood cells.




Figure 49. Result Screen Scatter Diagrams

RBC and PLT impedance-based measurements are represented by histograms. These diagrams show the number of cells on the Y axis and the cell size on the X axis. Smaller cells are displayed on the left side of the histogram, and larger cells are displayed on the right side. The height of the histogram at a given size represents the number of cells with that size. Greater histogram height indicates more cells.




Figure 50. Result Screen Histograms

Both the RBC and PLT histograms can be zoomed. Discriminators (thresholds) are displayed with red color. Platelets are located on the leftmost part of the RBC histogram to the left of the red discriminator. The PLT histogram is a magnified version of this region.

### 9.3.2 Warnings

The warnings area of the result screen shows flags and messages related to the evaluation of the sample. Clicking or tapping on the individual sections of the results screen displays a detailed warnings panel. Click or tap OK to close the detailed warnings panel.




Figure 51. Warnings Section of Results Screen

The 'Abacus 5' analyzer provides Warning flags, Morphological flags, and Interpretive messages in the Warnings panel. The following tables describe these flags and messages in detail. It also provides guidance on actions the operator can take to correct the problem described by the flag.



Warning flags		Warning flags are raised when the analyzer detects a problem in the analysis process that may cause incorrect or suspect results or the parameters in the valid blank measurement was too high.	
Flag	Meaning	Cause	Action
A	4 Diff Alarm	The 4DIFF populations are not clearly distinguishable on the scatter diagram.	If other parameters predict a 'normal' sample then perform a cleaning procedure and repeat the measurement. If 4DIFF results reported then treat them as 'reduced reliability'. If not all the parameters reported or the sample is predicted as pathological then perform a manual count on a stained smear.
B	WBC blank is high	WBC > WBC blank high limit	Check cleanliness of the reagents and the 'Abacus 5'. Perform cleaning of the 'Abacus 5'. Repeat blank measurement. See chapter: 7.2.5:Pneumatics start and blank measurement Treat WBC result as 'reduced reliability'.
b	RBC blank is high	RBC > RBC blank high limit	Check cleanliness of the reagents and the 'Abacus 5'. Perform cleaning of the 'Abacus 5'. Repeat blank measurement See chapter: 7.2.5 :Pneumatics start and blank measurement Treat RBC result as 'reduced reliability'.
C	WBC clogging	Drift in the probe voltage of the WBC capillary. The two main reasons for this drift are the misconnected reagent and the clogging of the capillary.	Check the reagent connections. Perform a cleaning procedure! See chapter: maintenance. Please repeat the measurement.
c	RBC/PLT clogging	Drift in the probe voltage of the RBC capillary. The two main reasons for this drift are the misconnected reagent and the clogging of the capillary.	Check the reagent connections. Perform a cleaning procedure! See chapter: maintenance. Please repeat the measurement.
D	MON - NEU alert	The MON and NEU populations are not clearly distinguishable on the scatter-gram.	If other parameters predict a 'normal' sample then perform a cleaning procedure and repeat the measurement. If MON, MON%, NEU, NEU% results reported then treat them as 'reduced reliability'. If not all the parameters reported or the sample is predicted as pathological then perform a manual count on a stained smear.
E	EO - NEU alert	The EO and NEU populations are not clearly distinguishable on the scatter-gram.	If other parameters predict a 'normal' sample then perform a cleaning procedure and repeat the measurement. If EO, EO%, NEU, NEU% results reported then treat them as 'reduced reliability'. If not all the parameters reported or the sample is predicted as pathological then perform a manual count on a stained smear.
F	Differential blank high	More than 100 cells detected during the '4Diff-blank' procedure.	Check cleanliness of the reagents and the 'Abacus 5'. Perform flow cell cleaning of the 'Abacus 5'. Repeat blank measurement. See chapter: 6.1.5 :Pneumatics start and blank measurement Treat 4 Diff results as 'reduced reliability'.

Warning flags		Warning flags are raised when the analyzer detects a problem in the analysis process that may cause incorrect or suspect results or the parameters in the valid blank measurement was too high.	
Flag	Meaning	Cause	Action
f	BASO blank high	More than 100 cells detected during the 'Baso-blank' procedure.	Check cleanliness of the reagents and the 'Abacus 5'. Perform flow cell cleaning of the 'Abacus 5'. Repeat blank measurement. See chapter: 6.1.5 :Pneumatics start and blank measurement Treat BASO result as 'reduced reliability'.
H	HGB blank high	Last accepted blank result: HGB $\geq 1.0$ g/dl	Check cleanliness of the reagents and the 'Abacus 5'. Perform cleaning of the 'Abacus 5'. Repeat blank measurement. See chapter: 7.2.5:Pneumatics start and blank measurement Treat HGB result as 'reduced reliability'.
I	PLT URI	The PLT-RBC gap is not clearly detectable on the PLT-RBC histogram: small MCV, fractured RBC-s, aggregated PLT-s (cold blood), side-effect of blood-transfusion.	Check the sample quality. Repeat the sample. If the problem persists then perform a manual count on a stained smear.
M	Out of RBC linearity range	The RBC coincidence is over the limit: too high RBC count.	Check the sample homogeneity. Repeat the sample with manual pre-dilution.
m	Close to RBC linearity range	The RBC coincidence is close to the limit: too high RBC count.	Treat RBC result as 'reduced reliability'. Check the sample homogeneity. Repeat the sample with manual pre-dilution.
O	High linearity range limit exceeded	If any of RBC, PLT or WBC results are higher than high linearity range limits.	Perform manual pre-dilution of sample and re-measure it.
o	Low linearity range limit exceeded	If any of RBC, PLT or WBC results are lower than Low linearity range limits.	None
p	PLT blank high	Last accepted blank result: PLT $\geq 15 \times 10^3/\mu\text{l}$	Check cleanliness of the reagents and the 'Abacus 5'. Perform cleaning of the 'Abacus 5'. Repeat blank measurement. See chapter: 7.2.5:Pneumatics start and blank measurement Treat PLT result as 'reduced reliability'.
Q	MON - LYM alert	The MON and LYM populations are not clearly distinguishable on the scatter-gram.	If other parameters predict a 'normal' sample then perform a cleaning procedure and repeat the measurement. If MON, MON%, LYM, LYM% results reported then treat them as 'reduced reliability'. If not all the parameters reported or the sample is predicted as pathological then perform a manual count on a stained smear.
S	WBC measurement statistics warning	The distribution of the WBC detection changes in time. It points to clogging, non-homogeneous sample, cold sample, partial coagulation of the sample.	Perform cleaning of the 'Abacus 5'. Check the homogeneity, the temperature and the coagulation of the sample. Repeat the sample.
s	RBC measurement statistics warning	The distribution of the RBC detection changes in time. It points to clogging, non-homogeneous sample, cold sample, partial coagulation of the sample.	Perform cleaning of the 'Abacus 5'. Check the homogeneity, the temperature and the coagulation of the sample. Repeat the sample.

Warning flags		Warning flags are raised when the analyzer detects a problem in the analysis process that may cause incorrect or suspect results or the parameters in the valid blank measurement was too high.	
Flag	Meaning	Cause	Action
T	Blood detector error	Blood Detector malfunction (received signals are all zeros or not received signals)	Check BD and its connection to the LSDACQ card. Check sampling tube connection to BD.
u	BASO high	The number of the BASO population is high or there are lyse-resistant cells in the sample.	Repeat the measurement or perform a manual count on a stained smear.
V	WBC vacuum warning	There is a (partial) clogging or leakage in the WBC part of measurement system. A faulty/ worn-out pump can create vacuum errors as well.	Perform cleaning of the 'Abacus 5'. Perform a self-test. If vacuum problems persist then call for service.
v	RBC vacuum warning	There is a (partial) clogging or leakage in the RBC part of measurement system. A faulty/ worn-out pump can create vacuum errors as well.	Perform cleaning of the 'Abacus 5'. Perform a self-test. If vacuum problems persist then call for service.
W	Sampling warning	Blood Detector detected a sample volume that was too small.	Check sample volume in sample vial. Check sampling process. Re-run sample.
X	Differential error	Not enough cells detected during the 4DIFF WBC differentiating.	Perform a flow cell cleaning procedure! Repeat the measurement. If the problem persists then ask for service.
x	BASO error	Not enough cells detected during the BASO WBC differentiating.	If the total-WBC count is below $2 \times 10^3/\mu\text{l}$ then perform a manual count on a stained smear. Perform a flow cell cleaning procedure! Repeat the measurement. If the problem persists then ask for service.
Y	Differential percentage error	Algorithm error while distinguishing the WBC-s in the 4DIFF procedure.	Repeat the measurement. If the problem persists on a given sample then perform a manual count on a stained smear.
y	BASO percentage error	Algorithm error while distinguishing the WBC-s in the BASO procedure.	Repeat the measurement. If the problem persists on a given sample then perform a manual count on a stained smear.
Z	WBC noise high	During WBC count there is too much (more than 10% of the total) detected pulse on the first 10 channel of the 255 channel full range. This may be caused by lyse resistant RBC, electronic noise, contaminated WBC chamber or contaminated reagents.	Perform Clean process. Perform Hard Clean process.

Table 13. Warning Flags

Interpretive flag	Interpretive flags infer the presence of distributional abnormalities such as leukocytosis.
Flag	Meaning
Leukopenia	WBC --
Leukocytosis	WBC + or WBC ++
Neutropenia	NEU --
Neutrocytosis	NEU ++
Lymphopenia	LYM --

Interpretive flag	Interpretive flags infer the presence of distributional abnormalities such as leukocytosis.
Flag	Meaning
Lymphocytosis	LYM ++
Monocytosis	MON ++
Eosinophilia	EOS ++
Basophilia	BASO ++
Anemia	RBC - or RBC --
Polycythemia	RBC + or RBC ++
Microcytic RBC	MCV - or MCV --
Macrocytic RBC	MCV + or MCV ++
Hypochromic	MCHC - or MCHC --
Hyperchromic	MCHC + or MCHC ++
Anisocytosis	RDW + or RDW ++
Thrombocytopenia	PLT --
Thrombocytosis	PLT + or PLT ++
Microcytic PLT	MPV - or MPV --
Macrocytic PLT	MPV + or MPV ++

Table 14. Interpretive Flags

Morphological flag		Morphological flags are raised when the presence of morphologically abnormal cell types such as immature granulocytes are suspected.	
Flag	Meaning	Cause	Action
G	Immature granulocytes	The count in the region of (low - high angle scatter of Abacus 5) where the Immature granulocyte typically located is > 3% of the total WBC count.	Perform a manual count on a stained smear.
L	Atypical lymphocytes	The count in the region of (low - high angle scatter of Abacus 5) where the atypical lymphocytes typically located is > 1% of the total WBC count.	Perform a manual count on a stained smear

Table 15. Morphological Flags

## 10 Database Functions

The 'Abacus 5' database stores all measured data including sample results, QC, and patient data. Each sample result stored includes the complete parameter list, histograms, flags, and identifying data. The 'Abacus 5' has a data capacity of up to 100,000 total records.

To activate Database panel, clicking or tapping the Database quick link on the top of the screen will take you to the database view from any screen.




Figure 52. Accessing the Database

The database can be accessed at any time, even when processing samples in the background. You can return to the measurement panel anytime to view ongoing sample measurements.

### 10.1 Database Overview

The 'Abacus 5' analyzer database provides a powerful, easy-to-use method of storing, accessing, and managing patient and sample result information. The database of the 'Abacus 5' analyzer is capable of storing up to 100,000 measurement, control and patient records.

The database presents information in a table view of rows and columns. Each row in the database panel represents one sample measurement. Rows are also called 'records.' Each column in the database consists of sample identification information or measurement parameter results.

The database information can be scrolled using navigation buttons along the right and the bottom of the database view, or by moving a vertical scroll bar. Database records can be selected individually or in groups. Detailed result views are available for selected database records. Selected records can be printed, sent to an LIS and exported into a tab file. Various sort options are provided for quickly finding results.

Database records can be exported for storage and backup. Exported database information can be later imported for review. Imported records do not become part of your 'Abacus 5' analyzer's database. They are only imported for viewing, but can be sorted, selected, and managed the same way as with regular data stored inside the 'Abacus 5' analyzer database.



The ‘Abacus 5’ database displays two numerical counters separated by a slash under the left side of the database table. The left number is the number of selected rows, and the right number is the total number of rows in the database.

### 10.2 Scrolling the Database View

The limited size of the interactive database display acts like a viewport into the large amount of information inside your database. You can scroll the database view to bring desired information into view. The ‘Abacus 5’ analyzer provides several powerful and easy-to-use scrolling methods. The numbered items below match the items in Figure 53 below.

- 1. The scrollbar: tapping, holding and moving the scrollbar labeled ‘1’ with the tip of your finger on the touch screen display will quickly scroll information in the vertical direction. Clicking and dragging with an optional external mouse will achieve the same effect.
- 2. Tapping or clicking the buttons labeled ‘2’ will advance the display by one line or one column.
- 3. Tapping or clicking the buttons labeled ‘3’ will advance the display by one page of information at one time.
- 4. Tapping or clicking the buttons labeled ‘4’ will advance the display to the first or last item in the list.




Figure 53. Scrolling and Selecting

### 10.3 Sorting Database Information

The database is sorted by Result ID by default. Clicking or tapping the following database column headings changes the sort order of the database records:

Sample ID	Sample ID entered by the operator or read from a sample tube barcode label
Patient	Identifier that associates a result record with a patient records
Analysis Time	Time and date this sample tube was analyzed
Result ID	Unique sequential identifier assigned by the ‘Abacus 5’ to each records

Table 16. Database Sort Criteria

Clicking or tapping the same column heading again changes the ascending or descending order of values in the column.

## 10.4 Manual Selection of Database Records

Clicking or tapping on a row in the database display single-selects a row. Clicking or tapping it again de-selects the row. Selecting another row will de-select previously selected rows in single selection. Selected rows are highlighted with white text on a black background.

Checking the 'Multiselect' check box below the left side of the database table allows the selection of multiple rows. Tapping or clicking once selects a row, but previous selections remain selected unless you click or tap it again to de-select it. Un-checking the Multiselect check box will de-select all selected records.




Figure 54. Multiselect and Multiple Selection

Multiple selections are useful for viewing, managing and quickly viewing statistics of only a few selected records.

## 10.5 Automatic Selection of Database Records

When the total number of records in a database becomes very large, manual selection of records can become difficult. The 'Select by' function in the Manage records panel provides the means to select a group of records automatically according to their identifying information. See section 10.8.1 for more information on the Select by database function.

## 10.6 Viewing Detailed Results

Multiple selection is also useful for viewing the detailed results of a few records rather than navigating through the database. After selecting several records, clicking or tapping the 'Details' button on the database panel will display the detailed results for the first record in the selection. Clicking or tapping the 'Up' or 'Down' buttons on the details screen displays the details of the next record in the selection.

## 10.7 Statistics

After selecting a group of records, click or tap the 'Statistics' button on the database panel to view the coefficients of variation (CV) associated with each of the 24 parameters. This function is useful to quickly measure the imprecision of multiple runs of the same sample.

## 10.8 Managing Database Records

'Abacus 5' analyzer database management functions are accessed by clicking or tapping the 'Manage records' button to bring up the Manage records panel.




**Figure 55. Manage Records Panel**

All data manipulation actions that change the data stored in the 'Abacus 5' database are logged in the 'Abacus 5' log. If the multi-user mode is enabled, the name of the operator logged on at the time the data manipulation occurs will also be stored in the log.

### 10.8.1 Select By

Click or tap the 'Select by' button on the Manage Records panel to display the Select by panel. The Select by panel allows you to select a group of records with dates and Result IDs that match criteria entered into the panel.



**Figure 56. Select By Panel**

The Select by panel provides various ways to select records. Only one of the four selection methods can be chosen.

Select day	Clicking or tapping the 'Change' button displays a virtual on-screen date keyboard. All records analyzed on the date entered will be selected.
Select interval	The two 'Change' buttons allow entry of a start and end date. All records on the start date, the end date, and all dates in between will be selected.
Select all	All records in the database will be selected.
Select Result ID	Entering one Result ID will select a record with the same Result ID. Entering two Result IDs will select a range of records including the start and end Result IDs.

Clicking or tapping the 'Run select' will create a multiple selection of all database records that match the criteria (if any). The 'Cancel' button returns to the database view and the current selection is not changed.

### 10.8.2 Importing

Clicking or tapping the 'Import' button displays a blank, modified version of the database panel.




Figure 57. Database Importing Panel


The database importing panel provides the following options:

- The 'Import' button brings up a panel to select previously exported database record files to load and view.
- 'Select by' performs the select by database function to select a group of imported records according to criteria entered by the operator.
- 'Details' views the detailed result information for all or a selection of imported records.
- 'Clear' removes selected records from the imported records table, but does not permanently remove the cleared records from the stored database file.
- 'Back' ends viewing of imported database information and returns to the 'Abacus 5' stored database.

### 10.8.3 Export

The database export panel allows you to select a directory location to store records in a selection for backup or reference purposes. To back up all records in the database, first use the select by function to select all records in the database and export them using the export function.

Clicking or tapping the 'Export' button displays the database export panel.



**Figure 58. Directory Panel for Data Storage**

Introduce a USB storage device (USB hard drive, memory stick, etc.) into an available USB port on the back panel of the 'Abacus 5' analyzer. Select the directory to store the exported data and press the 'Ok' button to start the export.

The file format of exported data is an internal format, and can only be read and interpreted by the 'Abacus 5' analyzer using the Import function. Each database record is stored as an individual file with an '.rp' extension.

### 10.8.4 Send to LIS

The 'Send to LIS' button immediately transmits selected database records to the external LIS according to your LIS settings.

### 10.8.5 Save Tab File

The Save tab file function is used to archive database data to a tab separated text file. The 'Save tab file' button displays a directory panel like Figure 58, except that the file saved in the selected location is a tab separated text file rather than a database backup file. Tab separated text files can be loaded and viewed by a text editor or imported into spreadsheet software.

The filename created by the Save tab file function is generated automatically in the following format: TAB\_YYYYMMDDHHMMSS.txt with the year, month, day, hour, minute, and second of the time when the operation started.

To archive database data using this function, insert a USB data storage device to one of the USB inputs on the back panel of the 'Abacus 5' analyzer. Select the file location of your USB data storage device and click or tap the 'OK' button to begin the export.

See section 0 for detailed information on the format of the tab separated text file.

#### 10.8.6 Save Raw Data

The save raw data function saves the raw data information associated with a selection of database records. Raw data files are only suitable for internal use by Diatron staff. Your service engineer may ask you to save raw data files to help diagnose a service issue.

A confirmation password is required to save raw data. The default factory-set confirmation password is **555**. This password can be changed only by a Diatron certified service engineer.

Insert a USB data storage device to one of the USB inputs on the back panel of the 'Abacus 5' analyzer. Click or tap the 'Save raw data' button on the Manage records panel, then enter the correct confirmation password. Click or tap the 'OK' button to display a directory panel like Figure 58. Select the file location of your USB data storage device and click or tap the 'OK' button to begin saving raw data.

#### 10.8.7 Delete

The Delete function is used to delete records from a selection permanently from the 'Abacus 5' analyzer database. To permanently delete records, select a group of records and click or tap the 'Delete' button on the Manage records panel.

To prevent accidental deletion of data, the 'Abacus 5' requests the same confirmation password used in the save raw data function. After entering the correct password and clicking or tapping the 'Ok' button, the records will be permanently deleted.

The password to allow deleting records is: **555**.



## 11 Calibration

The 'Abacus 5' analyzer must be calibrated in order to provide correct, reliable results. Calibration should be done at the following times:

- At initial installation
- After replacing any component critical to the dilution or measurement process
- Any time quality control shows a systematic error or is outside predefined limits
- At periodic time intervals determined by laboratory regulations

The 'Abacus 5' analyzer calibration process consists of running a commercial calibrator material or a human whole blood sample with known parameter values multiple times. The known parameter values and the average values of multiple runs are used to calculate calibration factors. The new calibration factors take effect once they are accepted by the operator.

If you are using a human whole blood as a calibrator, the calibration target values must be measured on a reference analyzer according to applicable regulations, standards, or laboratory procedures.

An automated calibration procedure is provided that guides the operator through the calibration process and performs the calibration factor calculation. Calibration factors can also be modified directly without running a calibration procedure. Only the parameters that are traceable to reference methods are calibrated: WBC, RBC, PLT, HGB, MCV, and MPV.

To view the calibration functions of the 'Abacus 5' analyzer, select the 'Calibration' icon on the main menu to display the calibration panel. The Calibration panel displays an option to calibrate the 'Abacus 5' analyzer or to view prior calibrations.



Figure 59. Calibration Options

### 11.1 Calibrating the 'Abacus 5'

Selecting the 'Calibrate' option on the calibration options panel displays the calibration panel.






Figure 60. Calibration Panel and Calibration Run Result Panel

The limits of the calibration target values are linked to the normal range of the 'Human' sample mode.



Commercial calibrator material and sample tubes should be handled as a potentially biohazardous material. All applicable laws, regulations, and laboratory practices should be followed in the handling and disposal of calibrator materials.

To perform the automated calibration procedure, ensure that the 'Abacus 5' pneumatic system is started up and that a valid blank measurement has been accepted. Perform the following steps:

- Prepare a sample tube with commercial calibrator material or human whole blood. Ensure that there is sufficient sample volume for three or seven consecutive runs.
- Place the calibrator sample tube in the sample rotor.
- Select 'Three measures' or 'Seven measures' according to your laboratory's standard practices.
- Select 'Calibrator' for a commercial calibrator material or 'Human' for whole human blood calibration.
- Enter the calibration target values provided with your commercial calibrator material or the previously determined parameter values for human blood calibration.
- Click or tap the 'Next' button to begin the calibration process.
- The analyzer displays a blank result screen and begins sample measurement of the calibration material. Wait for the screen to be populated with results when sample measurement is complete.
- Click or tap the 'Accept' button if the measurement run is acceptable.
- Click or tap the 'Start next' button to begin the next calibration run. Selecting 'Start next' without accepting a calibration run discards the run. Selecting 'Abort' will abort the entire calibration procedure and discards all calibration results.
- The calibration procedure ends when three or seven accepted calibration runs are completed. The calibration factors are calculated and presented in the calibrator factor. The previous calibration factors and the coefficient of variation (CV) of the sample runs for each parameter are also displayed.
- Select 'Accept' to replace the current calibration factors with the accepted factors. Selecting 'Discard' eliminates the entire calibration process.

			CV%
WBC	1.19	( 1.00 )	1.94
RBC	1.00	( 1.00 )	0.58
PLT	1.00	( 1.00 )	0.69
HGB	1.00	( 1.00 )	1.15
MCV	1.00	( 1.00 )	0.67
MPV	1.00	( 1.00 )	5.07
RDWcv	0.92	( 0.92 )	0.67

Discard Accept

AS state None 100% 100% 100% 0% hu-HU 17:05

**Figure 61. Calibration Factor Panel**


The automatic calibration process can be skipped by selecting the 'Manual' calibration mode on the calibration panel and clicking or tapping the 'Next' button proceeds directly to the calibration factor panel for manual entry of calibration values. Manually entered calibration values can be accepted in the same way as calculated factors from calibrator sample runs, and are treated by the 'Abacus 5' in the same way.

### 11.1.1 View Calibrations

The view calibrations feature is helpful for reviewing the calibration history for the analyzer. Increasing or decreasing calibration value trends over time may indicate a problem that requires service intervention.

To view the calibration history click or tap the 'View calibration' function on the calibration panel. The calibration history is presented in a table view that can be navigated in the same way as the database panel.

Any calibration can be deleted. Deleting the last calibration causes the 'Abacus 5' analyzer to revert to the last accepted calibration.



**Figure 62. View Calibration Panel**

To delete a calibration, click or tap the 'Delete' button. The system will ask you for a confirmation password to avoid accidental deletion of calibration data. The confirmation password is set at the factory to **555**, and can be changed by your Diatron certified service engineer.

## 12 Quality Control

Daily use of quality control (QC) ensures that the 'Abacus 5' analyzer is operating optimally. Stable recovery of standard QC material parameters assures the operator that the 'Abacus 5' is operating consistently on a day-to-day basis and is functioning correctly.

By analyzing control materials on a regular basis, day-to-day reproducibility and general condition of the analyzer can be monitored. Target values and acceptable (tolerance) ranges for each parameter can be specified for an unlimited number of QC materials.



Control materials must be used before their expiration date. Always observe the control manufacturer's instructions for storage and use. Control materials must be well mixed before use.

Diatron recommends the use of Diacon 5 controls from Diatron.

The 'Abacus 5' analyzer does not limit the number of QC material lots or the number of QC measurement runs that are stored in the control database. Each QC measurement run is saved to the selected QC material lot.

QC lots can be commercial control materials or a whole human blood sample whose parameter values and ranges are determined by other means. The advantage of a commercial control material is that its shelf life is much longer than human blood, allowing control of the analyzer for a longer period of time.


Controls in the 'Abacus 5' are processed in manual mode only. To access the QC panel, select the 'QC' icon on the main menu.



Figure 63. QC Panel

### 12.1 Set QC Reference

Before taking any QC measurements, target values and acceptable ranges must be specified for a QC lot of control material. Click or tap the 'Set QC reference' button to display the Set QC panel.




The screenshot shows the 'Set QC Reference' panel. At the top, there are icons for 'Measure', 'AS', 'Database', 'Print', and 'Menu'. Below these are input fields for 'QC name', 'LOT No.', and 'Expiration date'. The main area contains a grid of parameters for blood counts, each with a target value and tolerance range input field. The parameters are: WBC, NEU%, LYM%, RBC, HGB, HCT, MCV, PLT, PCT, MPV, MON%, EO%, BAS%, MCH, MCHC, RDWsd, RDWcv, PDWsd, and PDWcv. At the bottom, there are buttons for 'Load QC reference', 'Save reference', and 'Back'. A status bar at the very bottom shows 'AS state' with four colored indicators (green, yellow, orange, red) and values (65, 5, 96, 11), and a timestamp 'hu-HU 16:29'.

**Figure 64. Set QC Reference Panel**

The 'QC name' is the name that will be displayed for selecting the control to associate with a QC measurement. Enter the identifying information, and the target values and tolerance ranges for any parameters that are to be controlled with the 'Abacus 5' QC system. To exclude a parameter from the QC procedure, leave its target value and tolerance range field empty. Press the 'Save reference' button after data entry is completed.

Tapping the 'Load QC reference' button will allow filling QC material data electronically – using a file provided by the manufacturer of the analyzer. A dialog will open. Browse to the correct folder, and select the file by tapping the checkbox in front of the filename.

The file contains information for all three levels of the control material (Low, Normal and High). All three levels will be loaded into the system and will become selectable in the QC menu.



The screenshot shows the 'Load QC Reference' dialog. It has a 'Selected path' field showing 'E:\diatron\DCheck\_U011.as' and an 'Actual path' field showing 'E:\diatron'. Below these is a list of files with checkboxes. The file 'DCheck\_U011.as' is selected. At the bottom, there are buttons for 'Level up', 'Create directory', 'Cancel', and 'Ok'. There are also up and down arrow buttons.

**Figure 65. Load QC Reference Panel**

## 12.2 QC Measure

To run a QC measurement, please ensure that a valid blank measurement has been run and accepted. Associate the QC measurement to the correct control material lot by select the matching QC name in the selector below the 'QC measure' button on the QC panel. Place the control material sample tube in the sample rotor and click or tap 'QC measure'.

The QC measurement results will be added to the QC database and associated with matching QC material lot.

## 12.3 View QC References

Click or tap the 'View QC References' button to see a list of QC lot references stored in the 'Abacus 5' analyzer control database. Target values and ranges can be viewed for any stored QC references.

## 12.4 View QC Data

Selecting 'View QC data' displays QC runs stored in the QC database in a table format that operates the same way as the 'Abacus 5' analyzer database. QC measurement results will have sequential ID numbers. Option buttons are provided to switch back and forth between viewing QC data and viewing QC diagrams.

## 12.5 View QC Diagrams

Click or tap 'View QC diagrams' to see the QC diagrams for control measurements.




Figure 66. View QC Diagrams Panel

Means, standard deviations (StDev) and coefficients of variation (CVar) are calculated based on the QC analyses. The dotted lines delineate acceptable ranges on QC diagrams. Use the selector control at the top of the panel to view graphical control information for other lots of control materials.

## 12.6 View X-B Data

Abacus 5 offers X-B (a.k.a: X-bar,  $\bar{X}$ ) to facilitate trend analysis of samples run.

The X-B moving average is an algorithm which is capable of tracking changes of the accuracy-alternating of the measurements. The algorithm observes the alteration of the derived values which are generally stable.

For example the MCV, MCH and MCHC values show the proportion between the blood components. These values fall within narrow limits in human blood therefore these values are able to show the drifting of the device.

The X-B algorithm works with batches. There are 20 samples in the batches. The average deviation of the corresponding batch is calculated as follows:

The X-B algorithm is illustrated below:

$$X_{(B,i)} = X_{(B,i-1)} + \text{SGN} \left\{ \frac{\sum_{j=1}^N \text{SGN} [X_{(j,i)} - X_{(B,i-1)}] * \sqrt{[X_{(j,i)} - X_{(B,i-1)}]^2}}{N} \right\}$$

Where:	i	=	present batch
	j	=	present specimen
	$X_{(B,i)}$	=	ith X-B value (current batch)
	$X_{(B,i-1)}$	=	(i-1)th value (previous batch)
	SGN	=	the arithmetic sign or number in parenthesis
	N	=	the number of samples in the batch
	*	=	symbol used to represent multiplication

Divide the average deviation of the corresponding batch by the number of the items in the batch (N = 20). Add the square of the resulting deviation to the mean. Then the result represents a point on the diagram. The value of the mean will be changed with all counted average deviation from the corresponding batch.

The average of the first batch (X) equals the target value. All the average deviations belonging to the corresponding batch are added to the mean with correct signs.

The samples are used in calculation chronologically therefore the copulated deviation points show the deviation from the mean line.

## 12.7 View X-B Diagrams

View X-B diagrams displays X-B charts for MCV, MCH, and MCHC.

View X-B data function displays X-B information in a table format similar to view QC data.

## 13 Patients

The 'Abacus 5' analyzer allows the operator to define and store patients in the patient database during manual measurement. Patients can also be predefined while creating a list for one of the automatic processing modes.

Select the Patients icon on the main menu to display the Patients panel.



Figure 67. Patients Panel

The Patients panel displays patient information in the familiar table format of the 'Abacus 5' analyzer database, and operates the same way as the database panel. Use the 'New' function to define new patients, and 'Edit' to change existing patient information. 'Details' displays patient data in a read-only mode.

The Patient panel can also be accessed by clicking or tapping the Name or Patient ID field in the Measurement panel. Two additional function buttons are available if the Patient panel is accessed from the Measurement panel.

The 'Select ID' function associates the currently selected patient to the sample measurement and returns to the Measurement panel. The 'Cancel' function returns to the measurement panel without selecting any patient to associate to the sample measurement.

You can select either Name or Patient ID to be displayed as the patient identifier in the Measurement panel by changing the settings in the 'Patient's displayed data' in the Main Menu/Settings/Customize panel.



Measure AS Database Print Menu

Name

Title

Patient Unique No.

Gender

Date of birth

Nationality

Phone

E-mail

Country

City

Address

ZIP

AS state  70% 80% 100% 15% en-US 04 : 24

**Figure 68. Edit Patient Panel**

The 'Save' function saves any changes made. To discard any changes made, simply click or tap the 'Back' button.

## 14 Multi-User Mode

The multi-user mode provides the clinical laboratory manager the means to control access to the 'Abacus 5' analyzer by ensuring that only authorized users can operate the analyzer. With the multi-user mode, the laboratory manager can assign some users the authority to perform critical system functions such as maintenance and calibration, and restrict other users to routine activities such as running samples.

The multi-user mode also provides audit tracking of activities by recording the user ID of the currently logged on user in the event log. Events that are normally tracked in the event log such as maintenance and calibration will now also store the user ID of the user that initiated the event. Logon and logoff become new events that are tracked in the system log when the multi-user mode is enabled.

The multi-user mode is disabled by default, and can only be enabled by your Diatron certified service engineer. If you want to use the multi-user mode, ask your Diatron certified service engineer during the analyzer installation (or at any time) to create the first administrator user with the user ID and password you provide. This allows you administrator access to the 'Abacus 5' analyzer and permits you to create additional users.

### 14.1 Types of Users

The multi-user mode provides two types of users. Administrator users have access to all system functions except service functions, and a regular user has access to a limited set of system functions such as running samples. Your Diatron service engineer creates the first administrator user account for the laboratory manager. As an administrator user, you can create other administrator user accounts or regular user accounts.

All administrator users have the same rights and privileges of the original administrator user, including the ability to create other administrator users. Use caution in assigning administrator user privileges and restrict the number of administrator users to trusted, well trained users that have been fully trained on the operation and maintenance of the 'Abacus 5' analyzer.

The following table shows the menu items that are restricted for each user type.

Item	Administrator	Regular user
Measure	No restrictions	No restrictions
Auto Sampler	No restrictions	No restrictions
Database	No restrictions	No access to 'Delete' button in 'Manage Records'
Print	No restrictions	No restrictions
Menu – Patient	No restrictions	No restrictions
Menu – Exit	No restrictions	No restrictions
Menu – Settings	No restrictions	No access except 'Customize' button
Menu – Maintenance	No restrictions	No restrictions
Menu – Calibrations	No restrictions	No access to 'Calibrate' button
Menu – QC	No restrictions	No access to 'Set QC Reference' button No access to 'Delete' button in 'View QC Data' No access to 'Reject' button in 'View X-B data' No access to 'Undo Rejection' button in 'View Rejected' in 'View X-B Data'
Menu – Diagnostics	No restrictions	No restrictions

Table 17. Restrictions by User Type

### 14.2 Creating a New User

When the analyzer is started with the multi-user mode enabled, the following logon panel is displayed.

The Logon Panel is a light gray rectangular window. At the top, it has a label 'UserID' above a yellow rectangular input field. Below this is a label 'Password' above a white rectangular input field. At the bottom, there are three rounded rectangular buttons: 'Add user by administrator', 'Add user by service person', and 'Login'.

**Figure 69. Logon Panel**

Enter your user ID and password, then press or click the 'Add user by administrator' button. This will bring up the add user panel.

The Add User Panel is a light gray rectangular window. At the top, there is a list area with a header 'labmanager' and a large empty space below it. To the right of the list area are two square buttons with up and down arrows. Below the list area are two buttons: 'Select' and 'Close'. At the bottom, there are three input fields labeled 'UserID', 'Password', and 'Re-enter password'. To the right of these fields is a checkbox labeled 'Administrator'. At the bottom right is a 'Create' button.

**Figure 70. Add User Panel**

To create a new user, enter a new user ID in the 'UserID' text field, then enter the new user password into the 'Password' text field, then again in the 'Re-enter password' text field. If you want to create a administrator regular user, click or tap the 'Administrator' check box. To create a regular user, leave the 'Administrator' check box unchecked. Click or tap the 'Create' button when you have finished entering the user information. This will display a panel indicating the operation is complete.


The Add User Complete Panel is a light gray rectangular window. It has a label 'Complete' at the top left. At the bottom right is a rounded rectangular button labeled 'Ok'.

**Figure 71. Add User Complete Panel**

Click or tap the 'Ok' button to acknowledge. You can continue to create users, or press or click the 'Close' button when you are finished creating users. Pressing the 'Close' button will display the logon panel again.

## 14.3 Deleting or Changing Users

Anyone with a valid administrator account can delete administrator or regular delete users. Regular users do not have privileges to create or delete users. To delete a user, log on with a valid administrator account user ID and password, then click or tap the 'Menu' button at the top of the screen, then click or tap the 'Settings' menu item.




**Figure 72. Settings Panel (Administrator User)**

Click or tap the 'User' button to access the user management panel. Select the user you want to delete and click or tap the 'Delete user' button. Click or tap the 'Yes' button when asked to confirm the deletion of the selected user.



**Figure 73. Settings Panel (Administrator User)**

You can also change the access privileges of a user. Unchecking the 'Administrator' check box when an administrator user is selected will change the user to a regular user. Similarly, checking the 'Administrator' check box when a regular user is selected will change the user to an administrator user. Either way, you are asked to confirm the change. Click or tap 'Yes' button to confirm or 'No' to cancel the change.



**Figure 74. Confirm User Delete Panel**

## 15 Settings

The 'Abacus 5' analyzer can be tailored to your laboratory's needs by changing the settings. Click or tap the 'Settings' icon on the main menu to access the Settings panel.



**Figure 75. Settings Panel**


The options available on the Settings panel are:

- Customize
- Laboratory
- External Devices
- System
- Units
- Printer
- Profile Limits
- X-B
- User

Clicking or tapping one of the buttons in the Settings panel displays the selected panel where settings changes can be made. Each panel has a 'Back' and 'Save' button. The 'Save' button saves any changes made and returns to the Settings panel, and the 'Back' button discards any changes and returns to the Settings panel.

### 15.1 Customize Settings

The Customize panel allows you to customize the 'Abacus 5' analyzer user interface language and other preferences about the analyzer operation. Click or tap the 'Customize' button on the Settings panel to access the Customize panel.



**Figure 76. Customize Settings Panel**

On the Customize panel, you can change the following settings:

- Skins – By changing the ‘skin’ the appearance of the menus, GUI elements will change. The ‘skins’ has no effect on the behavior of the ‘Abacus 5’
- Language – select the language that the ‘Abacus 5’ analyzer user interface uses to communicate with the operator. Only English (United States) is available in the US version of the ‘Abacus 5’ analyzer.
- Limit style – select Number, Graphical1, and Graphical2 display modes to change the display method for parameters outside the normal ranges on the result screens and printed reports.
- Sound volume – modify the sound volume up or down
- On screen keyboard active – Select whether to use the on-screen virtual keyboards for displaying information, or an external keyboard.
- Patient’s displayed data – choose Name or Patient ID as the patient identifier displayed in the Measurement panel.
- Demo mode

## 15.2 Laboratory Settings

The Laboratory settings panel allows entry of seven lines of text up to 50 characters of descriptive information about the laboratory such as laboratory name, address, etc. All non-empty lines are displayed at the top of every result printout.


## 15.3 External Devices

The External Devices panel allows entry of external LIS settings. The ‘Abacus 5’ analyzer supports external LIS communications over a serial communications link using the Diatron version 3.1 protocol. It also supports an Ethernet communications link over the HL7 version 2.5 protocol.

The External Devices panel lets you change the following settings:

- Sending port baud rate – select 9600 or 115200 baud for serial LIS connection.

- Automatic LIS – The ‘Abacus 5’ will transmit every measurement result as soon as it is available if this check box is checked.
- LIS – the ‘Abacus 5’ analyzer will use an Ethernet connection rather than a serial connection if this check box is checked.
- IP – the IP address of the external LIS computer should be entered here if an Ethernet LIS connection is used. IP addresses consist of four numbers from 0 to 255 separated by periods. See your IT administrator for more information about the IP address of your host LIS computer.
- Port – the port number of the external LIS computer should be entered with here with no commas if an Ethernet LIS connection is used. Ports are a single number between 0 and 65535. See your IT administrator for more information about the port number of your host LIS computer.
- Bidirectional LIS – check this check box if the LIS computer will download lists for one of the list modes for automatic sample processing.



**Figure 77. External Devices Settings Panel**

NOTE: When using bidirectional LIS settings, the ‘Abacus 5’ analyzer “listens” for an external device (or LIS computer) to send it instructions. Consult your information technology staff regarding the configuration of the host computer and the ‘Abacus 5’ analyzer for proper LIS operation. Your information technology staff person may need to consult your Diatron certified service engineer to ensure the Internet Protocol (IP) settings of the ‘Abacus 5’ analyzer are set properly.

The ‘Abacus 5’ analyzer “listens” on port 6600 for incoming LIS instructions.

## 15.4 System Settings

System settings change various options about the operation and behavior of the ‘Abacus 5’ analyzer.

The System settings panel lets you change the following settings:

- Database columns order – Which column (parameter) of the measurement database is displayed/ hidden
- Waste container volume – select 10L, 20L, or None. Selecting 10L or 20L changes the waste container tracking display behavior in the reagent status area at the bottom of the screen. Selecting None assumes waste is being drained directly into the laboratory drain and disables waste container tracking.



- Database display limit – select All or Last month to change the number of database records displayed in the database panel. Last month only displays the last 30 days of measurement results, and All displays all measurement results stored in the 'Abacus 5' database.
- Use only Sarstedt Monovette tube from sample rotor – Check this check box to change the sampling depth for manually presented tubes at the sample rotor to accommodate Sarstedt Monovette tubes. Uncheck this check box to manually process standard sample tubes. This setting does not affect automatic sample processing with the optional Autosampler, which is able to recognize tube types and change sampling depth automatically.
- Standby time – This defines the time period of inactivity before the analyzer automatically leaves the ready state and drains the pneumatic components to prepare for overnight standby. If this time is set to zero, the analyzer will never enter the standby state automatically.
- Offline rinsing frequency – this defines the time period between automatic fluidic operations to ensure that fresh reagents are available in critical pneumatic components so that the analyzer can quickly return to the ready state from standby. If this time is set to zero, the analyzer will never perform automatic fluidic operations while in the standby state, but may increase the amount of time needed to return to the ready state from standby.
- Screen saver – this defines the time period of inactivity before the analyzer automatically starts the screen saver. Entering a time of zero prevents the analyzer from every starting the screen saver.
- Special flags (G, L) – Checking the check box disables the generation, storage, LIS transmission and display of the G and L flags on the screen, printout.




Figure 78. System Settings Panel

## 15.5 Units

The Units panel allows the user to choose which units are displayed for certain parameters.

The following options are available on the Units panel:

- HGB unit – choose g/dL or g/L to change the displayed units for the HGB and MCHC parameters.
- Count unit – Choose cells/ $\mu\text{L}$  or cells/L to change the displayed units for the parameters in the table below.

Parameter	Cells/ $\mu\text{L}$	Cells/L
WBC	$10^3/\mu\text{L}$	$10^9/\text{L}$
LYM	$10^3/\mu\text{L}$	$10^9/\text{L}$
NEU	$10^3/\mu\text{L}$	$10^9/\text{L}$
MON	$10^3/\mu\text{L}$	$10^9/\text{L}$
EO	$10^3/\mu\text{L}$	$10^9/\text{L}$
BAS	$10^3/\mu\text{L}$	$10^9/\text{L}$
PLT	$10^3/\mu\text{L}$	$10^9/\text{L}$
RBC	$10^6/\mu\text{L}$	$10^{12}/\text{L}$

Table 18. Count Unit Parameters

- PCT/HCT unit – Choose Abs/Percentage/ to change the displayed units for PCT/HCT.
- MCV unit – Choose fL/ $\mu\text{m}^3$  to change the displayed units for MCV




Figure 79. Units Settings Panel


## 15.6 Printer Settings

The Printer Settings panel lets you choose which installed printer to use, view the status of the printer, and change the print behavior of the analyzer during operation.

The following options are available on the Printer Settings panel:

- Printer – select an installed printer to use.

- Printer status – indicates the status of the currently selected printer.
- Color printing – check this check box to print in color, or uncheck to print in black and white.
- Double sided printing – check this check box to print on both sides of the paper, or uncheck to print on one side only. This option is only available if your installed printer supports double-sided printing.
- Items in queue – indicates how many printouts are currently waiting in the printer queue.
- Cancel all jobs – Click or tap this button to cancel printing and clear all items waiting in the printer queue.
- Pre-printed paper header height – select from 0 to 80 pixels.
- Automatic print – check this check box to print every sample or control as it is run.
- Logo visible – check this box to print the Diatron logo on every printed report.
- Warning flags are visible – check this box to print warning flags on every printed report
- Diagnostic flags are visible – check this box to print diagnostic flags on every printed report



**Figure 80. Printer Settings Panel**

## 15.7 Profile Limits Settings

The Profile Limits settings are used to view and modify the normal ranges associated with the 'Abacus 5' analyzer's five sample modes: Human, Male, Female, Alternate 1, and Alternate 2. NOTE: Alternate 1 and Alternate 2 are sample modes that are different than Human, Male and Female, and have user-definable normal ranges.

The control at the top of the panel selects one of the five sample modes. The panel is then populated with the lower and upper limit values for each of the 24 parameters. Modify the values as needed, then click or tap 'Save' to save the new values and return to the Settings panel, or select Back to discard any changes and return to the Settings panel.

Parameter	Lower Limit	Upper Limit	Unit
WBC	3.00	15.00	10 <sup>9</sup> /μL
NEU	1.50	7.00	10 <sup>9</sup> /μL
LYM	1.00	3.70	10 <sup>9</sup> /μL
MON	0.00	0.70	10 <sup>9</sup> /μL
EO	0.00	0.50	10 <sup>9</sup> /μL
BAS	0.00	0.15	10 <sup>9</sup> /μL
NEU%	37.00	72.00	%
LYM%	21.00	50.00	%
MON%	0.00	14.00	%
EO%	0.00	6.00	%
BAS%	0.00	1.00	%
RBC	3.50	5.50	10 <sup>6</sup> /μL
HGB	120	174	g/L
HCT	26.0	50.0	%
MCV	86.0	110.0	fL
MCH	25.0	38.0	pg
MCHC	300	350	g/L
RDWs	0.0	0.0	fL
RDWc	0.0	16.0	%
PLT	50	400	10 <sup>9</sup> /μL
PCT	0.13	0.43	%
MPV	9.0	13.0	fL
PDWs	0.0	0.0	fL
PDWc	0.0	0.0	%


**Figure 81. Profile Limits Settings Panel**

## 15.8 X-B Settings

The X-B panel displays the currently selected X-B constants. Change the values as needed, then click or tap 'Save' to save the new values and return to the Settings panel, or select Back to discard any changes and return to the Settings panel.

The following X-B settings are available:

- Lower and upper limit – A parameter is usable in X-B statistics if its value is between these criteria
- Target –The start value of the first batch in the X-B algorithm
- Action limit(%) – Maximum percent deviation of the running average values from the target value. The user needs to check the X-B diagrams regularly as no action is performed by the software if the action limit is reached



**Figure 82. X-B Settings Panel**

## 15.9 User Settings

The User Settings panel allows administrators to delete users, or assign administrator rights.




**Figure 83. User Settings Panel**

See section 1 for more details on creating, deleting and modifying users.

## 15.10 Date and Time Adjustment

Double-click or tap the time display on the bottom right side of the screen to display a time/date settings panel. Change the time and date if needed, then click or tap the 'OK' button to save the new settings. Click or tap 'Cancel' to exit the Time/Date settings without saving.



**Figure 84. Time/Date Settings**

The 'Abacus 5' analyzer does not support automatic time zone or daylight savings time changes. Please change the time and date manually when necessary.



# 16 Instrument Diagnostics

The instrument diagnostics section allows the operator to perform diagnostic procedures, check the operational history of the instrument, set or check the status of the reagents and view stored information about the system.

Click or tap the 'Diagnostics' button on the main menu to view the available options in the Diagnostics panel.



Figure 85. Diagnostics Panel

## 16.1 Analyzer Self Test






Figure 86. Self Test Panel

Initiating a self test will test the key functions and components of the analyzer. The process takes approximately one minute. The individually tested items and their acceptable ranges can be found in the following table. Clicking or tapping any of the individual test items provides detailed information including the item readings, limits and comments.

Clicking or tapping the 'Load last selftest' button lets you view the results of the previous self test for comparison.

Big and small buffer times for 0 to 1,000m (0 to 1,380 ft.) Altitude			
Big buffer time		Small buffer time	
Generate	3,000 to 12,000 ms	Generate 1	1,200 to 2,800 ms
Release	3,000 to 6,000 ms	Generate 2	1,200 to 2,800 ms
		Release	700 to 1,300 ms
Big and small buffer times for 1,000 to 2,000m (3,280 to 6,560 ft.) Altitude			
Big buffer time		Small buffer time	
Generate	3,000 to 16,800 ms	Generate 1	1,200 to 3,640 ms
Release	3,000 to 7,200 ms	Generate 2	1,200 to 3,640 ms
		Release	700 to 1,560 ms
Big and small buffer times for greater than 2,000m (6,560 ft.) Altitude			
Big buffer time		Small buffer time	
Generate	3,000 to 26,400 ms	Generate 1	1,200 to 5,040 ms
Release	3,000 to 10,200 ms	Generate 2	1,200 to 5,040 ms
		Release	700 to 2,210 ms
Big buffer drift		Small buffer drift	
Maximum	540 to 560 mbar	Maximum	225 to 235 mbar
Minimum	530 to 560 mbar	Minimum	215 to 235 mbar
Drift	-5 to +15 mbar	Drift	-5 to +15 mbar
Max vacuum P1		Max vacuum P2	
Maximum	550 to 1,000 mbar	Maximum	500 to 1,000 mbar
Time	0 to 6,000 ms	Time	0 to 6,000 ms
Pump Status		Null pressures	
Pump1	1 to 1	Sheath	-20 to +20 mbar
Pump2	1 to 1	Capillary	-20 to +20 mbar
		Chamber	-20 to +20 mbar
Noise/pulse		HGB head	
pls/8sec	0 to 20 pulses	HGB dark	0 to 3,000 pulses
20000 pulses	19,990 to 20,050 pulses	HGB light	3,000 to 60,000 pulses
TCU temperature		Laser parameters	
Actual	Reference $\pm 0.2^{\circ}\text{C}$ (0.35°F)	Reference	$\pm 10\%$
Sink	0 to 70°C (32 to 158°F)	Laser Power	32 to 128
Laser DC Level		Reagent lock	
Laser off	0.00 to 0.05 mV	Status	UNLOCKED / LOCKED / NA
Laser on	0.05 to 0.65 mV	Available measurements	No specification
Electrode		Battery/Power Supply	
Voltage	45 to 55 V	Battery Voltage	0.0 to 3.5 V
Current	620 to 680 $\mu\text{A}$	+12V	11.4 to 13.0 V


Offset	-5.0 to 5.0 mV	-12V	-13.0 to -11.4 V
Temperature			
WBC Pre-heater	None (LSDACQ v3.0)		
	34.0 to 38.0°C (LSDACQ v3.1)		
LSDACQ	No specification		

**Table 19. Self Test Item Limits**

If any of the values falls outside the above limits, the software will list the errors on the self test report screen.

## 16.2 Log

The system log provides information for service or troubleshooting purposes. It lists events and possible problems related to operation or measurements. It also lists routine events such as login. If the multi-user mode is enabled, the user ID of the operator that is logged at the time an event occurs will also be stored with the log entry.



**Figure 87. Log Panel**

The selector at the top of the panel filters the log entries according to the following categories:

Filter	Filter Description
All	Displays all log entries, no filter
Error	Displays only error log entries
Pneumatic	Displays only pneumatic system errors
Software	Displays only software system errors
User	Displays only user-related log entries
System	Displays only system action log entries

Table 20. Log Filters

The 'Details' button provides more information about the currently selected log entry. The 'Selftest' button provides direct access to the self test panel in the event a maintenance procedure is needed in response to a log entry.

## 16.3 Reagent Status

The reagent status panel displays the level of each of the on-line reagents and waste containers. It provides the means to reset the levels of each of the reagents in the event of a reagent container change. It also provides the means to reset the waste container level in the event the waste is emptied. A 'Reset all' button is also provided in the event of a change of all reagent containers and emptying of waste.



Figure 88. Reagent Status Panel

The 'Abacus 5' analyzer has built-in reagent consumption calculators. The calculated level of each reagent is graphically and numerically indicated on the reagent status panel. The current reagent levels are also indicated in the status area at the bottom of the screen.

The 'Abacus 5' analyzer notifies the user when any of the reagents are running low or when the waste container is almost full. If you receive a reagent low notification, check the reagent container and replace it if necessary. Click or tap the 'Reset Diluent,' 'Reset Lyse-5P' or 'Reset Diff-5P' button reset the software counters to a full container level. Similarly, if you receive a notification that the waste is almost full, empty the waste container and click or tap the 'Reset Waste' button on the reagent status panel to reset the software counters to an empty waste container level.

If the reagent-locking feature is activated, please remove the old hardware key and insert the new hardware key provided in the Lyse-5P package before clicking or tapping pressing the 'Reset Lyse-5P' button. Clicking or tapping the 'Reset Lyse-5P' button activates the uploading of additional measurement counts from the hardware key.

16.4 Statistics

This menu item will display cumulative information such as the measurements and number of error events experienced during operation. Statistical data can be only reset by a Diatron certified service engineer.




Figure 89. Statistics Panel

16.5 Information

This menu item displays the instrument serial number and software versions of various system components. This data is important for maintenance purposes, and your service engineer will ask for this information prior to a service call.



**Figure 90. Information Panel**

## 17 Maintenance

To ensure reliable operation, regular maintenance is required. This chapter describes necessary procedures and steps required to keep your analyzer in a good operating condition.



Only Diatron certified service personnel that have received Abacus 5 service training can resolve major problems that occur during operation. Electrical and moving parts should always be serviced by trained service personnel.

There are, however some maintenance actions the operator should perform to keep the analyzer in good operating condition



Any elements inside the analyzer like tubes, valves, chambers, containers should be handled as a potentially biologically and chemical dangerous material. Observe local laws and regulations regarding handling and during disposal.

### 17.1 Opening the Front Panel

Several components that require regular cleaning and maintenance are behind the front panel.

Make sure that nothing is placed on top or in front of the analyzer. Grasp the lower sides of the front panel, gently pushing the sides. Pull the lower part towards you, lifting up the front panel.

Upon opening, a lever becomes visible. Make sure to tilt the front panel upwards so that you can push the lever into the secure position.




Figure 91. Opening the Front Panel

### 17.2 Closing the Front Panel



Do not push on the front panel while the security lever is in the “lock” position.




Figure 92. Closing the Front Panel

Gently lift the front panel so that the security support lever can be moved to the free position.

Gently lower the front panel. When it reaches its lowest position, gently push on the front side to click the lock-levers in place.

### 17.3 Removing the Side Panels

With the front panel open and secured in the upright position, thumb screws become visible that fix the side panels to the instrument: 2 in the front and 2 in the back. The screws are designed to stay on the analyzer to avoid losing them.

The system valves and measuring chambers are behind the left side panel. Open the left side panel if chamber cleaning is requested by the analyzer.

The sampling unit, wash head and sampling needle are behind the right side cover. Open the right side cover if the wash head needs to be cleaned or replaced.

When closing the side panels, always make sure to re-tighten the thumb screws.

	<p>After removing the covers potentially hazardous parts become accessible such as electronic boards, motors, moving parts, the sampling needle, chambers, tubes and valves.</p>
	<p>The sampling needle and other components inside the analyzer may cause injury or can be damaged if handled incorrectly. Use extreme caution when operating the analyzer with the front cover open. Running measurements with an open cover is not recommended due to the risk of injury. Always wear safety gloves while performing maintenance actions.</p>

### 17.4 User Maintainable Parts of the Analyzer

There are three user maintenance points inside the 'Abacus 5' analyzer:

- Shear valve: this critical component and is responsible for correct sampling and dilution.
- Wash head: this component keeps the sampling needle clean ensuring reliable piercing and sampling.
- Measuring chambers: contamination induced problems such as noise or high blank values can usually be resolved by cleaning the plastic measuring chambers.

Always wear protective gloves when working on internal parts of the analyzer.

## 17.5 Software Maintenance Functions

Maintenance functions can be accessed by clicking or tapping the “Maintenance” main menu item. This opens the maintenance function panel.



Figure 93. Software Maintenance Functions

Clicking or tapping any of the function buttons initiates the corresponding maintenance functions. The maintenance functions are grouped into the following categories.

<b>Cleaning</b>
Bleaching: cleans the diluent buffer and the connecting tubes with external cleaner reagent introduced through the sample rotor
Flow cell backwash: cleans the optical head and the flow cell by draining the sheath buffer through the optical head with a frequently used back flush through the flow cell, then refills and rinses the system
Shearvalve cleaning: drains the shear valve so the operator can open it to clean the ceramic valve, then refills and rinses the system
Rinse: washes the tubing system with diluent
Cleaning: cleans the system with external cleaner reagent introduced through the sample rotor
Hard cleaning: extended cleaning of the system with external cleaner reagent introduced through the sample rotor
<b>Draining</b>
Drain Diff-5P: drains the Diff-5P buffer
Drain sheath: drains the sheath buffer
Drain Lyse-5P: drains the Lyse-5P buffer
Drain diluent: drains the diluent buffer
Drain All: Drains all of the above buffer
<b>Empty Chamber</b>
Empty Mix chamber: <help needed with this one. What is the “mix” for?>
Empty RBC chamber: empties the RBC counting chamber
Empty WBC chamber: empties the WBC counting chamber
<b>Prime</b>
Prime Diff-5P: primes the Diff-5P buffer and tubing with Diff-5P reagent



Prime sheath: primes the sheath buffer and optical head with diluent
Prime Lyse-5P: primes the Lyse-5P buffer and tubing with Lyse-5P reagent
Prime diluent: primes the diluent buffers and tubing with diluent or primes the entire system with diluent (the operator is given a choice by the software maintenance function)
Prime all: primes all reagent buffers with their corresponding reagents (the operator is given a choice to perform a full system prime or prime individual buffers sequentially)
<i>Fill</i>
Fill: initial filling of the system with reagents. This must be performed PRIOR to first use of the instrument after installation or after the system has been drained for transportation or storage.
<i>Touchscreen</i>
Touchscreen calibration: calibrates the touch screen by prompting the operator to tap specific areas of the display.

Table 21. Software Maintenance Procedures

## 17.6 Cleaning the Shear Valve

The inner surfaces of the shear valve builds up salt deposits from the reagents. If these deposits are not regularly cleaned, it will cause a malfunction. To avoid this, it is recommended that the operator clean the shear valve after every 1200 samples.

The following materials are required to perform a shear valve cleaning:


- A clean, soft, damp cloth
- A clean, soft, dry cloth
- Lint-free, soft, dry-cloth
- Gloves
- Tweezers
- Toothpicks
- Water



The shear-valve is in contact with the sample blood. Wear gloves while cleaning the shear-valve. Handle the materials use to clean the shear-valve as potentially infectious material.


The shear valve is fixed with a thumb screw. No screwdriver or other tools are required.

Click or tap the 'Maintenance' button on the main menu to bring up the maintenance functions panel. Click or tap the 'Shearvalve cleaning' button in the Cleaning section to initiate the shear valve cleaning maintenance function.



**Figure 94. Shearvalve Cleaning Maintenance Function**

The 'Abacus 5' analyzer asks for confirmation to start the procedure. After clicking or tapping the 'OK' button the shear valve and the connecting tubing are emptied of liquids. A small amount of liquid may remain inside the tubing and the shear valve. When the procedure completes, the 'Abacus 5' prompts the operator to proceed with the shear valve cleaning with the following message:



**Figure 95. Shearvalve Cleaning Maintenance Function Prompt**

Don't click or tap the 'OK' button until the shear valve has been cleaned and re-assembled! Perform the shear valve cleaning according to the instructions below.

1. Open the front cover and secure it with the latch. Locate the shear valve in the center of the analyzer.
2. Unscrew and remove the axis thumbscrew that secures the shear valve. Clean the axis thumbscrew with water and wipe it dry.



Figure 96. Removing the Shear Valve Axis Thumbscrew

3. Slide off the upper disk of the shear valve. Due to the extremely smooth surface of the ceramic discs the may stick together somewhat. If the shear valve was not in use for a few days, then apply a few drops of water to the contact are between the upper and lower disks. The salts will dissolve in a short time and release the upper disk.




Figure 97. Disassembling the shear valve

4. Gently clean the mating surfaces of the disks of the shear valve, the housing of the valve and the tube connections.
5. Remove any salt build-up. Use tweezers to push in the dampened and the dry cloths. Apply a few drops of water to soften any hardened deposits.
6. Use a tooth pick to remove any salt crystals from narrow places. Do not use any sharp metal or hard objects which can scratch the smooth surfaces of the shear valve.
7. Clean the outside of the shear valve and the housing and mounting parts of the shear valve. Pay attention to clean the aligning surface.
8. Ensure that no lint or fibers remain on the smooth mating surfaces of the ceramic disks.




Figure 98. Cleaning the Shear Valve

9. After cleaning the shear valve, housing and surrounding area, gently put the two disks together.
10. Put the axis thumbscrew back into the center of the upper disk. The spring on the axis thumbscrew guarantees the necessary closing force for the two disks. Gently press down and rotate the axis thumbscrew clockwise until it clicks into the lower part.
11. Twist the axis thumbscrew until it stops. The mechanical design of the thumbscrew prevents over-tightening.




Figure 99. Reassembling the Shear Valve

12. Clean the area surrounding the shear valve again. You can let salt crystals and other small debris fall down. Sweep any particles laying on bottom of the 'Abacus 5' through the ventilation holes. Close the front door and remove the gloves.
13. Click/ tap the 'OK' button to inform the 'Abacus 5' software that you completed the cleaning of the shear valve. The shear valve cleaning maintenance function will check the movement and end positions of the shear valve.

## **17.7 Cleaning the Washing Head**

The washing head cleans the outer surface of the tip of the sampling needle with diluent. This is a critical function that protects the integrity of sample results by preventing contamination with previous blood samples. Salt build-up on the lower surface of the washing head may cause a malfunction during operation. In order to clean the washing head, it must be removed from the needle assembly.

The following materials are required to clean the washing head:

- A clean, soft, damp cloth
- Q-tips

- Gloves
- Water



Use extreme caution when working near the sampling needle. The sampling needle is sharp and can cause injury. Always wear gloves when performing maintenance on the 'Abacus 5' analyzer.

1. Open the front panel and secure it with the support bar.
2. The right side panel of the analyzer must be removed to access the wash head. Loosen the 2+2 screws on the front and rear of the side panel and remove the right side panel.
3. Locate the sampling needle. The washing head is positioned near the top of the sampling needle.
4. Using extreme caution, gently twist off the washing head from the needle and pull it downwards. Use a soft cloth dampened with water to clean the bottom of the wash head.




Figure 100. Removing the Washing Head

5. Use a soft cloth dampened with water or a Q-tip to clean the bottom of the wash head.




Figure 101. Cleaning the Washing Head

6. Using extreme caution with the sampling needle, gently push the wash head onto the needle. Push it up as far as it will go and lock it back in place by twisting it into the grooves of the holding rods.
7. Replace the side cover and the side cover screws. Close the front panel.

## 17.8 Software Cleaning Maintenance Functions

Daily cleaning with Diatro•Hypocleaner CC reagent during shutdown helps to keep the 'Abacus 5' analyzer clean and in good operating condition. However, blank measurements results may sometimes still be too high during daily use



or calibration. When this occurs, the 'Cleaning' and 'Hard cleaning' maintenance function are available to perform additional system cleaning on demand.

The 'Cleaning' and 'Hard cleaning' software maintenance functions are also available to perform deeper cleaning of the system tubing and pneumatic components.

If the blank measurements are still too high after the 'Cleaning' maintenance function, it will be necessary to perform the 'Hard cleaning' maintenance function for a more in-depth, thorough system cleaning. Running several blank measurements may be required after performing software maintenance cleaning functions to return blank values to within specifications.

If blank measurements do not return to normal after the 'Clean' or 'Hard clean' software maintenance functions and after running several blank measurements, please contact your Diatron certified service engineer for additional assistance.

### 17.8.1 Daily Cleaning

The Shutdown function at the end of the daily routine will request that a sample vial with 2-3 ml of Diatro•Hypocleaner CC reagent with no cap be placed into the sample rotor. The software shutdown function will clean the tubing in the measurement system.



Rubber sample caps are not designed for multiple punctures. The Diatro•Hypocleaner CC reagent will cause rubber vial caps to harden, releasing small particles that will clog the analyzer. Always remove the rubber cap when performing daily cleaning.

After this operation, you can leave the 'Abacus 5' analyzer powered on for the rest of the day. It will automatically enter a low power state if left unattended.

### 17.8.2 Cleaning the Measurement System

During normal operation, the 'Abacus 5' analyzer may report high values after a blank measurement. Repeating the blank measurement will usually bring the blank values down to within specification. If the blank values remain too high, it will be necessary to perform a cleaning maintenance function.

If the system consistently reports high blank values, the measuring chambers must be cleaned.

1. Remove the cap from a clean, dry empty sample tube. Fill the sample tube with 2 ml of deionized water and 1 ml of Diatro•Hypocleaner CC reagent. Do not replace the cap onto the sample tube.
2. Insert the sample tube with Diatro•Hypocleaner CC reagent into the sample rotor.
3. Click or tap the 'Cleaning' button on the maintenance functions panel in the 'Cleaning' group.
4. Follow the on-screen instructions.
5. Run a blank measurement and ensure the blank measurement values return to acceptable levels. Re-run the blank measurement a few times if necessary.

### 17.8.3 Extended Cleaning (Hard Cleaning)

If the blank values do not return to acceptable levels after performing the 'Cleaning' maintenance function, repeat the above steps except for clicking or tapping the 'Hard cleaning' button in step 3 rather than the 'Cleaning' button. This maintenance function will perform a longer, deep cleaning of the measurement system, buffers and chambers.

If the blank measurements do not return to normal after performing the 'Clean' and 'Hard clean' software maintenance functions and running several blank measurements, please contact your Diatron certified service engineer for additional assistance.

## 17.9 Replacing Reagents

The 'Abacus 5' analyzer is designed to be able continue running without the need to stop if a reagent container needs to be replaced. To support this feature, internal buffers keep sufficient reagents for one measurement cycle and warns the operator when a reagent is about to run out. An exclamation mark within a yellow triangle will be displayed by the reagent level indicator. Double clicking or tapping on the exclamation mark will open an explanation panel. The panel will suggest a specific solution, such as changing one of the reagent containers. Click or tap the 'Go to..' button on the explanation panel to navigate to the reagent status panel.




Figure 102. Reagent Low Indicator and Explanation Panel

To replace reagent containers, perform the following steps:

1. Prepare the new container, and place it near the place of operation
2. Open the container cap (remove cardboard flaps if necessary on 5L and 20L containers)
3. Unscrew the reagent line from the used container and place it into the new container
4. Switch the old and new containers
5. Dispose of used containers according to local regulations
6. Push the Reset button corresponding to the new reagent container to reset the reagent level software counters to 100%

	<p>Never allow the reagent line and dip tube to touch the ground or any other surface to avoid microbial contamination of the reagents and the internal pneumatic components of the 'Abacus 5' analyzer. Always transfer the reagent line from the old reagent container to the new container without coming in contact with any other surface.</p>
--	---



Always follow local laws and regulations governing the disposal of hazardous waste.



Figure 103. Reagent Replacement Panel





## 18 Reagent Locking

The reagent locking feature assures you that the 'Abacus 5' analyzer always produces quality results and operates optimally by ensuring that only genuine Diatron reagents are used in your 'Abacus 5' analyzer. The 'Abacus 5' analyzer measuring system and data analysis software are designed to work only with genuine Diatron reagents. Diatron cannot guarantee optimal performance of the 'Abacus 5' analyzer with incompatible reagents.

Each Diatro•Lyse-5P reagent container comes with a reagent key. The key is contained in an electrostatic plastic envelope located underneath the tear-out flap on the top of the reagent container cardboard box. Each key contains a license for 800 measurements.



Figure 104. Location of Reagent Key in Diatro•Lyse-5P Reagent Container

The 'Abacus 5' analyzer tracks the number of measurements performed, and warns the operator when the stored number of measurements counts down to a 50 measurements, and will stop performing sample measurements when the stored measurements are exhausted. Blank, Calibration and QC measurements are included in the overall measurement count.

The new hardware key included in your new Diatro•Lyse-5P reagent container is used to replenish the stored measurement count. To replenish the stored measurement count, replace the reagent key before replacing the Diatro•Lyse-5P reagent container.

Remove the old reagent key from the reagent lock connector located on the rear panel. Remove the new reagent key from its electrostatic envelope and insert it in the connector. Proceed with the replacement of the Diatro•Lyse-5P reagent container as usual. Clicking or tapping the 'Reset Lyse-5P' button will upload the new license for 800 measurements as well as reset the software counter for the Diatro•Lyse-5P to 100%.



Figure 105. Reagent Key Connector

The hardware-key looks similar to a USB memory device in appearance, but is not a USB device. Do not insert the reagent key into a USB device connector, Only insert it into the special connector pictured above. Used reagent keys cannot be re-used. Collect them and return them to your Diatron certified service engineer or your sales

representative for recycling, or dispose of them according to local laws or regulations governing the disposal of electronic waste materials.

## 19 The Daily Routine

In this chapter a typical daily routine is described. Your daily routine may differ from the one described here. This chapter also contains some tricks and tips. The description below assumes that the laboratory works during a single daily shift only.

The daily routine:

- Power-up and start-up the 'Abacus 5' analyzer (see [Chapter 6.2.2](#) for details)
- Check reagent container levels and start obtain new reagents from storage if necessary
- Perform the pneumatics initialization and blank measurement (tap the sampling tube icon)
- Collect the blood samples
- Ensure that the correct anticoagulant tubes (K<sub>3</sub>-EDTA) are used
- Ensure that the samples ages are within the 7 hours of processing limit
- Power on the Autosampler (if installed)
- Perform the QC procedure as scheduled
- Start the measurements in manual or automatic mode (see chapter 1 for details)
- Validate and release the measurement results according to your procedures
- Perform any necessary maintenance (see chapter 19 for details)
- Shut down the instrument (see chapter 7.3 for details)
- Clean the instrument (see chapter 5.3)

Notes:

Blank measurements:

- If you want to check whether the reagent path was contaminated during reagent container replacement, then perform 4 to 5 blank measurements before accepting the blank result. The reagents from the internal reservoirs should be fully consumed before the new reagent reaches the measurement system.
- If you receive the samples in batches and more than 2 hours have elapsed between the last sample of the previous batch and the first sample of the next batch, then it is a good practice to repeat the blank measurement.

Quality Control:

- Perform QC measurements according to your local regulations.
- Monitor the expiration date and stability after first use of the QC material.
- If there are more than one hematology instruments in your laboratory then it is also a good practice to measure the same human sample across all the instruments to compare their performance. Use multiple parallel runs to reduce the effect of random error (CV%).
- We suggest performing at least one blank measurement before starting QC measurements.

Calibration:

- Calibrate the 'Abacus 5' according to laboratory regulations.
- Calibration should be performed with every new control material lot, or if QC results or comparison of results to other instruments shows significant deviation in measurement results.

- If you perform a whole blood measurement calibration, it is recommended to run a blank measurement before the calibration procedure (see chapter 1).

Cleaning: Perform cleaning or hard cleaning if:

- High blank results are experienced
- Clogging related flags are consistently appearing with measurement results
- After service actions or service level maintenance
- Measurement results or QC results are drifting (perform hard cleaning)
- Clog formation in the tubing (first remove any plugging or contamination)
- Reagents incorrectly connected to the wrong place



## 20 Troubleshooting


Your 'Abacus 5' provides all the necessary software maintenance functions required to keep the analyzer operating in optimum working order. Despite excellent maintenance, the analyzer can still encounter problems where operator intervention is required. The 'Abacus 5' analyzer software informs the operator of situations when further help is required. These are called error messages.

The Abacus 5' analyzer provides additional information in the form of numerical codes that help the operator to identify the root cause of the error. These error codes are important for service personnel. They can interpret the messages and take necessary actions.

If you receive error messages with codes, please write them down. If possible, describe what you were doing before the error appeared. This important information helps your service engineer to help you when necessary. The software will store and re-display every error message, but it cannot describe the operation you were performing at the time.

When an error comes up, the software usually attempt to offer a solution, or will try to resolve the problem and retry the current operation. If this recovery process fails, or the problem continues to occur, please contact your service engineer.

### 20.1 Software error messages

	<p>If Windows® XP® Embedded operation system errors appear, please accept the error message and repeat the operation. If additional such error messages occur, please restart the 'Abacus 5' analyzer. At the next service visit, please let the service engineer know about the problem.</p>
--	---

**Low level error messages** provide status information about the data acquisition system and pneumatic system. Each error message includes the affected part and the area of the equipment. Check this area and part to find any irregularities like clogging, mechanical problems or blocked elements. Try to repeat the operation. The system will try to recover from the failure. If the error continues to occur, please contact your service personnel.

### 20.2 Pneumatic error messages

The analyzer software is designed to be error tolerant and recover from with minor problems. However, if some physical obstructions or extreme operating conditions cause some mechanical parts to behave unexpectedly, the system will give you an error message such as this:




Figure 106. Pneumatical Error Message Example

Clicking or tapping the 'OK' button returns and retries the current operation. Clicking or tapping the 'Details' button provide youo with further informaion abou the error. It will describe additional conditions of the failure. The numbers in the second row are important for the service engineer.




Figure 107. Detailed Error Message

## 20.3 Mechanical Problems

Most mechanical problems arise because of some physical blockage of the moving parts. Problems caused by severe blocks are usually accompanied by a grinding noise. The noise is a result of a motor not being able to move, but this does not usually result in mechanical damage to moving parts.

Typical failure situations can be traced to salt build-ups around the wash head or the shear valve. Incorrect operation during cleaning of the system, especially after long inactive periods of time can cause valves or tubes to get filled with dry salt deposits. These problems can be the source of dilutor motor noise.

In extreme tube or valve clogging situations the pressure build-up inside the dilutors can cause some tubes to be released from their connections. This symptom is not normal, but will not cause damage to the system. However, it is recommended to contact your service engineer to determine the cause and resolve the problems.

### 20.3.1 Sample Rotor (SR) Failures

#### *20.3.1.1 SR Gives Grinding Noise And / Or SW Displays SR Error Messages*

- The front panel is not closed properly
- The opening is not aligned with the sample rotor door
- There is liquid in the SR, if you can see traces of salt around the SR, contact your service engineer
- Check the washing head for contamination (blood and/or salt) as excessive salt build up can prevent the wash head from lifting up effectively and the sample rotor door can get stuck

#### *20.3.1.2 SR Error Appears During Initialization Process:*

- Check the washing head for contamination (blood and/or salt)
- Contact your service engineer

#### *20.3.1.3 The SR Does Not Turn Into The Analyzer Even With Open Front Panel*

- check the washing head for contamination (blood and salt)

### 20.3.2 Needle Mechanics, Vertical Motor (Mvert) Problems

#### *20.3.2.1 The Needle Carriage Keeps Dropping Back (Down) At Initialization*

- The wash head and surrounding area should be free of salt build-ups

- Salt build up, or thick salt layer at the bottom or on the inside can block the movement of the needle in the wash head (or the movement of the wash head around the needle)
- The through hole of the needle in the vertical carriage should be free of salt.
- Salt around the needle can damage the needle and can influence the sampling process, sampling quality or amount of sample
- The wash head position relative to the needle and to the SR is not correct
- The wash head comes down too far, and if lifted, leaves no room for the SR door to turn
- The vertical rod holding the wash head is not inserted correctly, or was not pushed up to the maximum, and the wash head is sitting in a low position

### 20.3.3 Shear Valve (SV) Related Errors

#### 20.3.3.1 SV Error At The First Startup

- the SV cannot turn:
  - Make sure that the pull tab has been removed
  - Check free movement by hand (rotate the shear valve)
- The SV is stuck:
  - Try to open the Shear Valve

#### 20.3.3.2 Grinding Noise After SV Cleaning, (After SV Reinstallation)

- The rotating part cannot reach the opto sensor position – check if the thumb screw is secured correctly and that there is no gap between the upper and lower discs

#### 20.3.3.3 SV Leakage

- The upper disc is not sitting well on the lower disc
- Open, and reseal the upper disc
- There must not be any opening left between the closing screw and the upper disc




Figure 108. Proper Seating of Shear Valve Thumbscrew

### 20.3.4 Dilutor Errors

- pinched, clogged tube around dilutor or at the rear reagent tubes
- physical obstruction (foreign material)



## 20.3.5 Priming Problems

### 20.3.5.1 *The Analyzer Would Not Prime Liquids*

- The respective reagent is out
- aspirating tube (in container):
  - has fallen off,
  - has a leak,
  - or is broken
- A reagent buffer is leaking, look for liquid below the analyzer
- Damaged tube in the system, look for leakage, or traces of liquid, and contact Service

## 20.3.6 Electronics Related Problems

### 20.3.6.1 *No Image on Display, No Backlight*

- Contact your service engineer

### 20.3.6.2 *Touch Sensitive Surface Not Working*

- Connect an external USB mouse:
  - This will allow clicking on the screen elements
  - Contact service for resolution of the problem

### 20.3.6.3 *Touch (Click) Is Inaccurate*

- The touch screen needs calibration
- Connect an external USB mouse
  - Perform the calibration under Main Menu / Maintenance / Touch screen calibration

### 20.3.6.4 *The Cursor Seems To Be Moving With Good Ratios, But In a Smaller Area*

- Calibrate the touch screen. (use an external USB mouse to access the function)
- If the problem persists, contact service

## 20.3.7 The Analyzer Does Not Power On

- Check proper power connections, and the state of the (small) switch at the power supply
- If the above fails, contact Service

## 20.3.8 I<sup>2</sup>C Errors Displayed At Startup

- Contact Service

## 20.4 Measurement results related problems

### 20.4.1 Fluctuating PLT background values

When the Abacus 5 is used under increased workload, running even 300 samples a day in 5 hours of continuous operation is possible. Despite built-in maintenance programs, the system may require intensive cleaning from time to time. If you find that blank PLT values are increasing and decreasing unexpectedly, please perform a 'Hard cleaning' maintenance function.

Go to Main Menu, and click or tap the 'Maintenance' button. Click or tap the 'Hard cleaning' button in the Cleaning section to initiate the procedure. You will need a syringe (10ml) to dose cleaning liquid into the analyzer's measuring chambers. The process is described below.

#### 20.4.2 Long, smeared population

A long, smeared population with this appearance shows that a key part of the measurement process is not operating as expected. This could be the TCU, the reagents or the sample mixing. Possible root causes are:

- The two lyse reagents are interchanged
- Check tube connections or container connections




Figure 109. Smeared Population



## 21 Accessories

The following items can be found along with the analyzer.

- Power cord
- Reagent tube kit
- Cleaning tube kit
- Waste container (20L)
- Cardboard box for waste container
- Sample tube adapter
- Operator's manual

## 22 Appendix

### 22.1 Reagent Consumption

Function	Reagent consumption per function (ml)		
	Avans Lyse 5P	Avans Diff-5P	Avans Dil-Diff
StartUp	13,6	1,7	128
Measure Blank	5,7	0,7	52
Measure 5 part	5,7	0,7	52
Measure Calibration	5,7	0,7	52
Measure QC	5,7	0,7	52
Standby	0	0	11,2
Wakeup	1	0	9
Cleaning	6,7	0,7	133
Hard cleaning	6,7	0,7	144,2
Flow cell cleaning	0	0	63
Shear Valve cleaning	9	1,2	100,2
Offline (overnight) rinsing	8	1	72
<sup>2</sup> Prime Diluent (full)	0	0	103
<sup>2</sup> Prime Lyse	7.5	0	2.5
<sup>2</sup> Prime Diff5P	0	4	0
<sup>2</sup> Prime all	22	11.5	103
<sup>1</sup> DrainDilu	0	0	120
<sup>1</sup> DrainLyse	60	0	0
<sup>1</sup> DrainDiff5P	0	60	0
<sup>1</sup> DrainAll	60	60	180
FillUp	22	12	103
Shutdown	7,7	0,7	153,2
Note 1 – If internal chambers are filled (volume of chambers: 50 ml each, diluent 2x50 + 50 sheath)			
Note 2 – Prime functions fill internal chambers. The following pneumatic function uses reagents from the inside buffer while the inside buffer is continuously refilled to a steady level.			

Table 22. Reagent Consumption

## 22.2 Display Ranges

Parameter	Range	Selectable Dimensions
WBC	0.00 – 120.00 * 10 <sup>3</sup> /μl or 10 <sup>9</sup> /l	cells/μl, cells/l
NEU	0.00 – 120.00 * 10 <sup>3</sup> /μl or 10 <sup>9</sup> /l	cells/μl, cells/l
LYM	0.00 – 120.00 * 10 <sup>3</sup> /μl or 10 <sup>9</sup> /l	cells/μl, cells/l
MON	0.00 – 120.00 * 10 <sup>3</sup> /μl or 10 <sup>9</sup> /l	cells/μl, cells/l
EO	0.00 – 120.00 * 10 <sup>3</sup> /μl or 10 <sup>9</sup> /l	cells/μl, cells/l
BAS	0.00 – 120.00 * 10 <sup>3</sup> /μl or 10 <sup>9</sup> /l	cells/μl, cells/l
NEU%	0 -100%	Not selectable
LYM%	0 -100%	Not selectable
MON%	0 -100%	Not selectable
EO%	0 -100%	Not selectable
BAS%	0 -10%	Not selectable
RBC	0.00 – 10.00 * 10 <sup>6</sup> /μl or 10 <sup>12</sup> /l	cells/μl, cells/l
HGB	0 – 30 g/dl or 0 – 300 g/l	g/dl, g/l
HCT	0 – 60 %	percentage
MCV	0 – 250 fl	Not selectable
MCH	(no restrictions) pg	Not selectable
MCHC	(no restrictions) g/dl or g/l	g/dl, g/l
RDWsd	(no restrictions) fl	Not selectable
RDWcv	(no restrictions) %	Not selectable
PLT	0 – 1200 * 10 <sup>3</sup> /μl or 10 <sup>9</sup> /l	cells/μl, cells/l
PCT	(no restrictions) %	Not selectable
MPV	3.0 – (no upper restriction) fl	Not selectable
PDWsd	(no restrictions) fl	Not selectable
PDWcv	(no restrictions) %	Not selectable

Table 23. Display Ranges

## 22.3 Fluidic System




Figure 110. Fluidic System Diagram

## 22.4 Printed Report Formats





Demo										12/07/2011 17:22	
Demo1					Demo4						
Demo2					Demo5						
Demo3					Demo6						
Sample ID		24			Patient ID		1				
Date		09/04/2010 11:50			Name		Demo				
Mode		Human			Date of birth		01/01/1999				
Parameter		Result			Limit		Gender		Male		
WBC	10.41	+	10 <sup>3</sup> /μL	5.00	-	10.00					
NEU	5.47		10 <sup>3</sup> /μL	2.00	-	7.50					
LYM	3.29		10 <sup>3</sup> /μL	1.30	-	4.00					
MON	1.38	+	10 <sup>3</sup> /μL	0.15	-	0.70					
EO	0.18		10 <sup>3</sup> /μL	0.00	-	0.50					
BAS	0.09		10 <sup>3</sup> /μL	0.00	-	0.15					
NEU%	52.5		%	40.0	-	75.0					
LYM%	31.6		%	21.0	-	40.0					
MON%	13.3	+	%	3.0	-	7.0					
EO%	1.7		%	0.0	-	5.0					
BAS%	0.9		%	0.0	-	1.5					
RBC	3.75	-	10 <sup>6</sup> /μL	4.00	-	5.50					
HGB	10.1	-	g/dL	12.0	-	17.4					
HCT	33.7	-	%	36.0	-	52.0					
MCV	89.7		fL	76.0	-	96.0					
MCH	26.8	-	pg	27.0	-	32.0					
MCHC	29.9	-	g/dL	30.0	-	35.0					
RDWsd	52.0	+	fL	20.0	-	42.0					
RDWcv	11.8		%	0.0	-	16.0					
PLT	365		10 <sup>3</sup> /μL	150	-	400					
PCT	0.49		%		-						
MPV	13.5		fL	8.0	-	15.0					
PDWsd	25.4		fL		-						
PDWcv	34.2		%		-						
Warning flags											
Morphological flags											
Interpretive flags		Leukocytosis?, Anemia?, Hypochromic?									


Figure 111. Detailed Printout Format

Lab header 2			Lab header 5					
Lab header 3			Lab header 6					
Lab header 4			Lab header 7					
Sample ID	WBC	NEU	LYM	MON	EO	BAS		
Mode		NEU%	LYM%	MON%	EO%	BAS%		
Patient ID:	RBC	HGB	HCT	MCV	MCH	MCHC	RDWsd	RDWcv
Date	PLT	PCT	MPV	PDWsd	PDWcv			
AS_10_A2	6.25	4.53	1.21	0.33	-- 0.09	0.09		
Human		72.5	19.3	5.3	1.5	1.4		
1	* 8.89	+ 230	** ---	* 79.7	* 25.8	** ---	* 62.2	* 23.0
14/06/2011 15:51	-- 22	0.02	9.2	24.9	40.1			
20418140022	+ 11.97	+ 8.68	1.99	+ 1.03	- 0.14	0.13		
Human		72.5	16.6	8.6	1.2	1.1		
1	4.76	- 112	40.3	84.6	- 23.5	++ 277	55.9	+ 19.4
14/06/2011 15:53	148	0.21	+ 13.9	37.4	38.3			
AS_10_A2	6.25	4.53	1.21	0.33	-- 0.09	0.09		
Human		72.5	19.3	5.3	1.5	1.4		
1	* 8.89	+ 230	** ---	* 79.7	* 25.8	** ---	* 62.2	* 23.0
14/06/2011 15:57	-- 22	0.02	9.2	24.9	40.1			
AS_8_A1	+ 15.36	+ 12.79	1.43	0.26	0.84	0.03		
Human		+ 83.3	- 9.3	- 1.7	5.5	0.2		
1	+ 7.75	137	48.5	- 62.5	- 17.6	++ 282	54.0	+ 27.1
14/06/2011 15:59	+ 473	0.45	9.6	27.6	39.4			
00025-110523	9.49	6.13	2.10	+ 0.85	0.33	0.08		
Human		64.6	22.1	9.0	3.5	0.8		
1	4.93	149	44.0	89.3	30.3	++ 339	54.5	+ 17.6
17/06/2011 11:08	247	0.30	+ 12.2	28.0	38.8			

Figure 112. Listed Printout Format



## 22.5 Parameter Calculation

Parameter	Description	Units	Calculation
WBC	Leukocyte (white blood cells) concentration	$10^3/\mu\text{l}$ or $10^9/\text{l}$	$\text{WBC} = \text{WBC}_{\text{CAL}} \times \text{WBC}_{\text{MEASURED}}$
NEU LYM MON EO BAS	Neutrophil, Lymphocyte, Monocyte, Eosinophil and Basophil concentration	$10^3/\mu\text{l}$ or $10^9/\text{l}$	Calculated from differential percentage and WBC $\text{NEU}\% \times \text{WBC}$ $\text{LYM}\% \times \text{WBC}$ $\text{MON}\% \times \text{WBC}$ $\text{EO}\% \times \text{WBC}$ $\text{BAS}\% \times \text{WBC}$
NEU% LYM% MON% EO%	Neutrophil, Lymphocyte, Monocyte and Eosinophil five part differential percentage	%	Relative value are derived from the optical 4 diff scatter populations
BAS%	Basophil five part differential percentage	%	Relative value is calculated from the basophil optical scatter populations
RBC	Erythrocyte (red blood cells) concentration	$10^6/\mu\text{l}$ or $10^{12}/\text{l}$	$\text{RBC} = \text{RBC}_{\text{CAL}} \times \text{RBC}_{\text{MEASURED}}$
HGB	Hemoglobin concentration	g/dl or g/l	Measured photometrically at a wavelength of 568 nm; blank measurement is performed per cycle on diluent and lyse mixture $\text{HGB} = \text{HGB}_{\text{CAL}} \times (\text{HGB}_{\text{MEASURED}} - \text{HGB}_{\text{BLANK}})$
HCT	Hematocrit	%	Calculated from RBC and MCV $\text{HCT} = \text{RBC} \times \text{MCV} \times 100$
MCV	Mean Corpuscular Volume	fl	Average volume of individual erythrocytes derived from the RBC histogram
MCH	Mean Corpuscular Hemoglobin	pg	Average hemoglobin content of erythrocytes, calculated from RBC and HGB values. $\text{MCH} = \text{HGB} / \text{RBC}$
MCHC	Mean Corpuscular Hemoglobin Concentration	g/dl or g/l	Calculated from HGB and HCT $\text{MCHC} = \text{HGB} / (\text{HCT} / 100)$
RDWsd	Red Cell Distribution Width	fl	Calculated from RDWcv and the distribution width of the erythrocyte histogram at 20% of peak $\text{RDWsd} = \text{RDWcv} \times (P_2 + P_1) / 0.56$
RDWcv	Red Cell Distribution Width CV	%	 <p>Derived from the distribution width of the erythrocyte histogram at 20% of peak  <math>\text{RDWcv} = \text{RDW}_{\text{CAL}} \times 0.56 \times (P_2 - P_1) / (P_2 + P_1)</math>  Factor of 0.56 converts measurement at 60% of peak to 20% of peak </p>
PLT	Platelet (thrombocyte) concentration	$10^3/\mu\text{l}$ or $10^9/\text{l}$	$\text{PLT} = \text{PLT}_{\text{CAL}} \times \text{PLT}_{\text{MEASURED}}$
PCT	Thrombocrit	%	Calculated from PLT and MPV $\text{PCT} = \text{PLT} \times \text{MPV} \times 100$
MPV	Mean Platelet Volume	fl	Average volume of individual platelets derived from the PLT histogram

Parameter	Description	Units	Calculation
PDWsd	Platelet Distribution Width	fl	Calculated from RDWcv and the distribution width of the platelet histogram at 20% of peak $PDWsd = PDWcv \times (P_2 + P_1) / 0.56$
PDWcv	Platelet Distribution Width CV	%	Derived from same histogram as above $PDWcv = PDW_{CAL} \times 0.56 \times (P_2 - P_1) / (P_2 + P_1)$ Factor of 0.56 converts measurement at 60% of peak to 20% of peak

Table 24. Parameter Calculation

## 22.6 Specifications

Item	Specification			
Dimensions	Analyzer Only	H	20 <sup>1</sup> / <sub>2</sub> inches	520 mm
		W	16 <sup>1</sup> / <sub>8</sub> inches	410 mm
		D	19 <sup>3</sup> / <sub>8</sub> inches	493 mm
	Analyzer with Autosampler	H	20 <sup>1</sup> / <sub>2</sub> inches	520 mm
		W	26 <sup>3</sup> / <sub>8</sub> inches	671 mm
		D	19 <sup>3</sup> / <sub>8</sub> inches	493 mm
Weight	Analyzer Only	79.4 lbs.		36 kg
	Analyzer with Autosampler	103 lbs.		46.7 kg
Sample volume	Closed and open vial mode: 110 μl			
Sample type	Human whole blood (K <sub>3</sub> -EDTA anticoagulant)			
Tube Identification	By means of the front panel keyboard (enter ID) By means of the barcode labels (manual and/or auto-sampler)			
Sampling method	Ceramic shear valve with 3 separate primary loops			
Measured parameters	CBC+5DIFF mode (24 parameters): WBC, LYM, MON, NEU, EOS, BAS, LYM%, MON%, NEU%, EOS%,BAS%,RBC, HCT, MCV, RDWsd, RDWcv HGB, MCH, MCHC, PLT, MPV, PDWsd, PDWcv, PCT			
Throughput	60 tests/hour			
Measurement method	Volumetric impedance change for WBC, RBC, PLT Light scattering BASO measurement Light scattering 4-diff measurement: LYM, MON, NEU, EOS Spectrophotometry for HGB			
Aperture diameter	WBC: 80 μm, RBC/PLT: 70 μm			
Aperture length	WBC: 80 μm, RBC/PLT: 70 μm			
HGB measurement	Light source: green LED with 568 nm wavelength Detector: light to frequency converter			
Optical measurement	Light source: semiconductor laser diode with 650 nm wavelength and 7mW Beam divergence: ~1.5° Class 1 Laser Product (Class IIIB laser module if the protective housing is opened, and interlocks defeated) IEC 60825-1: Ed. 2(2007) Quartz flow cell with hydro-dynamic focusing Detector: fiber optic coupled PIN Si photodiodes Internal safety interlock			
Auto-alignment system	Horizontal calibration of laser beam position Fine calibration with 5 μm polystyrene microsphere calibration material			
Reagents	Diatro•Dil-5P (20 liter) Diatro•Diff-5P (1 liter) Diatro•Lyse-5P (5 liter) Diatro•Hypocleaner CC (100 ml) (Emergency cleaner)			
Dilution ratios	WBC/BAS dilution - 1:170 RBC/PLT dilution - 1:21,250 4 Diff dilution - 1:50			
Sheath fluid	Diluent			
Control material	Diacon 5, Manufacturer: Diatron			
Quality Control	QC diagrams, QC database			

Item	Specification
Flagging	Pathological (diagnostic) flags Lab limits (normal ranges) Reagents alert (3 measurement prealert-online reagent replacement) Instrumental alerts, internal buffer for reagents
Calibration	Manual and SW supported automatic mode
Languages available	English menu and support for other languages
Software upgrade	Via USB
Data storage capacity	100,000 records including flags, scatter- and histograms
Data processing	Intel Atom 1.6 GHz
Data store	Windows® XP® Embedded
Display	800 x 600 color graphic LCD, portrait layout
External printing	Via USB port, any Windows® XP® compatible printer
External keyboard	Via PS/2 or USB
Bar-Code reader	Optional Manual bar-code reader via USB Built in Bar-Code in the Auto-sampler
Peripheral ports	4 x USB 2.0, Ethernet, PS/2
Power requirements	Power supply input: 100-127VAC/200-240 VAC; 47Hz to 63Hz Power Consumption: maximum 400 VA
Main fuse	F 10A H 250V
Operating temperature	15-30°C (59-86°F)
Maximum relative humidity	80%

Table 25. Specifications

## 22.7 Performance Characteristics

### 22.7.1 Precision

Device repeatability and reproducibility are defined as a standard deviation (SD) or as a coefficient of variation (CV), whichever is greater. Specifications are provided for primary and derived parameters only.

Parameter	Repeatability		Reproducibility/Precision		Units
	SD	CV%	SD	CV%	
WBC	< 0.18	< 2.7 %	< 0.20	< 3.4%	10 <sup>3</sup> /μl
NEU%	< 3.50		< 3.50		%
LYM%	< 3.10		< 3.10		%
MON%	< 2.00		< 2.00		%
EOS%	< 2.00		< 2.00		%
BAS%	< 0.50		< 0.50		%
RBC	< 0.11	< 1.7 %	< 0.13	< 2.0%	10 <sup>6</sup> /μl
HGB	< 0.20	< 2.0 %	< 0.22	< 2.4%	g/dl
MCV	< 1.0	< 1.7 %	< 1.20	< 2.0%	fL
RDWcv	< 0.4	< 2.5%	< 0.45	< 3.0%	%
PLT	< 23	< 6.0 %	< 27	< 7.0%	10 <sup>3</sup> /μl
MPV	< 0.45	< 8.7 %	< 0.50	< 10.0%	fL

Table 26. Precision Performance

### 22.7.2 Accuracy

Device accuracy is defined as a maximum predicted difference (bias) between this analyzer and an alternate method. The difference must be less than the absolute or percent difference criteria in the table below over an expected range. Specifications are provided for primary and derived parameters only.

Parameter	Criteria		Range		Units
	Absolute	Percent	Low	High	
WBC	0.30	6.00%	0.00	100.00	10 <sup>3</sup> /μl
NEU%	3.00	10.00%	0.00	100.00	%
LYM%	3.00	10.00%	0.00	100.00	%
MON%	3.00	10.00%	0.00	40.00	%
EOS%	1.50	10.00%	0.00	30.00	%
BAS%	0.50	40.00%	0.00	5.00	%
RBC	0.15	6.00%	0.00	8.00	10 <sup>6</sup> /μl
HGB	0.30	6.00%	1.00	25.00	g/dl
MCV	1.00	6.00%	50.00	120.00	fL
RDWcv	0.50	6.00%	8.00	20.00	%
PLT	15.00	8.00%	0.00	1000.00	10 <sup>3</sup> /μl
MPV	0.80	10.00%	5.00	20.00	fL

Table 27. Accuracy Performance

### 22.7.3 Linearity

Linearity is defined for the primary measured parameters.

Parameter	Coefficient of Determination (R <sup>2</sup> )	Nonlinearity Absolute Error	Nonlinearity Relative Error	Linearity Range Low	Linearity Range High	Units
WBC	> 0.95	< 0.80	< 3.0 %	0.2	100.0	10 <sup>3</sup> /μl
RBC	> 0.95	< 0.20	< 3.0 %	0.36	7.19	10 <sup>6</sup> /μl
HGB	> 0.95	< 0.27	< 3.0 %	1.1	22.2	g/dl
PLT	> 0.95	< 35	< 3.0 %	15	1000	10 <sup>3</sup> /μl

Table 28. Linearity Performance

## 22.7.4 Carryover

Parameter	Allowable Carryover	Units
WBC	< 1.00%	10 <sup>3</sup> /μl
RBC	< 0.50%	10 <sup>6</sup> /μl
HGB	< 0.80%	g/dl
PLT	< 1.00%	10 <sup>3</sup> /μl

Table 29. Carryover Performance

## 22.7.5 Sample Stability

Sample stability is characterized by the absolute and percent difference from the baseline values measured after 30 minutes of venipuncture for the parameters WBC, LYM%, MID%, GRA%, RBC, HGB, HCT, MCV, RDW and PLT. For MPV, the baseline shall be 2 hours after venipuncture. Either the absolute or percent difference must be under either the absolute or the percent difference criteria indicated in the table in section 22.7.2. The minimum expected stability is 7 hours from venipuncture.

## 22.7.6 Mode to Mode

Parameter	n	Correlation	Intercept	Slope	Mean		Units
					Closed Vial	Open Vial	
WBC	40	0.9997	-0.055	1.008	13.152	13.200	10 <sup>3</sup> /μl
NEU%	40	0.9954	-0.573	1.006	58.025	57.868	%
LYM%	40	0.9963	-0.025	1.004	30.696	30.661	%
MON%	40	0.9750	0.202	0.970	8.780	8.767	%
EOS%	39	0.8953	0.323	0.927	1.356	1.659	%
BAS%	40	0.6944	0.454	0.547	0.959	0.980	%
RBC	40	0.9989	-0.014	1.003	4.274	4.276	10 <sup>6</sup> /μl
HGB	40	0.9956	-0.105	1.016	12.104	12.191	g/dl
MCV	40	0.9957	-1.675	1.020	86.330	86.404	fl
RDWcv	39	0.9925	-0.212	1.015	16.263	16.327	%
PLT	39	0.9976	4.401	0.997	246.177	249.886	10 <sup>3</sup> /μl
MPV	40	0.9315	0.290	0.953	7.830	7.763	fl

Table 30. Closed vs. Open Vial Mode Performance

Parameter	n	Correlation	Intercept	Slope	Mean		Units
					Closed Vial	Open Vial	
WBC	40	0.9998	-0.050	1.015	13.152	13.055	10 <sup>3</sup> /μl
NEU%	40	0.9967	-0.391	1.006	58.025	57.794	%
LYM%	40	0.9957	0.104	1.013	30.696	30.263	%
MON%	40	0.9654	0.551	0.887	8.780	9.314	%
EOS%	39	0.9224	0.219	0.917	1.356	1.541	%
BAS%	40	0.7494	0.400	0.587	0.959	0.986	%
RBC	40	0.9989	-0.056	1.017	4.274	4.259	10 <sup>6</sup> /μl
HGB	40	0.9985	0.054	1.007	12.104	12.055	g/dl
MCV	40	0.9977	-0.489	1.006	86.330	86.410	fl
RDWcv	39	0.9923	-0.354	1.025	16.263	16.140	%
PLT	39	0.9963	4.795	0.997	246.177	245.949	10 <sup>3</sup> /μl
MPV	40	0.9236	0.330	0.947	7.830	7.880	fl

Table 31. Closed vs. Autosampler Mode Performance

## 22.7.7 Reference Ranges

Evaluation of reference ranges for the Abacus 5 hematology analyzer was conducted using a data set of 240 normal human whole blood samples. Of the 240 samples, 120 were from female patients and 120 were from male patients  $\geq 22$  years of age. Diatron recommends each laboratory establish its own reference ranges.

Parameter	Unit	Lower limit	Upper limit
WBC	$10^3/\mu\text{l}$	4.50	10.37
LYM	%	14.76	45.40
MON	%	2.91	12.1
NEU	%	42.90	78.10
EOS	%	0.10	7.00
BAS	%	0.15	1.60
LYM	$10^3/\mu\text{l}$	1.08	3.17
MON	$10^3/\mu\text{l}$	0.20	0.91
NEU	$10^3/\mu\text{l}$	2.43	7.42
EOS	$10^3/\mu\text{l}$	0.01	0.53
BAS	$10^3/\mu\text{l}$	0.01	0.13
RBC	$10^6/\mu\text{l}$	3.86 F-3.86 M-3.91	5.62 F-5.18 M-5.62
HGB	g/dl	11.8 F-11.8 M-12.0	16.9 F-15.1 M-16.9
HCT	%	35.5 F-35.5 M-36.2	48.0 F-46.5 M-48.0
MCV	fl	81.6	97.7
MCH	pg	26.8	33.8
MCHC	g/dl	31.1	35.5
RDWcv	%	12.8	16.8
PLT	$10^3/\text{l}$	151	361
MPV	fl	6.1	14.1

Table 32. Reference Ranges

## 22.7.8 Background Limits

Parameter	Unit	Blank high limit
WBC	$10^3 \text{ cells}/\mu\text{l}$	> 0.20
HGB	g/dl	> 1.0
RBC	$10^6 \text{ cells}/\mu\text{l}$	> 0.05
PLT	$10^3 \text{ cells}/\mu\text{l}$	> 15

Table 33. Background Limits

## 22.7.9 Analytical Measurement Range

Parameter	AMR		Units
	Low	High	
WBC	0.21	100.00	$10^3 \text{ cells}/\mu\text{l}$
NEU%	0.0	100.0	%
LYM%	0.0	100.0	%
MON%	0.0	100.0	%
EOS%	0.0	100.0	%
BAS%	0.0	10.0	%
NEU	0.00	100.00	$10^3 \text{ cells}/\mu\text{l}$
LYM	0.00	100.00	$10^3 \text{ cells}/\mu\text{l}$
MON	0.00	100.00	$10^3 \text{ cells}/\mu\text{l}$
EOS	0.00	100.00	$10^3 \text{ cells}/\mu\text{l}$
BAS	0.00	10.00	$10^3 \text{ cells}/\mu\text{l}$
RBC	0.36	7.19	$10^6 \text{ cells}/\mu\text{l}$
HGB	1.1	22.2	g/dl
HCT	NA	NA	%
MCV	50	120	fl
MCH	NA	NA	pg

Parameter	AMR		Units
	Low	High	
MCHC	NA	NA	g/dl
RDWcv	8	40	%
RDWsd	10	250	fl
PLT	16	1000	10 <sup>3</sup> cells/ $\mu$ l
MPV	4.0	30.0	fl

Table 34. Analytical Measuring Range

#### 22.7.10 Interfering Substances

The following substances can interfere with parameter measurement and alternate measurement procedures may be required.

Parameter	Interference
WBC	> 5 NRBCs/100 WBCs , PLT clumps/large PLTs
RBC	WBC Count > 75.0 x10 <sup>3</sup> / $\mu$ L
MCV	WBC Count > 75.0 x10 <sup>3</sup> / $\mu$ L
PLT	PLT clumps/largePLTs
Hemoglobin	WBC count > 75.0 x10 <sup>3</sup> / $\mu$ L , Lipids > 280 mg/dL
Differential	> 5 NRBC's/100 WBC's , PLT clumps/largePLTs



## 22.8 Reagent System

### 22.8.1 Diluent

Diatro•Dil – 5P diluent      Isotonic solution, used to dilute whole blood and quantitative and qualitative determination of RBC, WBC, PLT and HGB concentration

### 22.8.2 Lyse Reagent 1

Diatro•Lyse - 5P      Reagent for stromatolysis of RBC and quantitative determination of WBC 5-part differentiation (LYM, MON, NEU, EOS, BAS) and HGB concentration measurement of human blood.

### 22.8.3 Lyse Reagent 2

Diatro•Diff - 5P      Quantitative determination of WBC, leukocyte four-part differentiation (LYM, MON, NEU, EOS).

### 22.8.4 Cleaner

Diatro•Hypocleaner CC      Capillaries, tubing and chambers, removing blood component precipitates.

## 22.9 Tab File Format

This file format is composed of a header line and consecutive lines containing records selected for saving. Each parameter is selected by the <TAB> (08h) character, allowing easy interfacing or importing into most data processing applications, like Microsoft® Excel. Scatter diagrams and histograms, flags are NOT saved in this format.

The header line contains the same column names you can see in the analyzer database.

Data values are saved as they are stored in the analyzer database, alphanumerical or numerical. Decimal separator is defined by the actual operating language of the software.

### Example “tab separated file”

```
Sample ID<TAB>Patient<TAB>Analysis  
time<TAB>Mode<TAB>WBC<TAB>NEU<TAB>LYM<TAB>MON<TAB>EO<TAB>BAS<TAB>NEU%<TAB>LYM%<T  
AB>MON%<TAB>EO%<TAB>BAS%<TAB>RBC<TAB>HGB<TAB>HCT<TAB>MCV<TAB>MCH<TAB>MCHC<TAB>RD  
Wsd<TAB>RDWcv<TAB>PLT<TAB>PCT<TAB>MPV<TAB>PDWsd<TAB>PDWcv<TAB>Warnings<TAB>Tray<  
TAB>User<TAB>SEQ<TAB>
```

```
27<TAB>1<TAB>4/9/2010 11:50:05  
AM<TAB>Human<TAB>4.85<TAB>2.66<TAB>1.32<TAB>0.48<TAB>0.33<TAB>0.06<TAB>54.8<TAB>  
27.3<TAB>9.9<TAB>6.8<TAB>1.2<TAB>3.28<TAB>95<TAB>31.6<TAB>96.2<TAB>29<TAB>302<TA  
B>50.2<TAB>11<TAB>279<TAB>0.48<TAB>17.2<TAB>33.4<TAB>34.2<TAB><TAB><TAB><TAB>7<T  
AB>
```

```
26<TAB>1<TAB>4/9/2010 11:50:25  
AM<TAB>Human<TAB>6.99<TAB>4.26<TAB>1.92<TAB>0.64<TAB>0.12<TAB>0.05<TAB>61<TAB>27  
.4<TAB>9.2<TAB>1.7<TAB>0.7<TAB>3.44<TAB>94<TAB>31.7<TAB>92.3<TAB>27.4<TAB>297<TA  
B>56.1<TAB>12.4<TAB>331<TAB>0.42<TAB>12.6<TAB>26.1<TAB>34.9<TAB><TAB><TAB><TAB>8  
<TAB>
```

```
25<TAB>1<TAB>4/9/2010 11:50:40  
AM<TAB>Human<TAB>6.85<TAB>3.8<TAB>2.16<TAB>0.6<TAB>0.27<TAB>0.02<TAB>55.5<TAB>31  
.5<TAB>8.8<TAB>3.9<TAB>0.3<TAB>4.19<TAB>107<TAB>36.7<TAB>87.5<TAB>25.5<TAB>292<T  
AB>49.4<TAB>11.8<TAB>378<TAB>0.49<TAB>12.9<TAB>25.2<TAB>34.8<TAB><TAB><TAB><TAB>  
9<TAB>
```

```
24<TAB>1<TAB>4/9/2010 11:50:56  
AM<TAB>Human<TAB>10.41<TAB>5.47<TAB>3.29<TAB>1.38<TAB>0.18<TAB>0.09<TAB>52.5<TAB>  
>31.6<TAB>13.3<TAB>1.7<TAB>0.9<TAB>3.75<TAB>101<TAB>33.7<TAB>89.7<TAB>26.8<TAB>2  
99<TAB>52<TAB>11.8<TAB>365<TAB>0.49<TAB>13.5<TAB>25.4<TAB>34.2<TAB><TAB><TAB><TA  
B>10<TAB>
```

## 23 Index

### A

aperture, 26, 32, 33

### B

BASO, 28, 70, 75, 79, 80, 81

blank, 27, 40, 46, 48, 49, 51, 56, 57, 59, 61, 64, 66, 72, 74, 78, 79, 86, 91, 96, 121, 128, 133, 134, 147

### C

calibration, 23, 36, 44, 47, 48, 56, 90, 91, 92, 93, 134, 139

cleaning, 1, 30, 44, 46, 50, 59, 70, 78, 79, 80, 121, 123, 124, 126, 128, 134, 137, 139, 140, 142

clogging, 78, 79, 80, 136, 137

### D

diluent, 147

Diluent, 24, 32, 142

dilutor, 137, 138

### E

EDTA, 54, 55, 133

EOS, 22, 24, 26, 27, 81, 143, 154

### F

flow cell, 27, 78, 79, 80

### H

hemoglobin, 147

HGB, 22, 24, 26, 33, 43, 48, 70, 79, 90, 108, 143, 147, 154

hypochlorite, 25

### I

impedance, 26, 32, 70, 75

### K

keyboard, 14, 16, 31, 36, 37, 40, 43, 57, 73, 86, 105

### L

laser, 22, 27, 28, 33

LYM, 22, 24, 26, 27, 74, 79, 80, 81, 108, 143, 154

Lyse, 24, 31, 32, 132, 142, 154

### M

maintenance, 12, 13, 20, 30, 36, 44, 46, 78, 118, 120, 121, 133, 134, 139

MON, 22, 24, 26, 27, 74, 78, 79, 81, 108, 143, 154

## N

needle, 12, 30, 33, 48, 54, 58, 70, 126, 127, 138  
NEU, 22, 24, 26, 27, 74, 78, 80, 108, 143, 154  
noise, 80, 121, 137  
normal ranges, 43, 56, 58, 64, 66, 73, 105, 109

## P

PLT, 22, 23, 24, 26, 43, 48, 70, 75, 76, 78, 79, 81, 90, 108, 139, 143, 147, 154  
pressure, 36, 137  
printer, 14, 41, 43, 69, 108, 109

## Q

QC, 22, 44, 47, 48, 56, 57, 59, 82, 94, 95, 96, 97, 132, 133, 134, 142

## R

RBC, 22, 24, 26, 32, 43, 44, 48, 70, 75, 76, 78, 79, 80, 81, 90, 108, 143, 147, 154  
reagent, 1, 13, 24, 25, 26, 31, 32, 39, 42, 44, 46, 50, 51, 52, 53, 78, 106, 117, 129, 132, 133, 138, 139  
reagent consumption, 1, 51, 53, 117

## S

sample ID, 37, 40, 56, 57, 59, 60, 63, 66, 73  
sampling needle, 12, 30, 54, 121, 127  
shear valve, 30, 46, 48, 70, 124, 137, 138  
shear-valve, 123, 124, 126  
standby, 31, 36, 39, 43, 47, 53, 107

## T

touch screen, 22, 30, 36, 83, 139  
tubing, 25, 31, 46, 48, 51, 52, 134, 154

## U

USB, 22, 31, 36, 57, 87, 88, 132, 139

## W

wash head, 121, 126, 127, 137, 138  
waste, 12, 17, 20, 31, 32, 39, 43, 44, 46, 51, 52, 106, 130, 142  
waste container, 20, 39, 43, 46, 51, 52, 106, 142  
WBC, 22, 24, 26, 27, 33, 43, 44, 48, 70, 78, 79, 80, 81, 90, 108, 143, 147, 154

